



**SANILUDO  
MEDIC**



# INSTRUMENTACIÓN QUIRÚRGICA DE ODONTOLOGÍA



---

[www.saniludomedic.com](http://www.saniludomedic.com)

(+593) 994766871 - (+593) 986212125

[fcastro@saniludomedic.com](mailto:fcastro@saniludomedic.com)

Urb. Milann Ferrara MZ.6 V.17



*Together*

*Conquer the world*





## CONTENTS

SR. NO.	PRODUCT CODE	PRODUCT NAME	PAGE NO.
1	T-7101	INFUSION SET UMAFLOW (LUER LOCK)	3
2	T-7102	INFUSION SET UMAFLOW (TOP)	3
3	T-7103	INFUSION SET UMAFLOW (SAFE-AN-FLOW)	4
4	T-7117	INFUSION SET UMAFLOW (PAPER POUCH)	4
5	T-7104	INFUSION SET UMAFLOW (1-UP)	5
6	T-7105	INFUSION SET UMAFLOW (MICRODRIP)	5
7	T-7106	INFUSION SET UMAFLOW (MICRODRIP)	6
8	T-7107	MEASURE VOLUME SET UMADRIP 150 ML.	6
9	T-7108	MEASURE VOLUME SET UMADRIP 110 ML.	7
10	T-7109	BLOOD ADMINISTRATION SET (UMAFLOW)	7
11	T-7110	BLOOD ADMINISTRATION SET (UMAFLOW)	8
12	T-7112	SCALP VEIN SET (UMAFLOW)	8
13	T-7114	HIGH PRESSURE MONITORING LINE (UMAFLOW)	9
14	T-7116	EXTENSION LINE WITH 3-WAY STOP COCK	9
15	T-7111	PERITONEAL DIALYSIS SET - UMAFLOW (Y SET)	10
16	U-7208	URINE COLLECTING BAG (TINY)	10
17	U-7209	URINE COLLECTING BAG (UMAFLOW)	11
18	U-7202	URINE COLLECTING BAG (UROFILL)	11
19	U-7203	URINE COLLECTING BAG (UROSOFT)	12
20	U-7204	URINE COLLECTING BAG (ROYAL-19)	12
21	U-7212	URINE COLLECTING BAG WITH BOTTOM OUTLET (ROYAL-5)	13
22	U-7211	URINE BAG-LEG BAG SET (ROYALCARE)	13
23	U-7210	URINE COLLECTING BAG WITH BOTTOM OUTLET (URO MEASURE)	14
24	U-7205	NELATON CATHETER (UMAFLOW)	14
25	U-7206	MALE EXTERNAL CATHETER (URICATH)	15
26	U-7207	T.U.R. SET (UMAFLOW)	15
27	G-7301	RYLE'S TUBE (UMAFLOW)	16
28	G-7304	RYLE'S TUBE (UMALINE)	16
29	G-7302	INFANT FEEDING TUBE (UMAFLOW)	17
30	G-7303	ABDOMINAL DRAINAGE KIT (UMA ADK)	17
31	G-7305	COLOSTOMY KIT (UMAFLOW)	18
32	A-7401	SUCTION CATHETER (UMAFLOW)	18
33	A-7407	SUCTION CATHETER - UMAFLOW (K-60,K-61)	19

SR. NO.	PRODUCT CODE	PRODUCT NAME	PAGE NO.
34	A-7409	SUCTION CATHETER THUMB CONTROL (UMAFLOW)	19
35	A-7402	URETHERAL CATHETER (UMAFLOW)	20
36	A-7408	URETHERAL CATHETER-UMAFLOW (K-90,K-91)	20
37	A-7403	NASAL OXYGEN CATHETER (UMAFLOW)	21
38	A-7404	TWIN BORE NASAL OXYGEN SET (UMAFLOW)	21
39	A-7405	OXYGEN FACE MASK (UMAFLOW)	22
40	A-7415	NEBULIZER KIT (UMAFLOW)	22
41	A-7413	HIGH CONCENTRATION MASK (UMAFLOW)	23
42	A-7414	ASEPTO PUMP (UMAFLOW)	23
43	A-7406	GUEDEL AIRWAY (UMAFLOW)	24
44	A-7410	TRACHEOSTOMY TUBE - UMAFLOW (PLAIN)	24
45	A-7411	ENDOTRACHEAL TUBE - UMAFLOW (PLAIN) ORAL/NASAL	25
46	A-7412	RESPIRATORY EXERCISER (UMACISER)	25
47	S-7501	CHEST DRAINAGE CATHETER (UMAFLOW)	26
48	S-7504	WATER SEALED CHEST DRAINAGE BAG (UMADRAIN)	26
49	S-7502	CLOSED WOUND SUCTION UNIT (UMA VAC SET)	27
50	S-7503	CLOSED WOUND SUCTION UNIT (MINOR VAC SET)	27
51	S-7505	VACCUM SUCK SUCTION SET (UMAFLOW)	28
52	S-7506	CONNECTING TUBE (UMAFLOW)	28
53	S-7507	KARMAN CANNULA (UMAFLOW)	29
54	S-7509	UMBILICAL CANNULA (UMAFLOW)	29
55	S-7510	C.V.P. MANOMETER (UMALOW)	30
56	M-7601	INFANT MUCUS EXTRACTOR (UMAFLOW)	30
57	M-7602	UMBILICAL CORD CLAMP (UMAFLOW)	31
58	M-7603	DISPOSABLE EXAMINATION PLASTIC GLOVES	31
59	M-7604	STERILE LATEX SURGICAL GLOVES (SOFTTEX)	32
60	M-7610	NON STERILE LATEX SURGICAL GLOVES (SOFTTEX)	32
61	M-7605	CORRUGATED DRAINAGE SHEET (UMAFLOW)	33
62	M-7606	DISPOSABLE FACE MASK	33
63	M-7607	DISPOSABLE CAP	33
64	M-7608	VETERINARY NEEDLE	34





The company has in-house quality assurance, quality control and testing laboratory to inspect at different stage in production cycle beginning with raw material to final product.



Saniludo has in-house highly equipped R & D section which helps to improve the quality standards of products constantly & to add new innovative product range as well.







## **INFUSION SET (UMAFLOW) T-7101**

(WITH "Y" CONNECTION & LUER LOCK)

- ✓ Confirms IS 12655
- ✓ Clear, transparent, flexible chamber with filter
- ✓ Soft & kink resistant medical grade PVC tube
- ✓ Big roller controller for accurate flow rate
- ✓ "Y" type injection port ( latex free )  
for extra medication
- ✓ Luer lock connection for extra safety
- ✓ Needle size : 21G
- ✓ Approximately 20 drops per ml  $\pm 0.1$  ml
- ✓ Individually packed & sterile



## **INFUSION SET (UMAFLOW) T-7102**

- ✓ Confirms IS 12655
- ✓ Clear, transparent, flexible chamber with filter
- ✓ Soft & kink resistant medical grade PVC tube
- ✓ Efficient roller controller for accurate flow rate
- ✓ Self sealing tube latex for extra medication
- ✓ Needle size : 21G
- ✓ Approximately 20 drops per ml  $\pm 0.1$  ml
- ✓ Individually packed & sterile



## **INFUSION SET (UMAFLOW) T-7103**

- ✓ Confirms IS 12655
- ✓ Clear, transparent, flexible chamber with filter
- ✓ Soft & kink resistant medical grade PVC tube
- ✓ Big roller controller for accurate flow rate
- ✓ Self sealing bulb latex for extra medication
- ✓ Needle size : 21G
- ✓ Approximately 20 drops per ml  $\pm 0.1$  ml
- ✓ Individually packed & sterile

## **INFUSION SET (UMAFLOW) T-7117**

- ✓ Confirms IS 12655
- ✓ Clear, transparent, flexible chamber with filter
- ✓ Soft & kink resistant medical grade PVC tube
- ✓ Efficient roller controller for accurate flow rate
- ✓ Bulb latex for extra medication
- ✓ Needle size : 21G
- ✓ Approximately 20 drops per ml  $\pm 0.1$  ml
- ✓ Individually packed in paper pouch & sterile







## **INFUSION SET (UMAFLOW) T-7104**

- ✓ Confirms IS 12655
- ✓ Clear, transparent, flexible chamber with air vented spike & air filter
- ✓ Soft & kink resistant medical grade PVC tube
- ✓ Big roller controller for accurate flow rate
- ✓ Bulb latex for extra medication
- ✓ Needle size : 21G
- ✓ Approximately 20 drops per ml  $\pm$  0.1 ml
- ✓ Individually packed & sterile



## **INFUSION SET UMAFLOW (MICRODRIP) T-7105**

- ✓ Confirms IS 12655
- ✓ Clear, transparent, flexible chamber with micro dropper & filter
- ✓ Soft & kink resistant medical grade PVC tube
- ✓ Roller controller for accurate flow rate
- ✓ Self sealing tube latex for extra medication
- ✓ Needle size : 23G
- ✓ Approximately 60 drops per ml  $\pm$  0.1 ml
- ✓ Individually packed & sterile





## **INFUSION SET UMAFLOW (MICRODRIP) T-7106**

- ✓ Confirms IS 12655
- ✓ Clear, transparent, flexible chamber with micro dropper & filter
- ✓ Soft & kink resistant medical grade PVC tube
- ✓ Big roller controller for accurate flow rate
- ✓ Bulb latex for extra medication
- ✓ Needle size : 23G
- ✓ Approximately 60 drops per ml  $\pm 0.1$  ml
- ✓ Individually packed in poly pouch & sterile



## **MEASURE VOLUME SET (UMADRIP) 150ML T-7107**

- ✓ Clear, transparent, flexible burette type chamber with micro dropper with capacity of 150 ml
- ✓ Chamber having retentive air inlet & injection site
- ✓ Automatic shutoff valve prevent air embolism
- ✓ Soft & kink resistant medical grade PVC tube
- ✓ Bulb latex for extra medication
- ✓ Needle size : 23 G
- ✓ Approximately 60 drops per ml  $\pm 0.1$  ml
- ✓ Individually packed in poly pouch & sterile







## **MEASURE VOLUME SET (UMADRIP) 110 ML T-7108**

- Clear, transparent, flexible burette type chamber with micro dropper with capacity of 110 ml
- Chamber having retentive air inlet & injection site
- Automatic shutoff valve prevent air embolism
- Soft & kink resistant medical grade PVC tube
- Bulb latex for extra medication
- Needle size : 23 G
- Approximately 60 drops per ml  $\pm$  0.1 ml
- Individually packed in poly pouch & sterile



## **BLOOD ADMINISTRATION SET (UMAFLOW) T-7109**

- Confirms IS 9824
- Clear, transparent, flexible drip chamber with filter to prevent particulate matter
- Soft & kink resistant medical grade PVC tube
- Efficient roller controller for accuracy
- Self sealing tube latex for extra medication
- Needle size : 18 G
- Individually packed & sterile





## **BLOOD ADMINISTRATION SET (UMAFLOW) T-7110**

- ✓ Confirms IS 9824
- ✓ Clear, transparent, flexible double drip chamber with filter to prevent particulate matter
- ✓ Soft & kink resistant medical grade PVC tube
- ✓ Efficient roller controller for accuracy
- ✓ Bulb latex for extra medication
- ✓ Needle size : 18G
- ✓ Individually packed in poly pouch & sterile

## **SCALP VEIN SET (UMAFLOW) T-7112**

- ✓ Designed to provide rapid venous access & greater patient comfort during the infusion
- ✓ Tube manufactured from the non-toxic & kink resistant PVC
- ✓ Thin walled siliconized steel cannula
- ✓ Provided with soft PVC wings for easy handling
- ✓ Each needle fitted with international colour code butterfly to easy detection of gauge of needle
- ✓ Individually packed in blister pack & sterile



SIZE	18G	20G	21G	22G	23G	24G
COLOUR	PINK	YELLOW	GREEN	BLACK	BLUE	RED





## HIGH PRESSURE MONITORING LINE (UMAFLOW) T-7114

- Small bore high pressure extension line suitable for high pressure monitoring
- With male & female connector
- Kink resistant medical grade PVC tube for pressure up to 54 bar
- Tubing diameter : I.D. 1.5 mm approx, O.D. 3.0 mm approx
- Individually packed in poly pouch & sterile

Tubing Length	50cm	100cm	150cm	200cm
---------------	------	-------	-------	-------



## EXTENSION LINE WITH 3-WAY STOP COCK (UMAFLOW) T-7116

- Extension line with three way stop cock at one end & male luer lock connector at other end
- Soft & kink resistant medical grade PVC tube
- Extension tube's diameter : I.D. 3.0 mm approx, O.D. 4.2 mm approx
- Individually packed in poly pouch & sterile

Tubing Length	10cm	25cm	50cm	100cm	150cm	200cm
---------------	------	------	------	-------	-------	-------






## PERITONEAL DIALYSIS SET (Y' SET) T-7111

- ❧ Specially designed "Y" shaped transfusion set for administration of peritoneal dialysis solution
- ❧ Provided with two upper clamps to facilitate the alternative change of bottle solution
- ❧ Lower clamps help easy import drainage of solution
- ❧ Individually packed in paper pouch & sterile

## URINE COLLECTING BAG (TINY) U-7208

(FOR PEDIATRIC USE)

- ❧ Manufactured from clinical grade plastic sheet
- ❧ Provided with adhesive property
- ❧ Capacity : 100ml approx
- ❧ Individually packed in paper pouch & sterile







## URINE COLLECTING BAG (UMAFLOW) U-7209

- Manufactured from clinical grade non-toxic, transparent PVC sheet
- Clear, transparent, soft, kink resistant medical grade PVC tube with adequate length provided
- Provided non-return valve to prevent backflow of fluid & bacterial migration
- Capacity : 2 liter (2000ml approx)
- Individually packed & sterile



## URINE COLLECTING BAG (UROFILL) U-7202

- Manufactured from clinical grade non-toxic, transparent PVC sheet
- Clear, transparent, soft & kink resistant medical grade PVC tube with adequate length provided
- Provided non-return valve to prevent backflow of fluid & bacterial migration
- Capacity : 2 liter (2000ml approx)
- Individually packed in poly pouch & sterile



## URINE COLLECTING BAG (UROSOFT) U-7203

- Manufactured from clinical grade non-toxic, transparent PVC sheet
- Clear, transparent, soft & kink resistant medical grade PVC tube with adequate length provided
- Provided non-return valve to prevent back flow off fluid & bacterial migration
- Provided with air filter at top of the bag
- Capacity : 2 liter (2000ml approx)
- Individually packed in poly pouch & sterile

## URINE COLLECTING BAG (ROYAL-19) U-7204

- Manufactured from clinical grade non-toxic, transparent, extra thick PVC sheet
- Clear, transparent, soft & kink resistant medical grade PVC tube with adequate length provided
- Provided non-return valve to prevent backflow of fluid & bacterial migration
- Provided with moulded hanger to carrying & hanging
- Capacity : 2 liter (2000ml approx)
- Individually packed in poly pouch & sterile







## URINE COLLECTING BAG (ROYAL-5) U-7212

(WITH BOTTOM OUTLET)

- Manufactured from clinical grade non-toxic, transparent extra thick PVC sheet
- Clear, transparent, soft & kink resistant medical grade PVC tube with adequate length provided
- Provided non-return valve to prevent backflow of fluid & bacterial migration
- Provided with moulded hanger to carrying & hanging
- Provided with bottom outlet for emptying of the bag
- Capacity : 2 liter (2000ml approx)
- Individually packed in poly pouch & sterile



## URINE BAG-LEG BAG SET (ROYALCARE) U-7211

- Used for short term & long term urine collection
- Provided freedom to the patient to move freely as the leg
- Available with leg straps
- Capacity : 800ml approx
- Individually packed & sterile



## URINE COLLECTING BAG WITH BOTTOM OUT LET (URO MEASURE) U-7210

(URINE COLLECTING BAG WITH MEASURE VOLUME METER FOR HOURLY MEASUREMENT OF URINE OUTPUT)

- Urine collecting bag marked with ml. for accurate measurement of urine output
- Close circuit system to eliminate the risk of contamination
- Soft, kink resistant, medical grade PVC tube with adequate length provided
- With bottom outlet for drainage
- Connecting tubing length 100 cm
- Capacity of volume meter : 250ml
- Capacity of urine bag : 2 liter (2000ml Approx)
- Individually packed in box & sterile

## NELATON CATHETER (UMAFLOW) U-7205

(DESIGNED FOR SHORT-TERM BLADDER CATHETERISATION THROUGH URETHRA)

- Manufactured from non-toxic, non-irritant medical grade PVC
- Frozen surface tubing for easy advancement
- Cone shaped distal end for atraumatic catheterization
- Two lateral eyes for efficient drainage
- Radio opaque line through out the catheter length
- Funnel shaped connector for easy connection to urine bag at proximal end
- Length : 40 cm
- Individually packed in paper pouch & sterile



SIZE	FG-8	FG-10	FG-12	FG-14	FG-16	FG-18	FG-20
COLOUR	BLUE	BLACK	WHITE	GREEN	ORANGE	RED	YELLOW





## **MALE EXTERNAL CATHETER (URICATH) U-7206**

(PENILE SHEATH)

- External catheter used for urine incontinence for day & night use for male patients
- Manufactured from pure latex provided with soft adhesive coated strip for proper fixing on the penis
- Proximal end is designed for safe & easy connection to urine container
- Individually packed in paper pouch & sterile



SIZE	SMALL	MEDIUM	LARGE	EX.LARGE
MM	20mm	25mm	30mm	35mm



## **T.U.R. SET (UMAFLOW) U-7207**

( "Y" SHAPE SET FOR ENDOSCOPIC IRRIGATION)

- Manufactured from non-toxic, non-irritant medical grade PVC
- "Y" shaped set for endoscopy irrigation during trans-urethral resection of prostate gland
- "G" clamp helps quick & smooth change of bottles
- Proximal end fitted with any kind of endoscope
- Individually packed in poly pouch & sterile



**RYLE'S TUBE  
(UMAFLOW) G-7301**

(FOR NASOGASTRIC INTRODUCTION  
FOR NUTRITION AND ASPIRATION OF  
INTESTINAL SECRETION)

- Ø Cone shaped distal end with radio opaque material for accurate placement
- Ø Four lateral eyes
- Ø Length marking on the tube
- Ø Soft, kink resistant medical grade PVC tube
- Ø Length : 105 cm
- Ø Individually packed in poly pouch & sterile

SIZE	FG-8	FG-10	FG-12	FG-14	FG-16	FG-18	FG-20
------	------	-------	-------	-------	-------	-------	-------

**RYLE'S TUBE  
(UMALINE) G-7304**

(NASOGASTRIC TUBE)

- Ø Cone shaped distal end with radio opaque material for accurate placement
- Ø Four lateral eyes
- Ø Length marking on the tube
- Ø Soft, kink resistant medical grade PVC tube
- Ø Colour code funnel end luer connector for easy identification of the sizes
- Ø Length : 105 cm
- Ø Individually packed in poly pouch & sterile



SIZE	FG-8	FG-10	FG-12	FG-14	FG-16	FG-18	FG-20
COLOUR	BLUE	BLACK	WHITE	GREEN	ORANGE	RED	YELLOW





## INFANT FEEDING TUBE (UMAFLOW) G-7302

(INFANT FEEDING TUBE FOR PEDIATRIC & NEONATAL USE)

- o Distal end is closed with two lateral eyes
- o Radio opaque line runs through the tube marked at 20 cm for accurate placement
- o Soft, kink resistant medical grade PVC tube
- o Available in female flexible connector
- o Length : 52 cm
- o Individually packed in paper pouch & sterile



SIZE	FG-5	FG-6	FG-7	FG-8	FG-9	FG-10
COLOUR	GRAY	GREEN	WHITE	BLUE	ORANGE	BLACK



## ABDOMINAL DRAINAGE KIT (UMA ADK) G-7303

UNIT CONSISTS OF

1. Extra soft drainage catheter
2. Collection bag of 2 liter (2000ml approx)
  - o Manufactured for postoperative abdominal drainage
  - o Soft & smooth catheter with large atraumatic eyes for efficient drainage
  - o Radio opaque line runs through out the catheter to help x-ray visualization
  - o 2 liter (2000ml approx) collecting bag to collect the abdominal fluid
  - o Specially designed hanger for carrying facilities
  - o Individually packed in poly pouch & sterile

SIZE	FG-16	FG-20	FG-24	FG-28	FG-30	FG-32	FG-36
------	-------	-------	-------	-------	-------	-------	-------



## COLOSTOMY KIT (UMAFLOW) G-7305

(COLON OSTOMY KIT)

- Colostomy kit used for long term for perfect relief & comfort
- Flexible light weight & given belt to fit the body well
- 15mm 'U' foam padding is provide between the sleeve & body for better comfort
- Provided with extra bags



## SUCTION CATHETER (UMAFLOW) A-7401

(TO REMOVE SECRETION FROM MOUTH, OROPHARYNX, TRACHEA & BRONCHIAL TUBE)

- Opened distal end with one lateral eye
- Soft, kink resistant medical grade PVC tube
- Proximal end with funnel type colour connector
- Length : 53 cm
- Individually packed in paper pouch & sterile



SIZE	FG-6	FG-8	FG-10	FG-12	FG-14	FG-16	FG-18	FG-20
COLOUR	GREEN	BLUE	BLACK	WHITE	GREEN	ORANGE	RED	YELLOW





## SUCTION CATHETER UMAFLOW (K-60, K- 61) A-7407

(TO REMOVE SECRETION FROM MOUTH,  
OROPHARYNX, TRACHEA & BRONCHIAL TUBE)

- Opened distal end with one lateral eye
- Soft, kink resistant medical grade PVC tube
- Proximal end with funnel type colour connector
- Length : 53 cm
- Individually packed in paper pouch & sterile



PRODUCT	K-60	K-61
SIZE	FG-14	FG-10



## SUCTION CATHETER THUMB CONTROL (UMAFLOW) A-7409

(REMOVE SECRETION FROM MOUTH,  
OROPHARYNX TRACHEA &  
BRONCHIAL TUBE)

- Opened distal end with one lateral eye
- "T" type vacuum control valve with thumb control is provided
- Frozen surface, kink resistant, medical grade PVC tube
- Length : 53 cm
- Individually packed in paper pouch & sterile

SIZE	FG-6	FG-8	FG-10	FG-12	FG-14	FG-16	FG-18	FG-20
COLOUR	GREEN	BLUE	BLACK	WHITE	GREEN	ORANGE	RED	YELLOW



## **URETHRAL CATHETER (UMAFLOW) A-7402**

(SHORT-TERM URINE DRAIN CATHETER)

- Use for short-term bladder catheterization through urethra
- Distal end is closed & proximal end has female colour code connector
- Soft, kink resistant, medical grade PVC tube
- Length : 40 cm
- Individually packed in paper pouch & sterile



SIZE	FG-8	FG-10	FG-12	FG-14	FG-16	FG-18	FG-20
COLOUR	BLUE	BLACK	WHITE	GREEN	ORANGE	RED	YELLOW

## **URETHRAL CATHETER (UMAFLOW(K-90, K-91) A-7408**

(SHORT-TERM URINE DRAIN CATHETER)

- Use for short-term bladder catheterization through urethra
- Distal end is closed & proximal end has female colour code connector
- Soft, kink resistant, medical grade PVC tube
- Length : 40 cm
- Individually packed in poly pouch & sterile



PRODUCT	K-90	K-91
SIZE	FG-14	FG-10





## **NASAL OXYGEN CATHETER (UMAFLOW) A-7403**

(CATHETER FOR OXYGEN ADMINISTRATION)

- Open distal end with multiple lateral eyes for even dispersion of oxygen & prevention of oxygen burns
- Straight pack
- Catheter length : 40 cm
- Individually packed in paper pouch & sterile



SIZE	FG-8	FG-10	FG-12	FG-14	FG-16	FG-18	FG-20
COLOUR	BLUE	BLACK	WHITE	GREEN	ORANGE	RED	YELLOW



## **TWIN BORE NASAL OXYGEN SET (UMAFLOW) A-7404**

- Suitable for easy & efficient administration of oxygen
- Star-lumen tube to avoid accidental blockage
- Soft, kink resistant, medical grade PVC tube
- Soft, twin prong nasal tip to ensure equal distribution of oxygen through both air passages
- Tube length : 2 meter (200 cm)
- Size : Adult, Child
- Individually packed in poly pouch



## **OXYGEN FACE MASK (UMAFLOW) A-7405**

(SUITABLE FOR OXYGEN THERAPY)

- ✓ Oxygen face mask for long-term use
- ✓ Clear & soft mask for patient comfort
- ✓ Adjustable elastic strip & integrated nose clip for proper positioning of mask
- ✓ Long star-lumen tube length : 2 meter (200 cm)
- ✓ Size : Adult, Child
- ✓ Individually packed in poly pouch



## **NEBULIZER KIT (UMAFLOW) A-7415**

KIT CONSISTS OF

1. Nebulizer mask
  2. Nebulizer chamber
  3. 200 cm star-lumen tube
- ✓ Soft, clear aerosol mask is ideal for long term therapy
  - ✓ Gentle rolled edges with nose clip for extra comfort
  - ✓ Rotating type connector for patient comfort
  - ✓ Specially designed nebulizer chamber with nebulizer rate of approximate 3cc
  - ✓ Individually packed in poly pouch





## **HIGH CONCENTRATION MASK (UMAFLOW) A-7413**

### UNIT CONSISTS OF

1. Mask with elastic strap & non-return valve
  2. Two meter star-lumen tube
  3. 600ml reservoir
- ☞ Specially designed for high concentration oxygen therapy
  - ☞ Clear & soft mask for patient comfort
  - ☞ Non-return valve provided at inhale & exhale point to prevent CO<sub>2</sub> Retention
  - ☞ Light weight reservoir improve oxygen efficiency & patient comfort
  - ☞ Individually packed in poly pack



## **ASEPTO PUMP (UMAFLOW) A-7414**

- ☞ Irrigation & feeding syringe from smooth tip facilitates easy & safe connection to any catheter for gastric tube
- ☞ Transparent bulb facilitate easy pumping
- ☞ Capacity : 60ml
- ☞ Individually packed in poly pouch & sterile





## **GUEDEL AIRWAY (UMAFLOW) A-7406**

(ORO-PHARYNGEAL AIRWAY)

- ✓ Suitable for maintain an unobstructed oro-pharyngeal airway during general anesthesia and in unconscious patients
- ✓ Bite block to prevent biting of tongue and air occlusion.
- ✓ Rounded atraumatic edges
- ✓ Individually packed & sterile

SIZE	00	0	1	2	3	4
COLOUR	BLUE	BLACK	WHITE	GREEN	ORANGE	RED

## **TRACHEOSTOMY TUBE UMAFLOW (PLAIN) A-7410**

(WITHOUT CUFFED)

- ✓ Made up of soft, smooth, medical grade PVC
- ✓ Lightweight & non-irritant material avoid inflammation around trachea
- ✓ Soft, flexible, flange help easy fixation
- ✓ Individually packed in blister pack & sterile



SIZE	4.0	4.5	5	6	7	7.5	8	9	10
------	-----	-----	---	---	---	-----	---	---	----





## ENDOTRACHEAL TUBE UMAFLOW (PLAIN) ORAL/ NASAL A-7411

(WITHOUT CUFFED)

- Nasoral tip suitable for nasal / oral intubation
- Manufactured from thermo sensitive material readily conforming body curvature stable at body temperature
- Radio opaque line throughout the length to facilitate exact location of the tube
- Proximal end fitted with standard 15mm connector
- Individually packed in blister pack & sterile



SIZE	2	2.5	3.0	3.5	4.0	4.5	5.0	5.5	6.0	6.5
------	---	-----	-----	-----	-----	-----	-----	-----	-----	-----



## RESPIRATORY EXERCISER (UMACISER) A-7412

UNIT CONSISTS OF

1. Transparent main body
  2. Small three balls
  3. Connecting tube with mouth piece
- Respiratory lung exerciser
  - Advance system to exercise respiration through inspiration
  - Designed with three chambers to provide smooth exercise to patient stepwise
  - Unit is composed of three-step chamber with different colour balls
  - By aspiration three ball will rise to maximum
  - Individually packed in box



## CHEST DRAINAGE CATHETER (UMAFLOW) S-7501

(INTRACOSTAL DRAINAGE CATHETER)

- Manufactured from soft, kink resistant medical grade PVC tube
- Atraumatic, soft, rounded tip with two lateral eyes for efficient drainage
- Individually packed in paper pouch & sterile

SIZE	FG-16	FG-20	FG-24	FG-28	FG-30	FG-32	FG-36
------	-------	-------	-------	-------	-------	-------	-------

## WATER SEALED CHEST DRAINAGE BAG (UMADRAIN) S-7504

- Use to collect drainage fluid from chest cavity
- PVC drainage bag with 1000ml capacity suitable for thoracic applications
- Soft, kink resistant, medical grade PVC tube with adequate length to ensure uninterrupted flow
- Special hanger given to carry the bag
- Individually packed & sterile








## **CLOSED WOUND SUCTION UNIT (UMA VAC SET) S-7502**

### UNIT CONSISTS OF

1. Bellow unit with connector
  2. Connecting tube with clamp & "Y" connector
  3. Curved trocar with catheter
  4. Spare, perforated drain catheter
- ✓ Suitable for closed wound drainage under negative pressure post operatively
  - ✓ With option to operate one or two catheter simultaneously
  - ✓ Calibrated bellow chamber
  - ✓ Bellow unit capacity : 800ml
  - ✓ Individually packed in box & sterile

SIZE	FG-10	FG-12	FG-14	FG-16	FG-18
COLOUR	BLACK	WHITE	GREEN	ORANGE	RED



## **CLOSED WOUND DRAINAGE SET (MINOR VAC SET) S-7503**

### UNIT CONSISTS OF

1. Bellow unit capacity 50 ml
  2. Curved needle with two catheters
- ✓ Catheter have radio opaque line & lateral eyes for efficient trauma free performance
  - ✓ Easy to operate by single person
  - ✓ Individually packed in box & sterile

SIZE	FG-6	FG-8
COLOUR	GREEN	BLUE



## VACCUM SUCK SUCTION SET (UMAFLOW) S-7505

(YANKAUR SUCTION SET)

- ✓ Suitable for pre-operative removal of secretion of body fluid
- ✓ Yankaur handle with four lateral eyes at distal end
- ✓ Yankaur handle with 2.5 meter long tubing & finger tip vaccum control facilities
- ✓ Proximal end can fitted with any type of suction apparatus
- ✓ Individually packed in poly pouch & sterile



## CONNECTING TUBE (UMAFLOW) S-7506

- ✓ Tubing is manufactured from non-toxic, non-irritant, medical grade PVC
- ✓ Tube is extra soft & flexible compare to latex tube
- ✓ Tubing length : 1.5 meter
- ✓ Individually packed in paper pouch & sterile







## **KARMAN CANNULA (UMAFLOW) S-7507**

(MTP CANNULA)

(WITH ADAPTOR & WITHOUT ADAPTOR)

- ♻ Specially designed for aseptic medical termination of pregnancy
- ♻ Coned shaped distal end with two large lateral eyes
- ♻ Proximal end fitted with MTP syringe or suction apparatus
- ♻ Individually packed in paper pouch & sterile

SIZE	5mm	6mm	7mm	8mm	10mm
------	-----	-----	-----	-----	------



## **UMBILICAL CANNULA (UMAFLOW) S-7509**

- ♻ Made up of soft, non-toxic PVC
- ♻ Designed for intermittent or continues access to umbilical artery's of newborn or pre mature infants
- ♻ Length : 40 cm
- ♻ Individually packed in paper pouch & sterile

SIZE	FG-6	FG-7	FG-8	FG-9
COLOUR	GREEN	WHITE	BLUE	ORANGE



## **C V P MANOMETER SCALE (UMAFLOW) S-7510**



- o C V P Manometer scale used to monitoring central venous pressure during infusion in I C U, I C C U, operation theater & trauma wards
- o Three way tap attach to C V P scale & scale is marked with 4 cm to 34 cm
- o 120 cm long tube with male & female connector for proper connection
- o Floating ball is provided to take reading of the pressure
- o Moulded clamps with latex rubber strings are given for fixation on I V Stand
- o Bulb latex given for extra medication
- o Individually packed in peel off pack & sterile

## **INFANT MUCUS EXTRACTOR (UMAFLOW) M- 7601**

- o Specially designed for aspiration of secretion from oropharynx in newly born babies
- o Chamber of mucus extractor is transparent & graduated
- o Suitable for mucus specimen for microbial examination
- o Individually packed in poly pouch & sterile







## **UMBILICAL CORD CLAMP (UMAFLOW) M-7602**

- ✓ Suitable for clamping the umbilical cord after the childbirth
- ✓ Security lock to prevent accidental opening after clamping
- ✓ Grooved clamping area to prevent slipping of the umbilical cord
- ✓ Individually packed in blister pack & sterile



## **DISPOSABLE EXAMINATION PLASTIC GLOVES M-7603**

1. FEATHER TOUCH (POLY PACK)
2. PROTECTION (ATTACH WITH PAPER)
1. FEATHER TOUCH
  - ✓ Manufactured from thin, soft HDPE
  - ✓ Individually packed & Sterile
2. PROTECTION
  - ✓ Manufactured from ultra thin, soft HDPE attach with paper
  - ✓ Available in 25'S bag



## STERILE LATEX SURGICAL GLOVES (SOFTTEX) M-7604

- ✓ Made from pure natural rubber latex for greater finger dexterity
- ✓ Non-finger fatigue keeps the hand cool & comfortable
- ✓ Powder free surgical gloves
- ✓ For single use only
- ✓ Individually packed & sterile

SIZE	6	6½	7	7½	8
------	---	----	---	----	---

## NON STERILE LATEX SURGICAL GLOVES (SOFTTEX) M-7610

- ✓ Made from natural rubber latex for greater finger dexterity
- ✓ Non-finger fatigue keeps the hand cool & comfortable
- ✓ Powder free surgical gloves
- ✓ Available in 50'S box



SIZE	6	6½	7	7½	8
------	---	----	---	----	---






## **CORRUGATED DRAINAGE SHEET (UMAFLOW) M-7605**

(MULTICHANNEL DRAINAGE SYSTEM)

- ✧ Manufactured from non-toxic, non-irritant medical grade PVC
- ✧ Multichannel drainage system for efficient drainage during operation
- ✧ Provided with x-ray opaque line
- ✧ Size : 25mm × 250mm
- ✧ Individually packed & sterile



## **DISPOSABLE FACE MASK M-7606**

- ✧ Made from soft, non-woven cloth
- ✧ Double layer mask offering excellent bacterial filtration efficiency
- ✧ Individually packed & sterile




## **DISPOSABLE CAP M-7607**

- ✧ Specially cap designed to protect the patient against dropping of hair/dandruff of the doctors and other attending staff
- ✧ Provided with stitched interlocking adaptable elastic for better grip on the forehead & comfortable fitting
- ✧ Individually packed & sterile





## VETERINARY NEEDLE (UMAFLOW) M-7608

-  Veterinary needle used for veterinary use
-  Needles are available with plastic hub
-  Available in 10'S box



GAUGE	18	18	18	18	20	20	20
LENGTH	1/2"	3/4"	1"	1 1/2"	3/4"	1"	1 1/2"

