





Crafted to Excellence



SURTEX INSTRUMENTS

Crafted to Excellence

INTRODUCTION

• INTRODUCTION TO THE CATALOGUE

In this catalogue we are pleased to introduce the ophthalmic surgery instruments range available from Surtex Instruments. Following feedback from customers we have prepared this as a standalone catalogue to allow for quick referencing and searching of our products required on a routine basis. The specialist ranges of our instruments are covered in the following additional catalogues available upon request.

- ENT
- Cardiovascular & Neurosurgery
- Ophthalmic
- Gynaecology and Obstetrics
- Dental Instruments
- Sterilization & Holloware

• ABOUT SURTEX INSTRUMENTS

Surtex Instruments is a privately owned medical technology company dedicated to innovation and excellence in crafting of surgical instruments. The company started as a specialist in reconstructive surgery instruments and based on its success and innovative drive gradually expanded to cover a vast range of surgical instruments.

We now manufacture over 8,000 surgical instruments and offer innovative solutions in trauma, gynecology, cardiovascular, microsurgery and dental surgery. All our products are available to view online where you can download our entire range of catalogues or request a printed version.



• SURTEX QUALITY STANDARDS

In European Union medical devices are subject to a series of uniform laws and the EC directive places rigid demands on surgical instruments. Surtex Instruments is proud to claim that we have always surpassed the legal requirements to ensure that each instrument reaching the surgeons hand is indeed a testament to surgical craftsmanship.

We are **BS EN ISO 13485:2012/ISO 9001:2008 (SGS UK)** certified and **CE** approved to ensure compliance with the standards defined in line with medical devices directive **93/42/EEC**.



WHY IS SURTEX FOCUSED ON INSTRUMENTS QUALITY?

Surgeons are highly talented and skilled individuals. They operate under significant time constraints and are required to make split second decisions while performing complex operations. At Surtex we understand that in such a scenario the surgeon requires tools that are not just fit for the purpose but are reliable, dependable and crafted by a company that has a track record for relentless pursuit for excellence in the making of surgical instruments.



• AFFORDABLE QUALITY

Surgical instruments represent a significant proportion of hospital's total assets. Hospitals and health centers are often required to minimize costs of surgical instruments without impacting quality or patient's health. Decision makers tend to respond via group purchasing contracts & multiple quotes to shop around for the lowest price. However this approach has significant drawbacks as numerous instrument companies have focused towards competing on price alone with a disregard towards instruments quality while claiming their instruments to be guaranteed and free of defect.

The outcome for the buyer being a perceived notion of a quality instrument, which is actually little more than a knock-off and delivers only poor value for money.

At **Surtex Instruments** our success has been based on helping customers navigate successfully through the instrument procurement process by identifying real value for money through fair comparison. We believe that quality should be affordable and you should not pay too much for over priced brands and neither should you pay too much for a low quality instrument.

• LIFE TIME WARRANTY

We encourage our customers to calculate lifetime costs prior to making purchase decisions and investing in surgical instruments.

Surtex Instruments are guaranteed for life* against manufacturing defects of material and workmanship. This means we will repair or replace free of charge any surgical instrument that does not meet its functional requirements when used for its intended purpose and maintained properly.



A member of the ABHI
ABHI
Health Tech For Life
Association of
British HealthTech
Industries

*Conditions apply.

[OUR PRODUCT BRANDS]



A new dimension in cutting and dissecting performance. Our class leading integral self-sharpening technology standard on all Surtex scissors enabling precise cutting with minimum force even after pro-longed use.



These Surtex premium scissors feature finely ground micro serrations along one blade of the scissor preventing tissue slippage. The ultra sharp edge of the other blade cuts through the tissue precisely and atraumatically.



Instruments with this symbol benefit from tungsten carbide (TC) inserts expertly bonded on the cutting edge of the instrument. This material is harder than stainless steel and enables the instruments to last up to 6 times longer than standard instruments.



A true craftsmanship masterpiece and surgeons favourite. These are the best scissors in our range and combine the features of SuperEdge™ & UltraCut™ into a perfect pair of scissors that is precise, durable and reliable.



The only choice when the perfect grip matters. Utilising Creep Feed technology individual teeth are grinded into hard wear resistant tungsten carbide (TC). The result is super sharp teeth with well-defined point and inserts with impeccable gripping properties and instruments that last up to 6 times longer.



These instruments have jaws, which are coated with diamond dust resulting in improved gripping enabling superior manipulation of sutures, needles and tissues.



Our specially designed high strength fibre grip handles resulting in lightweight instruments. These ergonomic handles provide the surgeon with excellent feel, better grip, reduced shock, and less hand fatigue.



Surtex ENT diagnostic and anaesthesia instruments featuring direct bulb illumination of xenon grade bulbs.



This range features distal fibre optic to carry light from xenon bulb allowing cold illumination on operating site in our ENT diagnostics and anaesthesia instruments.



Premium grade ENT diagnostics and anaesthesia instruments combining features of distal fibre optic technology and ultra bright LED illumination ensuring the best in clear vision of operating site.

**CURRENT REFERENCE STANDARDS
FOR SURTEX PRODUCTS**

STANDARDS	CONTENTS OF THE STANDARD
BS EN ISO 9001:2008	Quality Management Systems certification standard.
BS EN ISO 13485:2003	Medical Devices – Quality Management Systems certification standard.
MDD 93/42/EEC, Annex II, Section 3.2	Medical Device Directive, 'EC Declaration of Conformity' (full quality assurance system).
MDD 93/42/EEC, Annex I	Essential Requirements for Medical Devices.
BS EN ISO 7153-1:2001 (BS 5194:Part 1:1991)	Surgical Instruments – Metallic materials part 1 - Stainless Steel
BS 5194:Part 2:1989	Specification for instruments with pivot joints (excluding cutting instruments)
BS 5194:Part 3:1985	Specification for Dissecting Forceps
BS 5194:Part 4:1985	Specification for scissors, shears, and other jointed cutting instruments
ISO 7151:1988	Surgical instruments – non-cutting, articulated instruments – general requirements & test methods
ISO 7740:2004	Scalpels with detachable blades - Fitting dimensions
ISO 7741:1986	Scissors and shears - General requirements and test methods
ISO 6508-1:2005 (DIN 50103-3:1995)	Specifies the method for Rockwell and Rockwell superficial hardness tests for metallic materials.
BS EN ISO 14971:2012	Medical Devices – Application of risk management to medical devices.
DIN 58298:2010-09	Medical instruments - Materials, finish and testing.
DIN 58299:1964-01	Serrations for surgical instruments; profile angles, groove distances.
DIN 58300:1982-02	Joints for surgical instruments.
HTM 2010	Recommendations for sterilization
AS/NZ 4187-2003	Cleaning, disinfecting & sterilising of re-usable surgical instruments

This information is available on the internet at
www.surtx-instruments.com

INSTRUCTIONS FOR REPROCESSING OF REUSABLE SURGICAL INSTRUMENTS

The following instructions are for all reusable medical devices supplied by Surtex Instruments Ltd, unless stated otherwise with the packaging of the product. These instructions are intended for use only by persons with the required specialist knowledge and training.

Warnings

- Follow instructions and warnings as issued by the manufacturers of any decontaminants, disinfectants and cleaning agents. Wherever possible avoid use of mineral acids and harsh, abrasive agents.
- No part of the process shall exceed 140 °C.
- Some sensitive materials (e.g. Aluminium) are damaged by high alkaline solutions (pH>10).
- Devices with long, narrow cannula, hinges and blind holes require particular attention during cleaning.

Note: When reprocessing medical devices, always handle with care, wearing protective clothing, gloves and eye-wear in accordance with local health & safety procedures.

Limitations on Reprocessing

- Repeated processing has minimal effect on these instruments.
- End of life normally determined by wear and damage in use.
- Any specific limitations on the number of reprocessing cycles shall be made available with instrument.

From Point of Use

Whenever possible, do not allow blood, debris or bodily fluids to dry on instruments. For best results and to prolong the life of the medical device reprocess immediately after use. If they cannot be reprocessed immediately, use an enzymatic foam spray cleaner to help prevent soil from drying.

Preparation for Decontamination

- Process all the instruments as soon as it is reasonably practical following use.
- Disassemble only where intended, without the use of tools unless specifically provided by the manufacturer. Instructions for disassembly are available with the device, where required.

Cleaning: Automated

- Use only either CE marked or validated washer-disinfector machines and low foaming, non-ionizing cleaning agents and detergents following the manufacturer's instructions for use, warnings, concentrations and recommended cycles.
- Load instruments carefully, with any box joints and hinges open and so that any fenestration in instruments can drain.
- Place heavy instruments with care in the bottom of containers, taking care not to overload wash baskets.
- Place instruments with concave surfaces facing down to prevent pooling of water.
- Where possible, use appropriate attachments to flush inside reamers and devices with lumens or cannula.
- Ensure that soft, high purity water which is controlled for bacterial endo-toxins is used in the final rinse stage.

Note: Automated cleaning may not be suitable for all lumens and cannula in which case clean manually with a water jet gun, if available, and an appropriate brush (and stilette if provided) that reaches the depth of the feature. After manually cleaning, pass all devices through an automatic cleaning cycle to achieve disinfection.

Note: These instruments have been validated using a washer-disinfector cycle validated to include two cold rinses at <35 °C, a detergent cycle and a rinse cycle both at >50°C, a disinfection cycle operating at a temperature of between 80 °C and 90 °C for a minimum holding time of 1 minute (actual holding time in excess of 2 minutes 50 seconds) and a 20 minute drying cycle. The detergent used was low foaming, non-ionizing spray wash detergent cleaner (max 12pH) and the rinse aid a neutral pH low foaming, non-ionic surfactant with isopropyl alcohol.

Cleaning: Manual

- Manual cleaning is not advised if an automatic washer-disinfector is available. If the equipment is not available, use the following process:-
 - Use a double sink system (wash/rinse) dedicated for instrument cleaning (not used for hand washing). Ensure that the water temperature does not exceed 35 of the re-processor to ensure that the reprocessing as actually performed using equipment, materials and personnel in the reprocessing facility achieve the desired result. This requires validation and routine monitoring of the process. Likewise any deviation by the re-processor from the instructions provided must be properly evaluated for effectiveness and potential adverse consequences.
 - In the first sink, keeping the instrument submerged, with an autoclavable brush, apply CE marked cleaning solution to all surfaces until all soil has been removed. Pay particular attention to serrations, teeth, ratchets and hinges, always brushing away from the body and avoiding splashing. Ensure rongeurs and hinged instruments are thoroughly cleaned in both open and closed positions.
 - In the second sink, rinse instruments thoroughly with soft, high purity water which is controlled for bacterial endotoxins, so that the water reaches all parts of the instrument, then carefully hand dry or use a drying cabinet.

Note: Manual Cleaning is NOT a disinfection process: When manual cleaning is used it may not be possible to disinfect the device prior to further handling.

Cleaning Inspection

- After cleaning visually inspect, visually inspect all surfaces, cannulations, ratchets, joints, holes and lumens for complete removal of soil and fluids. If ANY soil or fluid is still visible, return the instrument for repeat decontamination.

Maintenance

- Apply surgical grade lubricant to hinges, joints and moving parts as per the lubricant manufacturer's instructions.

Continued on next page

Inspection and Functional Testing

- Visually inspect and check: - all instruments for damage and wear; cutting edges are free from nicks and present a continuous edge; jaws and teeth align correctly; all articulated instruments have a smooth movement without excess play; locking mechanisms (such as ratchets) fasten securely and close easily; long, slender instruments are not distorted; any component parts fit and assemble correctly with mating components.
- Remove for repair and replacement any blunt, worn out, flaking, fractured or damaged instruments.

Note: If an instrument is returned to the manufacturer / supplier, the instrument must be decontaminated and sterilized and be accompanied with the relevant documented evidence.

Packaging

- All instruments to be packed following local protocol in accordance with ISO 11607-1 or AAMI/CSR technique.

Sterilization (USA)

- Autoclave should comply with the requirements of, and be validated and maintained in accordance with ANSI/AAMI ST79.
- Pre-Vacuum moist heat sterilization operating at temperature 132 °C bar for a minimum holding time of 4 minutes always following the instructions of the machine manufacturer.
- When sterilizing multiple instruments in one autoclave cycle, ensure that the sterilizer manufacturers stated maximum load is not exceeded.
- Ensure instruments are dry before sterilization.

Sterilization (Outside USA)

- Autoclave should comply with the requirements of, and be validated and maintained in accordance with EN 285, EN 13060, EN ISO 17665.
- Either CE marked or validated vacuum autoclave operating at 134-137 °C bar for a minimum holding time of 3 minutes always following the instructions of the machine manufacturer.
- When sterilizing multiple instruments in one autoclave cycle, ensure that the sterilizer manufacturers stated maximum load is not exceeded.

Storage

- Ensure instruments are dry before storage, and stored in dry, clean conditions at an ambient room temperature.

Additional Information

- Other forms of cleaning (i.e. ultrasonic) and sterilization (i.e. Low temperature steam and Formaldehyde, Ethylene Oxide and Gas Plasma) are available. However, always follow the instructions for use as issued by the manufacturer and always consult with them if in any doubt over the suitability of any process used.
- Cleaning and sterilizing guidelines are available in HTM 2030 and HTM 2010. Contact: The NHS Estates Stationers Office Publications Cent for details at www.tsionline.gov.uk. For further information contact: NHS Estates Information Center, Department of Health, 1 Trevelyan Square, Boar Lane, Leeds, LS1 6AE, UK or visit www.nhsestates.gov.uk.


Note: It is the responsibility of the re-processor to ensure that the reprocessing as actually performed using equipment, materials and personnel in the reprocessing facility achieve the desired result. This requires validation and routine monitoring of the process. Likewise any deviation by the re-processor from the instructions provided must be properly evaluated for effectiveness and potential adverse consequences.

TABLE OF CONTENTS (A - Z)

Cannulas	28 - 38
Chisels & Gouges	166
Curettes & Spoons	39 - 40, 72
Fascia Needles, Cautery Balls	9
Fixation Rings	82 - 83
Forceps	123 - 162
Foreign Body Instruments	16 - 17
Haemostats & Bulldog Clamps	163
Hooks	81, 83 - 86, 178
Knives	11 - 24
Lens Loops	64 - 66
Markers	86 - 89
Mallet	78
Micro Manipulators	40 - 62, 66 - 70, 73 - 80
Needle Holders	167 - 177
Probes & Dilators	90 - 92, 179
Punches	26
Retractors, Lid Plates	7 - 10
Rongeurs	165, 179
Rulers & Calipers	90, 178
Scalpel Handles & Blades	93 - 94
Scissors	95 - 122
Sclarel Markers & Depressors	7
Speculums	1 - 6
Sphere Introducers & Enucleation Instruments	71
Storing Cases	180
Towel Clamps	164
Trephines	25, 27
Ophthalmic Surgeries Sets	191-203



Eye Speculum
Lidsperrer
Blefarostatos
Blépharostats
Blefarostati



Barraquer Wire Speculum

OP-010-10

Pediatric

Wire Ø 0.9 mm


Blade Size 10 mm

OP-010-15

Adult

Wire Ø 0.9 mm

Blade Size 15 mm



OP-012-15

Barraquer Wire Speculum

Heavy Wire Ø 1.0 mm

Blade Size 15 mm



Barraquer Wire Speculum
(Solid Blades)

OP-020-10

Pediatric

Blade Size 10 mm

OP-020-15

Adult

Blade Size 15 mm



OP-022-20

Feaster Wire Speculum

Solid Blades

Blade Size 20 mm



Kratz-Barraquer Wire Speculum

OP-030-13

With Open Blades, Wire Ø 0.9 mm

Blade Size 13 mm

OP-030-15

With Open Blades, Wire Ø 0.9 mm

Blade Size 15 mm



OP-032-15

Kratz-Barraquer Wire Speculum

With open blades, Wire Ø 1.0 mm

Blade Size 15 mm

Eye Speculum
 Lidsperrer
 Blefarostatos
 Blépharostats
 Blefarostati

OP

OPHTHALMIC SURGERY INSTRUMENTS



OP-040-10
Frankel Infant Lid Speculum
 Adjustable Opening
 Blade Size 10 mm



OP-042-00
Alfonso Infant Lid Speculum
 With Finger Grip, Wire Ø 0.9 mm
 Blade Size 5 mm



OP-046-00
Infant Lid Speculum
 Wire Ø 0.9 mm
 Blade Size 4 mm



OP-050-09
Barraquer Lid Speculum (Infant)
 Wire Ø 0.9 mm
 Blade Size 9 mm



- Cook Eye Speculum**
- | | |
|--|-----------|
| OP-060-08 | Infant |
| Overall Length 52 mm, Blade Size 8 mm | |
| OP-060-11 | Pediatric |
| Overall Length 55 mm, Blade Size 11 mm | |
| OP-060-15 | Adult |
| Overall Length 55 mm, Blade Size 15 mm | |



Eye Speculum
Lidsperrer
Blefarostatos
Blépharostats
Blefarostati




Murdoch Eye Speculum

OP-064-12

Overall Length 55 mm
Blade Size 12 mm

OP-064-15

Overall Length 57 mm
Blade Size 15 mm



OP-070-15

Mellinger Eye Speculum

Overall Length 57 mm
Blade Size 15 mm



OP-072-02

Bangerter Eye Speculum

Overall Length 55 mm With Locking Screw
Blade Size 15 mm



Williams Eye Speculum

OP-074-11

Overall Length 69 mm, With Locking Screw
Blade Size 11 mm

OP-074-15

Overall Length 86 mm, With Locking Screw
Blade Size 15 mm

Eye Speculum
 Lidsperrer
 Blefarostatos
 Blépharostats
 Blefarostati


OP
OPHTHALMIC SURGERY INSTRUMENTS

OP-076-16
Lancaster Eye Speculum

Overall Length 69 mm, Adjustable Opening
 Blade Size 16 mm


OP-078-00
Lester-Burch Eye Speculum

Overall Length 82 mm, Adjustable Opening
 Blade Size 21 mm


OP-081-15
Tanna-Lieberman Eye Speculum

Overall Length 86 mm, Adjustable Opening
 Blade Size 15 mm


Lieberman Eye Speculum
OP-083-10

Pediatric
 Overall Length 73 mm, Adjustable Opening
 Blade Size 10 mm


OP-083-15 A
Lieberman Eye Speculum

With Aspirating Ports Overall Length 86 mm
 Adjustable Opening Blade Size 15 mm


OP-085-15
Lieberman Eye Speculum


Overall Length 86 mm Adjustable Opening
 Blade Size 15 mm

OP-083-15
Adult

Overall Length 86 mm, Adjustable Opening
 Blade Size 15 mm



Eye Speculum
Lidsperrer
Blefarostatos
Blépharostats
Blefarostati



Lieberman Eye Speculum
(Solid Blades)

OP-087-10

Pediatric
Overall Length 73 mm, Adjustable Opening
Blade Size 10 mm

OP-087-15

Adult
Overall Length 86 mm, Adjustable Opening
Blade Size 15 mm



OP-089-15 A
Double X-K Wire Eye Speculum
With Aspirating Ports Overall Length
86 mm, Adjustable Opening
Blade Size 15 mm




Castroviejo Eye Speculum

OP-091-14

Overall Length 75 mm
Adjustable Opening, 14 mm

OP-091-16

Overall Length 82 mm
Adjustable Opening, 16 mm



OP-089-15

Double X-K Wire Eye Speculum
Overall Length 86 mm, Adjustable Opening
Blade Size 15 mm




OP-095-04
Saunders Eye Speculum

For Premature Infants Overall Length 35 mm
With Locking Thumb Screw Blade Size 4 mm

Scleral Fixation Ring
 Scleral Fixation Ring
 Anillo de fijación escleral
 Anneau de fixation scléral
 Sclerale Fixation Anello

OP

OPHTHALMIC SURGERY INSTRUMENTS


OP-100-00
Flieringa Scleral Fixation Ring

Set of 9 From OP-100-14 up to OP-100-22


OP-100-14
 14 mm

OP-100-15
 15 mm

OP-100-16
 16 mm

OP-100-17
 17 mm

OP-100-18
 18 mm

OP-100-19
 19 mm

OP-100-20
 20 mm

OP-100-21
 21 mm

OP-100-22
 22 mm

McNeill-Goldman Scleral Fixation Ring
OP-111-15

Small, Ø 15 mm

OP-111-16

Medium Ø 16 mm

OP-111-17

Medium Ø 17 mm

Depressors, Eyelid Retractors
 Depressors, Augenlid Wundhaken
 Depresores, Retractor de párpados
 Dépresseurs, Écarteur de paupières
 Depressore, Palpebra Divaricatore



OP-151-01
Gass Scleral Marker
 1.5 mm Marking Tip



Schepens Scleral Depressor

OP-153-01

Tip Width 7 mm Thimble Ø 15 mm Tapered to Ø 13 mm

OP-153-02

Tip Width 7 mm Thimble Ø 16 mm Tapered to Ø 14 mm

OP-153-03

Tip Width 7 mm Thimble Ø 17 mm Tapered to Ø 15 mm



OP-160-14
Shocket Scleral Depressor
 Double End With Clip
 Overall Length 14 cm (5 1/2")



OP-163-00
Gonnin Amsler Marker
 Measuring Range 7 - 23 mm With 1 mm
 Increments
 Overall Length 13 cm (5 1/4")



OP-171-15
Schepens Orbital Retractor
 Blade Size 14 x 56 mm With 4.5 mm Wide Notch
 Overall Length 15 cm (6")

Eyelid Retractors
 Augenlid Wundhaken
 Retractor de párpados
 Écarteur de paupières
 Palpebra Divaricatore

OP

OPHTHALMIC SURGERY INSTRUMENTS



Sauter Orbital Spatula
 Overall Length 14 cm (5 ½")

OP-172-01

Blade Size 22 x 13 mm

OP-172-02

Blade Size 22 x 16 mm

OP-172-03

Blade Size 22 x 19 mm



Jaffe Wire Lid Retractor (Pair)

OP-173-13

Blade Width 13 mm

OP-173-15

Blade Width 15 mm



Desmarres Lid Retractor
 Overall Length 13 cm (5")

OP-175-11
 Thin Solid Blades
 Fig. 0, 11 mm

OP-175-15
 Thin Solid Blades
 Fig. 2, 15 mm



OP-175-13
 Thin Solid Blades
 Fig. 1, 13 mm



OP-175-17
 Thin Solid Blades
 Fig. 3, 17 mm



Scleral Shield

OP-180-01

20 x 28 mm, Small (Pair)

OP-180-02

22 x 30 mm, Large (Pair)

Eyelid Retractors
 Augenlid Wundhaken
 Retractor de párpados
 Écarteur de paupières
 Palpebra Divaricatore



OP-192-11
Jaeger Lid Plate
 20 x 35 / 24 x 35
 Overall Length 11 cm (4 1/4")



OP-197-14
Wright Fascia Needle
 Oval Hole, 1 x 6 mm
 Overall Length 14 cm (5 1/2")



OP-199-14
Wills Cautery
 With Copper Ball Ø 12 mm
 Overall Length 14 cm (5 1/2")



OP-201-12
Walton Conjunctiva Retractor
 Thin Solid Blade, 6 mm
 Overall Length 12.5 cm (5")



Helveston Tissue Retractor
 Thin Curved Blade
 Overall Length 12.5 cm (5")



OP-211-07
"Barbie"
 Blade Width 7 mm



OP-211-09
"Big Barbie"
 Blade Width 9 mm



OP-211-11
"Great Big Barbie"
 Blade Width 11 mm

Eyelid Retractors
Augenlid Wundhaken
Retractor de párpados
Écarteur de paupières
Palpebra Divaricatore

OP

OPHTHALMIC SURGERY INSTRUMENTS

**OP-220-00****Stevenson Lacrimal Sac Retractor**

3 x 3 Blunt Prongs, Adjustable Opening
Overall Length 85 mm

**OP-222-00****Stevenson Lacrimal Sac Retractor**

Solid Blades With Serrated Edge Adjustable
Opening Blade Size 13 mm
Overall Length 85 mm

**OP-226-00****Agricola Lacrimal Sac Retractor**

3 x 3 Blunt Prongs
Overall Length 14 mm



Keratome
Keratome
Queratomos
Kératotomes
Cheratomo



Dastoor Gouge & Lacrimal Sac Dissector
Overall Length 13 cm (5 1/4")

OP-228-01

With Sharp Gouge and Sharp Dissector

OP-228-02

With Sharp Gouge and Blunt Dissector



Bonn Scleral Knife
Rounded Tip
Overall Length 12 cm (4 3/4")

OP-239-40

Angled, 4 mm Wide Tip

OP-239-70

Angled, 7 mm Wide Tip



Jaeger Keratome
Overall Length 13 cm (5")

OP-240-01
Straight, Fig. 1
8.0 mm x 9.5 mm

OP-240-02
Straight, Fig. 2
9.5 mm x 11.0 mm

OP-240-03
Straight, Fig. 3
11.0 mm x 12.5 mm

OP-241-01
Angled, Fig. 1
8.0 mm x 9.5 mm

OP-241-02
Angled, Fig. 2
9.5 mm x 11.0 mm

OP-241-03
Angled, Fig. 3
11.0 mm x 12.5 mm

Corneal Knives
 Hornhaut-Messer
 Cuchillos corneales
 Couteaux cornée
 Coltelli corneali

OP

OPHTHALMIC SURGERY INSTRUMENTS


OP-250-11
Tooke Corneal Knife

Straight, Blade Size 3 x 18 mm Smooth Blade
 Overall Length 11 cm (4 1/4")


OP-251-01

Fig. 1, Blade Size 5.5 mm x 2.5 mm

OP-251-02

Fig. 2, Blade Size 4.5 mm x 2.0 mm


OP-253-12
DLEK Corneal Dissector

12 mm Long x 3 mm Wide With Semi-Sharp Spatulated Tip
 Overall Length 12 cm (4 3/4")


OP-255-12
DLEK Corneal Dissector

12 mm Long x 3 mm Wide With Semi-Sharp Spatulated Tip
 Overall Length 12 cm (4 3/4")



Corneal Knives
Hornhaut-Messer
Cuchillos corneales
Couteaux cornée
Coltelli corneali



OP-256-01

Fig. 1, 1.5 x 25 mm

OP-256-02

Fig. 2, 1.7 x 27 mm

OP-256-03

Fig. 3, 2 x 30 mm

OP-256-04

Fig. 4, 2.3 x 33 mm

OP-256-05

Fig. 1, 2.5 x 35 mm

Graefe Cataract Knife

Straight Blade and With Sharp Cutting Edge
Overall Length 12 cm (4 $\frac{3}{4}$ ")



OP-257-05

Fig. 1, 5 mm

OP-257-06

Fig. 2, 6 mm

OP-257-07

Fig. 3, 7 mm

Ziegler Iris Knife

Overall Length 12 cm (4 $\frac{3}{4}$ ")

OP-258-12

Straight, Blade Size 2.4 mm
Wide x 4.5 mm Long Blade



OP-259-12

Curved, Blade Size 2.4 mm
Wide x 4.5 mm Long Blade



Walter Cataract Needle

Overall Length 12 cm (4 $\frac{3}{4}$ ")

Corneal Knives
Hornhaut-Messer
Cuchillos corneales
Couteaux cornée
Coltelli corneali

OP

OPHTHALMIC SURGERY INSTRUMENTS

**OP-260-12****Gills Corneal Knife**

Straight, Rounded Cutting Edge Blade Size 3 mm
Overall Length 12.5 cm (5")

**OP-261-12****Gills Corneal Knife**

Slightly Curved, Rounded Cutting Edge Blade Size 3 mm
Overall Length 12.5 cm (5")

**OP-263-12****Pauifique Graft Knife**

Blade Size 1.6 mm x 1.75 mm Angled Cutting Edge
Overall Length 12 cm (4 3/4")

**OP-265-12****Fukasaku Hockey Knife**

3 mm Sharp Front Edge 6 mm Semi Sharp Bottom Edge
Overall Length 12 cm (4 3/4")

Corneal Knives
 Hornhaut-Messer
 Cuchillos corneales
 Couteaux cornée
 Coltelli corneali



OP-266-11

Barkan Goniotomy Knife

Blade Size 2.5 mm Long Blade
 Overall Length 11.5 cm (4 1/2")



OP-268-01

Micro Round Viscocanalostomy Knife


1.0 mm Rounded Sharp Cutting Edge
 Overall Length 11.5 cm (4 1/2")

Gouges and Spuds
 Hohlmeißel und meißel
 Gubias para Cinceles
 Gouges et Ciseaux-burins
 Sgorbie e Scalpelli


OP
OPHTHALMIC SURGERY INSTRUMENTS

OP-271-12
Golf Club Foreign Body Spud


Blade Size 2 x 0.8 mm Small Golf Club Shaped Tip
 With Sharp Edge
 Overall Length 12 cm (4 ¾")


OP-273-12
Francis Foreign Body Spud

Blade Size 2 x 0.8 mm Flag Shaped Tip With Sharp Edge
 Overall Length 12 cm (4 ¾")


OP-275-11
Ellis Foreign Body Spud

Light Curved 1.0 mm Rounded Tip
 Overall Length 11 cm (4 ¼")




Front View

Side View


OP-276-12
Walton Foreign Body Gouge

1.4 mm Rounded Tip
 Overall Length 12.5 cm (5")

Gouges and Spuds
 Hohlmeißel und meißel
 Gubias para Cinceles
 Gouges et Ciseaux-burins
 Sgorbie e Scalpelli



OP-280-18

Dix Foreign Spud & Gouge

3 X 5 mm Needle Tip On One Side 1.7 X 15 mm
 Gouge Tip On Other Side In Reversible Handle
 Overall Length 18 cm (7")



OP-285-00

Arruga Micro Burr

1.0 mm Round Burr
 Overall Length 12 cm (4 3/4")



OP-283-11

Wills Eye Hospital Foreign Body Spud & Burr

Delicate Curved
 11 cm (4 1/4")



OP-291-11

Foreign Body Eye Magnet and Loop

Overall Length 11 cm (4 1/4")

Diamond Knives
 Diamantmesser
 Bisturí de diamante
 Couteaux de diamant
 Coltelli di diamanti

OP
OPHTHALMIC SURGERY INSTRUMENTS

OP-300-12

 6-Facet Blade for Bi-directional
 Cutting Pre-calibrated Depth
 Settings of 0.6 mm

OP-302-12

 6-Facet Blade for Bi-directional
 Cutting Pre-calibrated Depth
 Settings of 0.5 mm, 0.55 mm
 and 0.6 mm

OP-304-12

 6-Facet Blade for Bi-directional
 Cutting 0 to 1.5 mm Blade Extension
 in 10 Micron Increments

LRI Diamond Step Knife
 Overall Length 12 cm (4 $\frac{3}{4}$ ")

LRI Diamond Micrometer Knife
 Overall Length 12 cm (4 $\frac{3}{4}$ ")

LRI Diamond Micrometer Knife
 Overall Length 12 cm (4 $\frac{3}{4}$ ")




Diamond Knives
 Diamantmesser
 Bisturí de diamante
 Couteaux de diamant
 Coltelli di diamanti




AK Diamond Knife
 Overall Length 11 cm (4 1/4")



Diamond Micrometer Knife
 Overall Length 13 cm (5")




OP-311-02
 Double Vertical and 45° Cutting Edge
 Each Depth Can Be Adjusted Without
 Steps The Scale is Divided Each
 Hundredth Millimeter




OP-311-03
 Triple Vertical and 15° Cutting Edge
 Each Depth Can Be Adjusted Without Steps
 The Scale is Divided Each Hundredth Millimeter

DIAMAX™


Diamond Knives
 Diamantmesser
 Bisturí de diamante
 Couteaux de diamant
 Coltelli di diamanti

OP
OPHTHALMIC SURGERY INSTRUMENTS



OP-316-01
 Single 1 mm Wide
 30° Cutting Edge




OP-316-02
 Single 1 mm Wide
 45° Cutting Edge



OP-316-03
 Single 1 mm Wide
 Tri-Facet Cutting Edge



OP-316-04
 1 mm Wide Spear
 Shaped Cutting Edge



OP-316-05
 1 mm Wide
 10-Facet Cutting Edge

Diamond Knife
 Overall Length 12.5 cm (5")




Diamond Knives
 Diamantmesser
 Bisturí de diamante
 Couteaux de diamant
 Coltelli di diamanti




Diamond Step Knife
 Overall Length 12.5 cm (5")



OP-320-03
 10-Facet Cutting Edge 1.0 mm Wide
 With Adjustable Cutting
 Depth of 0.15 / 0.2 / 0.25 / 0.3 /
 0.37 / 0.5 / 4.0 mm




OP-320-01
 Single 45° Cutting Edge 1.0 mm Wide
 With Adjustable Cutting
 Depth of 0 / 0.3 / 0.4 / 0.5 / 0.6 /
 0.8 / 1.0 / 2.5 mm



OP-320-02
 Trifacet Cutting Edge 1.0 mm Wide
 With Adjustable Cutting
 Depth of 0 / 0.3 / 0.4 / 0.5 / 0.6 /
 0.8 / 1.0 / 2.5 mm




Crescent Diamond Knife
 Overall Length 11.5 cm (4 1/2")




OP-331-25
 Crescent Shaped 9-Edge
 Vertical Blade 2.0 mm Wide
 Angled

DIAMAX™


Diamond Knives
 Diamantmesser
 Bisturí de diamante
 Couteaux de diamant
 Coltelli di diamanti

OP
OPHTHALMIC SURGERY INSTRUMENTS



OP-340-28
 Straight, 2.8 mm
 Spear Shaped Cutting
 Edge




OP-340-30
 Straight, 3.0 mm
 Spear Shaped Cutting
 Edge



OP-340-32
 Straight, 3.2 mm
 Spear Shaped Cutting
 Edge



OP-341-28
 Angled, 2.8 mm
 Spear Shaped Cutting
 Edge



OP-341-30
 Angled, 3.0 mm
 Spear Shaped Cutting
 Edge



OP-341-32
 Angled, 3.2 mm
 Spear Shaped Cutting
 Edge

Diamond Phaco Knife
 Straight, Spear Shaped
 Overall Length 13 cm (5")


Diamond Phaco Knife
 Angled, Spear Shaped
 Overall Length 13 cm (5")



Diamond Phaco Knife
Overall Length 13 cm (5")




OP-351-25
Angled, 2.5 mm
Spear Shaped 4 - Sided
Cutting Edge



OP-351-30
Angled, 3.0 mm
Spear Shaped 4 - Sided
Cutting Edge




OP-351-28
Angled, 2.8 mm
Spear Shaped 4 - Sided
Cutting Edge



OP-351-32
Angled, 3.2 mm
Spear Shaped 4 - Sided
Cutting Edge



Diamond Phaco Knife
Overall Length 12 cm (4 3/4")




OP-355-27
Angled, 2.75 / 3.2 mm Trapezoid
Shaped 4 - Sided Cutting Edge


DIAMAXTM


Diamond Knives
 Diamantmesser
 Bisturí de diamante
 Couteaux de diamant
 Coltelli di diamanti


OP

OPHTHALMIC SURGERY INSTRUMENTS


OP-357-18

 Angled, 1.8 / 2.2 mm Trapezoid
 Shaped 4 - Sided Cutting Edge

OP-357-20

 Angled, 2.0 / 2.3 mm Trapezoid
 Shaped 4 - Sided Cutting Edge

OP-359-18

 Angled, 1.8 mm
 Spear Shaped 4 - Sided
 Cutting Edge

OP-359-22

 Angled, 2.2 mm
 Spear Shaped 4 - Sided
 Cutting Edge

MICS Diamond Knife
 Overall Length 12 cm (4 ¾")

MICS Diamond Knife
 Overall Length 13 cm (5")

DIAMAXTM



Trephines
Trehinen
Trépanos
Trépans
Trepani



Elliot Scleral Trephine Blade
3.5 cm (1 1/2")

OP-370-10

Ø 1.0 mm

OP-370-15

Ø 1.5 mm

OP-370-20

Ø 2.0 mm



OP-370-90
Elliot Trephine Handle
4.5 cm (1 3/4")



Camellin LASEK Trephine
270° Semi-Sharp Edge 90° Cut Out for Flap Hinge
Overall Length 12 cm (4 3/4")



Camellin LASEK Trephine & Alcohol Retaining Well
Overall Length 12.5 cm (5")

OP-381-08

Ø 8 mm Trehine and Ø 8.5 mm Well

OP-381-08

Ø 8 mm

OP-381-09

Ø 9 mm

OP-381-10

Ø 10 mm

OP-385-08

Ø 9 mm Trehine and Ø 9.5 mm Well


OP-385-09

Ø 10 mm Trehine and Ø 10.5 mm Well

OP-385-10

Ø 10 mm Trehine and Ø 10.5 mm Well

Ophthalmic Punches
 Augenheilkunde Stanze
 Punzones de oftalmología
 Poinçons ophtalmiques
 Punzoni oftalmici



OP-391-00
Gass Corneoscleral Punch

 With 360° Rotatable Cutting Tip 1.5 mm x 0.75 mm Bite
 Overall Length 8.5 cm (3 1/2")

OP-395-00
Luntz-Dodick Trabeculectomy Punch

 With 360° Rotatable Cutting Tip 1 mm x 0.5 mm Bite
 Overall Length 11.5 cm (4 1/2")

OP-397-00
Kelly Descemet's Membrane Punch

 For Glaucoma 19 Gauge / 1.0 mm 0.9 mm x 0.45 mm Bite
 Overall Length 12.5 cm (5")

OP-399-00
Christensen Punctum Punch

 For Opening Stenotic Punctum 1.0 mm Wide Tip
 Overall Length 10 cm (4")

OP
OPHTHALMIC SURGERY INSTRUMENTS



Trephines
Trephinen
Trépanos
Trépans
Trepani



Corneal Trehpine Blade

OP-400-01	OP-400-02	OP-400-03	OP-400-04	OP-400-05	OP-400-06	OP-400-07	OP-400-08
Ø 6.0 mm	Ø 6.5 mm	Ø 6.75 mm	Ø 7.0 mm	Ø 7.25 mm	Ø 7.5 mm	Ø 7.75 mm	Ø 8.0 mm
OP-400-09	OP-400-10	OP-400-11	OP-400-12	OP-400-13	OP-400-14		
Ø 8.25 mm	Ø 8.5 mm	Ø 8.75 mm	Ø 9.0 mm	Ø 9.5 mm	Ø 14.0 mm		



Corneal Trehpine Handle

OP-401-01	OP-401-02	OP-401-03	OP-401-04	OP-401-05	OP-401-06	OP-401-07	OP-401-08
Ø 6.0 mm	Ø 6.5 mm	Ø 6.75 mm	Ø 7.0 mm	Ø 7.25 mm	Ø 7.5 mm	Ø 7.75 mm	Ø 8.0 mm
OP-401-09	OP-401-10	OP-401-11	OP-401-12	OP-401-13			
Ø 8.25 mm	Ø 8.5 mm	Ø 8.75 mm	Ø 9.0 mm	Ø 9.5 mm			



Cannulas
Kanülen
Cánuas
Canules
Cannule

OP

Lacrimal Cannula
OP-410-23

 Malleable Tip
Straight, 23 G

OP-410-27

 Malleable Tip
Straight, 27 G

Lacrimal Cannula
OP-411-23

 Malleable Tip
Curved, 23 G

OP-411-27

 Malleable Tip
Curved, 27 G

OP-412-23
McIntyre Lacrimal Cannula

 Straight, 23 G,
Blunt Tip Dual
0.3 mm Side Ports


OP-413-23
McIntyre Lacrimal Cannula

 Curved, 23 G,
Blunt Tip Dual
0.3 mm Side Ports

OP-415-23
West Lacrimal Cannula

 Blunt Tip, 23 G
0.3 mm Side Port

OP-417-27
DALK Air Injection Cannula

 Flattened, Vaulted
0.2 mm Posterior Port
Near Tip, 27 G

OP-419-27
Rosenwasser DALK Cannula

 Flattened "Cobra Shaped"
Tip 0.2 mm Posterior Port
27 G

Air Injection Cannula

OP-421-27

 Angled at
5 mm, 27 G

OP-421-30

 Angled at
5 mm, 30 G

OP-423-27

 Angled at
7 mm, 27 G

OP-423-30

 Angled at
7 mm, 30 G



Cannulas
Kanülen
Cáñulas
Canules
Cannule



Randolph Cyclodialysis Cannula
Flat Tip Angled at 13 mm



OP-427-19
19 G

OP-427-21
21 G



Viscoelastic Aspirating Cannula
Round Tip Angled at 10 mm

OP-429-19
10 mm, 19 G

OP-429-22
10 mm, 22 G



OP-433-25
Troutman Olive Tip Cannula
Delicate, 25 G

OP-437-27


Flat Tip, Angled at 7 mm
27 G

OP-437-30

Flat Tip, Angled at 7 mm
30 G


Rainin Air Injection Cannula
Flat Tip, Angled at 7 mm, 27 G

Cannulas
Kanülen
Cánuas
Canules
Cannule

OP
OPHTHALMIC SURGERY INSTRUMENTS

OP-439-25
LASEK Alcohol Cannula

Angled at 10 mm

Delicate Olive Shaped Tip, 25 G


OP-441-19
Stevens Sub Tenon's Anesthesia Cannula

Gently Curved, Flat Tip, 19 G


OP-443-23
Atkinson Retrobulbar Needle

Straight, 35 mm Long Sharp Tip 23 G


OP-445-23
Corbin Sub Tenon's Anesthesia Cannula

Smooth Rounded Tip With 0.3 mm

Side Port 23 G


OP-447-24
Fukasaku Sub Tenon's Anesthesia Cannula

Gently Curved Smooth Tip, 24 G

Cannulas
 Kanülen
 Cáñulas
 Canules
 Cannule

OP-451-23
Simcoe Irrigating - Aspirating Cannula

 15 mm Long With 0.3 mm Aspirating Side Port, 23 G
 Aspiration Through Tubing Hub and Irrigation Through
 Luer-Lock Hub


OP-453-23
Simcoe (Reversed) Irrigating - Aspirating Cannula

 15 mm Long With 0.3 mm Aspirating Side Port, 23 G
 Aspiration Through Luer-Lock Hub and Irrigation Through
 Tubing Hub

OP-455-23
Simcoe Irrigating - Aspirating Cannula

 Curved Right, 23 G
 U-Shaped, With 0.3 mm Anterior Aspiration Port
 Aspiration Through Tubing Hub and Irrigation
 Through Luer-Lock Hub

OP-457-23
Simcoe Irrigating - Aspirating Cannula

 Curved Left, 23 G
 U-Shaped, With 0.3 mm Anterior Aspiration Port
 Aspiration Through Tubing Hub and Irrigation
 Through Luer-Lock Hub

OP-459-23
Gills Irrigating - Aspirating Cannula

 23 G, Angled Tip, Front Openings
 With 0.5 mm Recessed Irrigating Port

OP-462-22
Kelman Cystotome

23 G, Sharp Cutting Edge With Sharp Point

Cannulas
Kanülen
Cánuas
Canules
Cannule

OP

OP-463-22
Kratz Cystotome

22 G, Sharp Cutting Edge With Sharp Point


OP-465-23
McIntyre Cystotome

23 G, Angled at 10 mm, Sharp Point


Jansen Capsule Polisher Cannula
OP-467-23

 23 G, Gently Curved Shaft With
Lower Half Diamond Dusted Olive Tip

OP-471-19
Bishop-Harmon Anterior Chamber Irrigator

 23 G, Angled Shaft With Flat Blunt Tip Including Adapter
and Blue Silicone Bulb

OP-475-25
Hydrodissection Cannula

Angled at 8 mm, Flat Tip, 25 G


OP-477-26
Seeley Hydrodissection Cannula

Angled at 12 mm, Flat Tip, 25 G

Cannulas
 Kanülen
 Cábulas
 Canules
 Cannule


Knolle Anterior Chamber Irrigation Cannula
 Angled at 2 mm With Front Opening Tip



OP-481-23
 23 G



OP-481-27
 27 G



OP-481-30
 30 G




OP-485-27
Sauter Hydrodissection Cannula
 Vertically Flattened Tip With Front Opening, 27 G



OP-487-19
Bracken Irrigation Cannula
 Smooth Beveled Tip With Front
 Opening, 19 G



OP-489-23
Irrigation Cannula
 Curved Shaft With Blunt Tip and
 Dual 0.4 mm Side Ports, 23 G



OP-491-23
Alio LASIK Irrigation Cannula
 With One Front and Two Side Ports
 23 G



OP-493-21
Irrigation Cannula
 Curved Shaft With Blunt Tip and Dual
 0.5 mm Side Ports, 21 G



Cannulas
Kanülen
Cánuas
Canules
Cannule

OP

**OP-495-26****Spatulated LASIK Cannula**

Flattened Spatula Shaped Tip, 26 G

OPHTHALMIC SURGERY INSTRUMENTS

**OP-501-23****Fukasaku LASIK Cannula**

Flattened Low Profile Tip With Front Opening, 23 G

**OP-503-25****Carter LASIK Cannula**



Flattened Low Profile Tip With Front Opening, 25 G

**OP-507-25****Banaji LASIK Cannula**

6 Irrigating Ports With Bullet Shaped Tip, 25 G




Cannulas
Kanülen
Cáñulas
Canules
Cannule



Cannulas
Kanülen
Cánuas
Canules
Cannule

OP
OPHTHALMIC SURGERY INSTRUMENTS


OP-527-23
Charleux Cortex Extractor
Curved Shaft With Beveled Irrigation
Port 23 G



OP-527-25
De La O Aspirating Cannula
Angled at 12 mm With 0.3 mm
Side Port, 25 G



OP-529-05
5 mm Long Tip



OP-529-07
7 mm Long Tip



OP-529-08
8 mm Long Tip



OP-529-10
10 mm Long Tip



OP-531-23
Fukasaku Hydro-Chop Cannula
Vertically Flattened Tip, Dual Side Ports, 23 G



OP-532-27
Chang Hydrodissection Cannula
Flat Tip With Beveled Opening, 27 G



OP-533-23
Drews Irrigation Vectis
One Front and Two Backward
Opening Ports, 23 G



OP-535-20
Sheets Irrigation Vectis
Three Front Opening Ports, 20 G

Cannulas
 Kanülen
 Cáñulas
 Canules
 Cannule


OP-537-23
Knolle-Pearce Irrigation Vectis
 Three Front Opening Ports, 23 G



OP-539-25
Simcoe Irrigation Lens Nucleus Loop
 Finely Serrated Tip With Two Front Opening
 Ports, Right, 25 G



OP-541-25
Simcoe Irrigation Lens Nucleus Loop
 Finely Serrated Tip With Two Front Opening
 Ports, Left, 25 G



Alfonso Nucleus Irrigating Spoon


OP-542-25
 2.5 mm Wide Spoon With Single
 0.3 mm Front Opening Port

OP-542-40
 4.0 mm Wide Spoon With Single
 0.3 mm Front Opening Port



OP-543-22
Drof Infusion Cannula
 18 Gauge Shaft With
 22 Gauge 3 mm Long Tip

Cannulas
Kanülen
Cábulas
Canules
Cannule

OP
OPHTHALMIC SURGERY INSTRUMENTS

OP-549-01 T
Irrigation Handpiece

Titanium

 For Bimanual Technique, 21 G, Dual 2 x 0.5 mm Side Ports
 Overall Length 10.5 cm (4")

OP-549-02 T
Aspiration Handpiece

Titanium

 For Bimanual Technique, 21 G, 2 x 0.5 mm Side Port
 Overall Length 10.5 cm (4")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Ophtalmologie Instruments
 Oftalmologia Strumenti



Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Ophtalmologie Instruments
 Oftalmologia Strumenti

OP-550-00
 Fig. 0, Ø 1.5 mm

OP-550-01
 Fig. 1 Ø 1.75 mm

OP-550-02
 Fig. 2, Ø 2.25 mm

OP-550-03
 Fig. 3, Ø 2.5 mm

OP-550-04
 Fig. 4, Ø 3.5 mm

Meyerhofer Chalazion Curette
 Overall Length 12.5 cm (5")

OP-552-10
 Fig. 0
 1 x 1.7 mm

OP-552-15
 Fig. 1
 1.5 x 2.5 mm

OP-552-20
 Fig. 2
 2.0 x 3.3 mm

OP-552-25
 Fig. 3
 2.5 x 4.0 mm

OP-552-30
 Fig. 4
 3.0 x 5.0 mm

OP-552-35
 Fig. 5
 3.5 x 5.7 mm

Hebra Chalazion Curette
 Overall Length 12.5 cm (5")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Optalmologie Instruments
 Oftalmologia Strumenti

OP

OPHTHALMIC SURGERY INSTRUMENTS


OP-571-12
Nightingale Capsule Polishing Curette

 Semi-Sharp Posterior Edge Ø 1.75 mm Open Ring
 Overall Length 12.5 cm (5")

OP-573-12
Shepherd Capsule Polishing Curette

 Semi-Sharp Anterior & Posterior Edges Ø 1.75 mm Open Ring
 Overall Length 12.5 cm (5")


Top View

Side View

OP-587-13
Blumenthal Conjunctiva Dissector

 Ø 2.0 mm Blunt Disc Shaped Tip
 Overall Length 13 cm (5")

OP-591-12
Buratto LASIK Flap Protector

 12.0 mm Wide Blade
 Overall Length 12.5 cm (5")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Ophtalmologie Instruments
 Oftalmologia Strumenti



OP-594-11
Buratto LASIK Flap Protector

Double Ended With 12.0 mm Wide Blades Angled at 15°
 Overall Length 11.5 cm (4 1/2")



OP-597-11
Mannis-Buratto LASIK Flap Protector

Double Ended With 12.0 mm Wide Blades Angled at 45°
 Overall Length 11.5 cm (4 1/2")



OP-599-13
Rowen LASIK Flap Protector/Spatula

12.0 mm Wide Protector Blade and 8.0 mm Long Spatula Blade
 Overall Length 13 cm (5")

Ophthalmology Instruments
Ophthalmologie Instrumente
Oftalmología Instrumentos
Ophtalmologie Instruments
Oftalmologia Strumenti

OP

OPHTHALMIC SURGERY INSTRUMENTS


OP-611-12

Sloane LASIK Micro Hoe
2.0 mm Wide Semi-Circular
Tip With Beveled Edges
Overall Length 12 cm (4 $\frac{3}{4}$)


OP-615-11

Sloane Epi Peeler
1 x 5 mm Spatula Shaped
Semi Sharp Tip Edge
Overall Length 11.5 cm (4 $\frac{1}{2}$)


OP-619-13

Sloane LASIK Micro Hoe & Epi Peeler
2.0 mm Wide Semi-Circular Tip With Beveled Edges
1 x 5 mm Spatula Shaped Semi Sharp Tip Edge
Overall Length 13.5 cm (5 $\frac{1}{4}$)

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Ophtalmologie Instruments
 Oftalmologia Strumenti



OP-623-11
Sloane LASIK Epi-Dissector
 7.0 mm Long Spatula Shaped Tip
 With Semi-Sharp Edge
 Overall Length 11.5 cm (4 ½")



OP-627-12
Sloane LASIK Flap Repository
 Gently Curved Bi-Convex Spatulated Tip
 With Beveled Smooth Edge
 Overall Length 12 cm (4 ¾")



OP-629-01
DALK Corneal Dissector
 12 mm Long Flat Vaulted Spatula
 With Beveled Tip
 Overall Length 12.5 cm (5")



OP-629-02
Rosenwasser DALK Trisector
 10 mm Long Triangular Shaft With Spatula Shaped Tip
 Overall Length 12.5 cm (5")



OP-629-03
Bonfadini DALK Spatula
 Vaulted 1.0 mm Wide 8 mm Long Blade Spatula
 Shaped Blade With Central Groove
 Overall Length 12.5 cm (5")



Ophthalmology Instruments
Ophthalmologie Instrumente
Oftalmología Instrumentos
Ophthalmologie Instruments
Oftalmologia Strumenti

OP

OPHTHALMIC SURGERY INSTRUMENTS

**OP-629-14****Dastoor Pupil and Iris Repositor**

Double Ended 1 mm Wide Notched Tip on One Side and Spatula Shaped Tip on Other side
Overall Length 13 cm (5")

**OP-630-01****Culler Iris Spatula**

Straight 1 mm Wide and 28 mm Long Spatula
Overall Length 13 cm (5")

**OP-630-02****Culler Iris Spatula**

Straight 2 mm Wide and 28 mm Long Spatula
Overall Length 13 cm (5")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Ophtalmologie Instruments
 Oftalmologia Strumenti



OP-632-13

Wheeler Double-End Spatula

Double Ended Slightly Curved Spatula Shaped Blades (2 x 33 mm Blade on One Side and 2 x 24 mm Blade on The Other Side)
 Overall Length 13 cm (5")



Barraquer Iris Spatula

Overall Length 12 cm (4 3/4")

OP-635-01

Vaulted 0.25 mm Wide and 15.0 mm Long Blade With Blunt Tip



OP-637-12

Gimbel Nucleus Spatula

30° Angled Shaft, 9 mm Long Ø 0.25 mm Round Blade
 Overall Length 12.5 cm (5")

OP-635-02

Vaulted 0.50 mm Wide and 15.0 mm Long Blade With Blunt Tip

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Optalmologie Instruments
 Oftalmologia Strumenti

OP

OPHTHALMIC SURGERY INSTRUMENTS


OP-639-12
Escaf Nucleus Manipulator

With 1 mm Wide and 1.4 mm Long Tip
 Angled at 45° - 13 mm Long from Bend
 to Tip Shaft
 Overall Length 12 cm (4 3/4")


OP-640-12
Knolle Nucleus Spatula

0.5 mm Wide, Flat Malleable Spatula
 Overall Length 12.5 cm (5")


OP-643-01
OP-643-02
Koch Nucleus Spatula

Overall Length 11.5 cm (4 1/2")


OP-643-01

Slightly Curved 0.8 mm Wide Duckbill Shaped Tip


OP-643-02

Slightly Curved 0.8 mm Wide Duckbill Shaped Tip
 With Notches

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Ophtalmologie Instruments
 Oftalmologia Strumenti


OP-645-11
Akahoshi Nucleus Sustainer

Ø 0.4 mm Ball Shaped Tip 1.8 mm
 Angled With 10 mm Long Shaft
 Overall Length 11.5 cm (4 ½")


OP-647-12
Escaf Nucleus Sustainer

Ø 0.5 mm Ball Shaped Tip With Inferior Textured Surface 1.8 mm Angled With 10 mm Long Shaft
 Overall Length 12 cm (4 ¾")


Lieberman Nucleus Manipulator

Overall Length 12 cm (4 ¾")

OP-649-01

1.5 mm Curved Right Tip With 10 mm Long Shaft from Bend to Tip

OP-649-02


1.5 mm Curved Left Tip With 10 mm Long Shaft from Bend to Tip

OP-649-01
OP-649-02


Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Optalmologie Instruments
 Oftalmologia Strumenti

OP

OPHTHALMIC SURGERY INSTRUMENTS


Lieberman Nucleus Manipulator and Spatula

Overall Length 14 cm (5 1/2")

OP-651-01

1.5 mm Curved Right Tip With 10 mm Long Shaft
 from Bend to Tip on One Side and 1 mm Wide Spatula
 Shaped Tip on Other Side

OP-651-02

1.5 mm Curved Left Tip With 10 mm Long Shaft
 from Bend to Tip on One Side and 1 mm Wide Spatula
 Shaped Tip on Other Side

Fukasaku-Lieberman Nucleus Manipulator and Hook

Overall Length 13.5 cm (5 1/4")

OP-653-01

For Right Hand Use 3.5 mm Long Micro Finger Chopper
 on One Side & 1.25 Long Blunt Tip Hook With 0.8 mm
 Long Wedge Shaped Inferior Edge & Notch on Front
 Side on Other Side

OP-653-02

For Right Hand Use 3.5 mm Long Micro Finger Chopper
 on One Side & 1.25 Long Blunt Tip Hook With 0.8 mm
 Long Wedge Shaped Inferior Edge & Notch on Front
 Side on Other Side

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Ophtalmologie Instruments
 Oftalmologia Strumenti


OP-657-13
Masket Phaco Spatula

Combination of Ø 1 mm
 Disc Shaped Spatula & 0.8 mm
 Long Rosen Type Phaco Splitter Tip
 Overall Length 13 cm (5")


OP-659-01

For Right Hand Use 3.5 mm Long Micro Finger Chopper on One Side & 0.9 mm Long Quick Chop Tip on Other Side

OP-659-02

For Left Hand Use 3.5 mm Long Micro Finger Chopper on One Side & 0.9 mm Long Quick Chop Tip on Other Side


OP-661-11
Minardi Phaco Chopper

30° Angled Shaft With 0.7 mm Long Quick Chop Tip Angle to Tip Length of 10 mm
 Overall Length 11.5 cm (4 ½")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Optalmologie Instruments
 Oftalmologia Strumenti

OP


OPHTHALMIC SURGERY INSTRUMENTS


OP-663-13
Hwang Phaco Chopper

Double 1 mm Long Vertical Chopper On One Side and 1 mm Wide Blunt Chisel Tipped Spatula On Other Side, Angle to Tip Length of 10 mm on Both Sides
 Overall Length 13.5 cm (5 1/4")


OP-665-11
Agarwal Karate Phaco Chopper

0.9 mm Long Tip With 0.6 mm Semi-Sharp Inferior Edge & Pointed Tip Angle to Tip Length of 11 mm
 Overall Length 11.5 cm (4 1/2")


OP-667-11
Chang Micro Finger Chopper

3.5 mm Long Curved Finger Shaped Tip With Defined Inner Edge for Horizontal Chopping, Angle to Tip Length of 10 mm
 Overall Length 11.5 cm (4 1/2")


Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Ophtalmologie Instruments
 Oftalmologia Strumenti



OP-669-12
Brierley Nucleus Splitter
 1 mm Long Splitter With Incision
 Sealing Shaft Angle to Tip Length of
 12 mm
 Overall Length 12 cm (4 $\frac{3}{4}$ ")



Dodick-Kammann Nucleus Chopper
 Overall Length 12 cm (4 $\frac{3}{4}$ ")
OP-670-15
 Straight Shaft With 1.5 mm Long Chopper
 Tip For Bimanual Nucleus Division
OP-670-20
 Straight Shaft With 2.0 mm Long Chopper
 Tip For Bimanual Nucleus Division



Dodick-Kammann Nucleus Chopper
 Overall Length 12 cm (4 $\frac{3}{4}$ ")
OP-671-15
 Angled Shaft With 1.5 mm Long Chopper
 Tip For Bimanual Nucleus Division, Angle to Tip Length
 of 12 mm
OP-671-20
 Angled Shaft With 2.0 mm Long Chopper
 Tip For Bimanual Nucleus Division, Angle to Tip Length
 of 12 mm

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Optalmologie Instruments
 Oftalmologico Strumenti

OP
OPHTHALMIC SURGERY INSTRUMENTS

OP-673-12
Escaf Nucleus Divider

1.25 mm Long & 0.5 mm Wide Wedge Shaped Tip for Vertical Chopping
 Overall Length 12 cm (4 $\frac{3}{4}$ ")


OP-675-12
Neuhann Nucleus Divider


0.5 mm Wide Strongly Curved Spatula Shaped Tip
 Overall Length 12 cm (4 $\frac{3}{4}$ ")


OP-677-11
Nagahara Nucleus Lens Spatula & Chopper

Chopper With 1.5 mm Long Tip With 0.9 mm Long Inferior Wedge and Blunt Tip on One Side and Spatula With Flat Inferior Edge on Other Side Angle to Tip Length of 10 mm on Both Sides
 Overall Length 11 cm (4 $\frac{1}{4}$ ")


OP-678-11
Nagahara Nucleus Chopper

Chopper With 1.5 mm Long Tip With 0.9 mm Long Inferior Wedge and Blunt Tip Angle to Tip Length of 10 mm
 Overall Length 11 cm (4 $\frac{1}{4}$)


OP-679-11
Fukasaku Snapper Hook

1.25 mm Long Blunt Tip With 0.8 mm Long Wedge Shaped Inferior Edge Angle to Tip Length of 9 mm
 Overall Length 11.5 cm (4 $\frac{1}{2}$ ")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Ophtalmologie Instruments
 Oftalmologia Strumenti


OP-681-11

Fukasaku Snapper Hook & Small Pupil
 1.25 mm Long Blunt Tip With 0.8 mm
 Long Wedge Shaped Inferior Edge &
 Notch on Front Angle to Tip Length of
 9 mm
 Overall Length 11.5 cm (4 1/2")


OP-683-11

Rosen Phaco Splitter
 0.8 mm Long Sharp 45°
 Inferior Cutting Edge Angle
 to Tip Length of 9 mm
 Overall Length 11.5 cm (4 1/2")


OP-685-11

Rosen Phaco Splitter
 0.8 mm Long Wedge
 Shaped 60° Inferior
 Cutting Edge With Blunt Tip
 Angle to Tip Length of 9 mm
 Overall Length 11.5 cm (4 1/2")


OP-687-13

Rowen Phaco Spatula
 Double Ended With Angled
 Shafts, 9 mm from Angle to
 Tip For Rosen Chopper on
 One Side and 7 mm Angle To
 Tip for Bechert Type Nucleus
 Rotator on Other Side
 Overall Length 13 cm (5")

Ophthalmology Instruments
Ophthalmologie Instrumente
Oftalmología Instrumentos
Ophthalmologie Instruments
Oftalmologia Strumenti

OP

OPHTHALMIC SURGERY INSTRUMENTS

**OP-689-01****Bechert Nucleus Rotator**

0.78 mm Horizontally Forked Tip With
Angled to Tip Length of 7 mm
Overall Length 11.5 cm (4 1/2")

**OP-689-02****Jaffe-Bechert Nucleus Rotator**

0.78 mm Horizontally Forked Tip With
Angled to Tip Length of 9 mm
Overall Length 11.5 cm (4 1/2")

**OP-691-11****Clayman Nucleus Rotator**

0.5 mm Wide Flat Spatula With V-Notched Tip
Angled to Tip Length of 10 mm
Overall Length 12 cm (4 3/4")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Ophtalmologie Instruments
 Oftalmologia Strumenti


OP-693-11

Drysdale Nucleus Manipulator
 0.8 mm x 2.4 mm Paddle
 Shaped Tip With Angle to Tip
 Length of 10 mm
 Overall Length 11 cm (4 1/4")


OP-695-11

Pisacano Nucleus Rotator
 1 mm Horizontal Flat Disc
 Shaped Tip With Angled to
 Tip Length of 11 mm
 Overall Length 11.5 cm (4 1/2")


OP-697-11

Pisacano Nucleus Rotator
 1 mm Vertical Flat Disc Shaped
 Tip With Angled to Tip Length
 of 11 mm
 Overall Length 11.5 cm (4 1/2")



OP-699-11

Siepser Endocapsular Controller
 Ø 0.7 mm Trumpet Shaped Tip
 With Angled to Tip Length of
 12 mm
 Overall Length 11.5 cm (4 1/2")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Optalmologie Instruments
 Oftalmologica Strumenti

OP

OPHTHALMIC SURGERY INSTRUMENTS


OP-700-12

Anis Nucleus Manipulator
 Straight Shaft & Curved Tip
 With Ø 0.5 mm Ball at Tip
 Overall Length 12 cm (4 ¾")

OP-701-12

Anis Nucleus Manipulator
 Angled Shaft And Curved Tip
 With Ø 0.5 mm Ball at Tip,
 Angle to Tip Length of 10 mm
 Overall Length 12 cm (4 ¾")

OP-703-12

Tennant Nucleus Rotator
 Ø 0.6 mm Ball Shaped Tip
 With Angled to Tip Length
 of 10 mm
 Overall Length 12 cm (4 ¾")

OP-705-11

Lippman Nucleus Spatula
 Ø 0.5 mm Shaft With
 Blunt Forked Tip Angle to Tip
 Length of 10 mm
 Overall Length 11.5 cm (4 ½")



Ophthalmology Instruments
Ophthalmologie Instrumente
Oftalmología Instrumentos
Ophtalmologie Instruments
Oftalmologia Strumenti


OP-707-11
Girard Synechia Spatula

Vaulted 0.5 mm Wide Blade
With Semi Sharp Tip
Angle to Tip Length of 10 mm
Overall Length 11.5 cm (4 1/2")



OP-709-12
Obstbaum Synechia Spatula

0.25 mm Wide Flat Spatula With
Semi Sharp Tip Angled to Tip Length
of 10 mm
Overall Length 12 cm (4 3/4")


OP-711-13
Moran LASIK Enhancement Spatula

Double Ended Fukasaku Spatula On
One Side & Miniature Paton Spatula
on Other Side
Overall Length 13.5 cm (5 1/4")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Optalmologie Instruments
 Oftalmologia Strumenti

OP
OPHTHALMIC SURGERY INSTRUMENTS

OP-713-14

Hersh IntraLase Flap Spatula
 Double Ended 2 x 6 mm and
 2 x 10 mm Spatula With
 Beveled Blade Edges
 Overall Length 14 cm (5 1/4")


OP-715-13

Binder Flap Separator & Elevator
 1 mm Wide Spatula With 0.7 mm Long
 Rounded Tip (0.1 mm x 0.3 mm), Angle
 to Tip Length of 10 mm
 Overall Length 13 cm (5")


OP-717-12

Fukasaku LASIK Spatula
 1.2 mm Wide Gently Curved
 Spatula With Round Anterior
 Edge and Flat Posterior Edge,
 Angled to Tip Length of 12 mm
 Overall Length 12 cm (4 3/4")


Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Ophtalmologie Instruments
 Oftalmologia Strumenti



OP-719-12

Guimaraes LASIK Spatula

1.2 mm Wide Gently Curved Bi-Convex
 Spatula With Semi-Sharp Edge
 Overall Length 12 cm (4 ¾")



OP-720-13

Katzen LASIK Flap Spatula

Double Ended With 2 mm Wide Bi-Convex
 Notched Spatula on One Side And Ø 0.2 mm
 Modified Sinskey Hook on Other Side
 Overall Length 13.5 cm (5 ¼")



OP-723-12

Maddox LASIK Spatula

Double Ended Spatula Combines a 2 mm Wide
 Flat And a Ø 0.25 mm Cylindrical Spatula With
 7 mm from Angle to Tip Length On Both Sides
 Overall Length 12.5 cm (5")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Optalmologie Instruments
 Oftalmologia Strumenti

OP

OPHTHALMIC SURGERY INSTRUMENTS



OP-725-11
LASIK Flap Manipulator

Delicate Ø 0.3 mm Olive Tip Manipulator
 With Angle to Tip Length of 10 mm
 Overall Length 11.5 cm (4 1/2")



OP-726-11
Hirschmann Lens Spatula
 0.6 mm Wide Spatula With 0.2 mm
 Notches on Both Sides
 Overall Length 11.5 cm (4 1/2")



OP-727-12
Jaffe Lens Spatula
 0.5 mm Thick Smooth Flat Spatula With Angled
 to Tip Length of 10 mm
 Overall Length 12 cm (4 3/4")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Ophtalmologie Instruments
 Oftalmologia Strumenti



OP-728-12

Straight Ø 0.25 mm Hourglass
 Shaped Tip With 0.5 mm Height



OP-729-12

Angled Ø 0.25 mm Hourglass
 Shaped Tip With 0.5 mm Height
 With Angled to Tip Length of
 10 mm

Lester IOL Manipulator
 Overall Length 12 cm (4 ¾")



OP-731-11

Lester IOL Manipulator
 Ø 0.25 mm Hourglass Shaped Tip
 With 0.5 mm Height With Angled to Tip
 Length of 10 mm Vaulted Shanks
 Overall Length 11.5 cm (4 ½")



OP-733-12

Hunkeler IOL Manipulator
 Ø 0.3 mm Hemi-Ball Shaped Tip With Angle
 to Tip Length of 10 mm
 Overall Length 12 cm (4 ¾")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Optalmologie Instruments
 Oftalmologia Strumenti

OP
OPHTHALMIC SURGERY INSTRUMENTS


OP-735-12
Alio-Kelman IOL Hook
 Ø 0.2 mm Tip With
 0.5 mm Height With Angle
 to Tip Length of 10 mm
 Overall Length 12 cm (4 ¾")



OP-738-12
Aker Lens Pusher
 0.8 mm Smooth, Fish
 Mouth Shaped Jaws
 Overall Length 12.5 cm (5")




OP-739-11
Bechert Lens Pusher
 0.78 mm Vertically Forked
 Tip With Angled to Tip
 Length of 10 mm
 Overall Length 11.5 cm (4 ½")



OP-741-13
Elschnig Cyclodialysis Spatula
 Blunt Tip 1 mm Wide Spatula
 With Angle To Tip Length of 10 mm
 Overall Length 11.5 cm (4 ½")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Ophtalmologie Instruments
 Oftalmologia Strumenti


OP-743-13
Castroviejo Synechia Spatula

Double Ended Spatula 0.5 mm Wide With 0.25 mm Thickness of Blades, 10 mm From Angle to Tip on One Side and 15 mm From Angle to Tip on Other Side, Semi Sharp Tips Overall Length 13 cm (5")


OP-745-25
Castroviejo Cyclodialysis Spatula

Overall Length 13 cm (5")

OP-745-50

Double Ended Spatula 0.25 mm Wide With 0.25 mm Thickness 10 mm From Angle to Tip on One Side and 15 mm From Angle to Tip on Other Side, Blunt Tip

OP-745-50

Double Ended Spatula 0.25 mm Wide With 0.5 mm Thickness 10 mm From Angle to Tip on One Side and 15 mm From Angle to Tip on Other Side, Blunt Tip


OP-749-12
Barrett Nucleus Expressor

Curved Shaft With Ø 2 mm Bulbous Tip Overall Length 12 cm (4 ¾")



Ophthalmology Instruments
Ophthalmologie Instrumente
Oftalmología Instrumentos
Ophthalmologie Instruments
Oftalmologia Strumenti

OP

OPHTHALMIC SURGERY INSTRUMENTS

**Lewis Lens Loop**

Overall Length 13 (5")

OP-750-01

Smooth 5 mm Wide Curved Loop Shaped Tip

OP-750-02

Smooth 6 mm Wide Curved Loop Shaped Tip

**OP-752-13****New Orleans Lens Loop**

Smooth 5 mm Wide Slightly Curved

Loop Shaped Tip

Overall Length 13.5 cm (5 1/4")

**OP-754-13****Wilder Lens Loop**

Serrated 4 mm Wide Curved Loop

Shaped Tip

Overall Length 13.5 cm (5 1/4")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Ophtalmologie Instruments
 Oftalmologia Strumenti


OP-757-15
Kirby Lens Loop and Muscle Hook

Double Ended With 5 mm Wide Loop Shaped Tip on One Side and 12 mm Muscle Hook on Other Side With 2.0 mm Bulbous Tip
 Overall Length 15 cm (6")


OP-758-14
Morrison Lens Loop and Probe

Double Ended With 4.5 mm Wide Loop Shaped Tip on One Side and Ø 0.8 mm Olive Tip Probe on Other Side
 Overall Length 14.5 cm (5 ¾")


OP-760-13
Kansas Lens Loop

Serrated 3.0 mm Wide Slightly Curved Loop Shaped Tip
 Overall Length 13.5 cm (5 ¼")

Ophthalmology Instruments
Ophthalmologie Instrumente
Oftalmología Instrumentos
Ophtalmologie Instruments
Oftalmologia Strumenti

OP

OPHTHALMIC SURGERY INSTRUMENTS

**OP-761-13****Keener-Arlt Lens Loop**

Serrated 3.0 mm Wide Curved

Loop Shaped Tip

Overall Length 13.5 cm (5 1/4")

**OP-765-10****Alfonso Nucleus Trisector**

Oval 2.5 mm Wide Semi Sharp Fenestrated

Blade To Be Used With OP-775-02

Overall Length 10.5 cm (4 1/4")

**OP-767-10****Bucher-Kansas Nucleus Trisector**

Triangular 2.5 mm Wide Semi Sharp

Fenestrated Blade To Be Used With OP-775-03

Overall Length 10.5 cm (4 1/4")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Ophtalmologie Instruments
 Oftalmologia Strumenti



OP-769-11

Kansas Nucleus Trisector

Triangular 4 mm Wide Semi Sharp Fenestrated
 Blade To Be Used With OP-773-12
 Overall Length 11 cm (4 1/4")



OP-771-11

Kansas Nucleus Bisector

10 mm Long Semi Sharp Blade Used
 With OP-773-12 For Small Manual Incision Surgery
 Overall Length 11.5 cm (4 1/4")



OP-773-12

Kansas Nucleus Vectis

3.0 mm Wide Spatula Shaped Blade
 Overall Length 12 cm (4 3/4")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Optalmologie Instruments
 Oftalmologia Strumenti

OP

OPHTHALMIC SURGERY INSTRUMENTS



Alfonso Cutting Platform
 Overall Length 12.5 cm (5")

OP-775-02

2.0 mm Wide Blade Used as Cutting
 Block For Nucleus Bisector and Trisector

OP-775-03

3.0 mm Wide Blade Used as Cutting
 Block For Nucleus Bisector and Trisector



OP-779-12
Giannetti Nucleus Capture Spoon
 Inclined, Arcuate 2.6 mm Wide Spoon
 Overall Length 12.5 cm (5")



OP-781-12
Giannetti Nucleus Capture Spatula
 2 mm Wide Blunt Tip Arcuate Nucleus
 Capture Spatula
 Overall Length 12.5 cm (5")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Ophtalmologie Instruments
 Oftalmologia Strumenti



OP-783-12
Giannetti Nucleus Capture Hook
 Straight Hook For Nucleus Manipulation
 Overall Length 12.5 cm (5")



OP-790-13
Culler Lens Spoon
 5 x 9 mm Spoon Shaped Tip
 Overall Length 13 cm (5")



OP-793-14
Paton Spatula & Spoon
 6 mm Wide and 15 mm Long Grooved Spoon Blade and
 2 x 31 mm Spatula Shaped Blade With Last 7 mm Angled Tip
 Overall Length 14 cm (5 1/2")


Ophthalmology Instruments
Ophthalmologie Instrumente
Oftalmología Instrumentos
Ophtalmologie Instruments
Oftalmologia Strumenti

OP

OPHTHALMIC SURGERY INSTRUMENTS

**OP-795-04**

Rosenwasser Lamella Donor Shovel
9 mm Long and 4 mm Wide Plate With Retaining
Rim on Two Sides
Overall Length 11.5 cm (4 1/2")

**OP-795-05**

Rosenwasser Lamella Donor Shovel
9 mm Long and 5 mm Wide Plate With Retaining
Rim on Three Sides
Overall Length 11.5 cm (4 1/2")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Ophtalmologie Instruments
 Oftalmologia Strumenti



OP-796-00
Carter Sphere Introducer
 With Spring Loaded Plunger To Hold and
 Introduce Eye Spheres
 Overall Length 12.5 cm (5")



OP-798-00
Foster Enucleation Snare
 With Ratcheting Wing-Nut Handle
 Overall Length 14.5 cm (5 1/2")



OP-798-90
Twisted Wires For Foster Enucleation Snare
 Pack of 12



OP-799-13
Eye Enucleation Forcep
 Curved, S Shaped, Longitudinally Serrated
 Overall Length 13.5 cm (5 1/4")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Optalmologie Instruments
 Oftalmologia Strumenti

OP

OPHTHALMIC SURGERY INSTRUMENTS


OP-800-01

Bunge Evisceration Spoon
 Flat and Smooth Ø 10.0 mm Spoon
 Shaped Blade
 Overall Length 14 cm (5 ½")

OP-800-02

Bunge Evisceration Spoon
 Flat and Smooth Ø 12.0 mm
 Spoon Shaped Blade
 Overall Length 14 cm (5 ½")

OP-802-15

Wells Enucleation Spoon
 21 mm Wide Spoon With
 5 x 10 mm Notch
 Overall Length 15 cm (6")

OP-803-00

Dastoor Erisophake (Bell)
 Overall Length 6 cm (2 ¼")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Ophtalmologie Instruments
 Oftalmologia Strumenti




OP-804-12
Graether Collar Button

Straight, Ø 0.6 mm Collar Button Shaped Micro Iris Retractor
 Overall Length 12.5 cm (5")



OP-805-12
Graether Collar Button

Angled, Ø 0.6 mm Collar Button Shaped Micro Iris Retractor With Angle to Tip Length of 10 mm
 Overall Length 12.5 cm (5")




OP-811-11
Z Hook Manipulator

Ø 0.2 mm Z Shaped Hook With Angle to Tip Length of 10 mm
 Overall Length 11.5 cm (4 ½")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Optalmologie Instruments
 Oftalmologia Strumenti

OP

OPHTHALMIC SURGERY INSTRUMENTS


OP-813-11

Sinskey (Reversed) Lens Manipulating Hook
 Ø 0.2 mm Blunt Tip With 0.35 mm Height
 With Angle to Tip Length of 10 mm
 Overall Length 11.5 cm (4 1/2")


OP-815-12

Bonn Micro Iris Hook
 Curved Hook With Ø 0.15 mm Blunt Tip
 Overall Length 12.5 cm (5")



OP-817-12

Shepard Iris Hook
 1.0 mm Wide Ø 0.2 mm Blunt Tip
 Overall Length 12.5 cm (5")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Ophtalmologie Instruments
 Oftalmologia Strumenti


OP-819-01


Jaffe Lens Manipulating Hook
 Ø 0.15 mm Hook and 0.35 mm Height,
 Flat Handle
 Overall Length 11.5 cm (4 ½")


OP-819-02

Clayman-Jaffe Lens Manipulating Hook
 Ø 0.15 mm Hook and 0.35 mm Height With
 Angle to Tip Length of 10 mm Flat Handle
 Overall Length 11.5 cm (4 ½")


OP-821-01

Sinskey Lens Manipulating Hook
 Ø 0.25 mm Blunt Tip With 0.35 mm
 Height, Flat Handle
 Overall Length 11.5 cm (4 ½")


OP-821-02

Harris-Sinskey Lens Manipulating Hook
 Ø 0.25 mm Blunt Tip With 0.35 mm Height With
 Angle to Tip Length of 10 mm Flat Handle
 Overall Length 11.5 cm (4 ½")


OP-823-12


Lewicky Lens Manipulating Hook
 Vaulted Shaft, Ø 0.2 mm Blunt Tip With Angle
 to Tip Length of 11 mm
 Overall Length 12 cm (4 ¾")



Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Optalmologie Instruments
 Oftalmologia Strumenti

OP

OPHTHALMIC SURGERY INSTRUMENTS


OP-825-11

Hirschmann Iris Hook
 Ø 0.2 mm Blunt Round Tip
 Angled to Tip Length of 10 mm
 Overall Length 11 cm (4 1/4")

OP-826-01

Sinskey Lens Manipulating Hook
 Ø 0.25 mm Blunt Tip With 0.35 mm
 Height
 Overall Length 11 cm (4 1/2")

OP-826-02

Sinskey Lens Manipulating Hook
 Ø 0.25 mm Blunt Tip With 0.35 mm
 Height With Angle to Tip Length of
 10 mm
 Overall Length 11 cm (4 1/2")


Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Ophtalmologie Instruments
 Oftalmologia Strumenti



OP-826-03
Sinskey II Lens Manipulating Hook
 Ø 0.2 mm Blunt Tip With 0.35 mm
 Height
 Overall Length 11.5 cm (4 1/2")



OP-826-04
Sinskey II Lens Manipulating Hook
 Ø 0.2 mm Blunt Tip With 0.35 mm
 Height With Angle to Tip Length of
 10 mm
 Overall Length 11.5 cm (4 1/2")



OP-827-01
Straight
 0.76 mm x 0.75 mm
 x 0.25 mm Rotator Tip



OP-827-02
Angled
 0.76 mm x 0.75 mm
 x 0.25 mm Rotator Tip

Kuglen Iris Hook and Lens Manipulator
 Overall Length 11.5 cm (4 1/2")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Optalmologie Instruments
 Oftalmologia Strumenti

OP

OPHTHALMIC SURGERY INSTRUMENTS


OP-829-01

Kuglen Iris Hook and Lens Manipulator
 Ø 0.76 mm x 0.75 mm x 0.25 mm Clover
 Leaf Model Tip
 Overall Length 11.5 cm (4 ½")


OP-831-12

Rentsch Boat Hook
 Ø 0.2 mm Double Hooks
 Overall Length 12 cm (4 ¾")


OP-833-12

Osher Y Hook
 0.67 mm Wide Vertical Fork Shaped
 Tips With Angled to Tip Length of 9 mm
 Overall Length 12 cm (4 ¾")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Ophtalmologie Instruments
 Oftalmologia Strumenti



OP-835-01
Fenzl Lens Manipulating Hook
 Ø 0.15 mm Hook and
 0.5 mm Height
 Overall Length 12 cm (4 ¾")



OP-835-02
Fenzl Lens Manipulating Hook
 Ø 0.15 mm Hook and
 0.5 mm Height With Angle to
 Tip Length of 9 mm
 Overall Length 12 cm (4 ¾")



OP-837-01
Straight
 Ø 0.15 mm Dual Fork
 Shaped Tip



OP-837-02
Angled
 Ø 0.15 mm Dual Fork
 Shaped Tip With Angle
 to Tip Length of 9 mm

Maltzman-Fenzl Lens Manipulating Hook
 Overall Length 12 cm (4 ¾")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Optalmologie Instruments
 Oftalmologia Strumenti

OP

OPHTHALMIC SURGERY INSTRUMENTS


OP-839-02

Jaffe-Maltzman Lens Manipulating Hook
 Angled, Ø 0.15 mm Dual Fork
 Shaped Tip With Angle to Tip Length of
 11 mm
 Overall Length 12 cm (4 ¾")


OP-841-01

Giannetti Push Pull Manipulator
 Straight
 0.76 mm x 0.75 mm x 0.25 mm Tip
 Overall Length 12.5 cm (5")


OP-843-02

Deitz ICL Manipulator
 0.8 mm Wide Vertical Fork
 Shaped Tips, 11 mm Long
 From Bend to Tip
 Overall Length 12 cm (4 ¾")

*ICL® is a registered trademark of
 STAAR Surgical Company*

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Ophtalmologie Instruments
 Oftalmologia Strumenti



OP-845-01
Fixation Hook
 Curved 2 Sharp
 Pronged Hook With
 1.5 mm Spread
 Overall Length 12.5 cm (5")



OP-845-02
Fixation Hook
 Curved 2 Sharp
 Pronged Hook With
 2.0 mm Spread
 Overall Length 12.5 cm (5")



OP-849-01
 Right
 (Clockwise)



OP-849-02
 Left
 (Counter-Clockwise)




OP-851-09
Silverman Fixation Hook
 Blunt Cone Shaped Tips With
 3.0 mm Spread
 Overall Length 9.5 cm (3 3/4")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Optalmologie Instruments
 Oftalmologia Strumenti


OP
OPHTHALMIC SURGERY INSTRUMENTS

OP-861-11


Lu-Mendez Fixation Ring and Guide
 I.D. Ø 11 mm, LRI Guide Marked in 10°
 Increments and Fixation Ring With 15
 Blunt Teeth
 Overall Length 11 cm (4 1/4")


OP-863-11

Gimbel Stabilization Ring
 Offset C Shaped Swivel Fixation
 Ring With Ø 13 mm With
 10 mm Opening and 9 Blunt Teeth
 Overall Length 11.5 cm (4 1/2")


OP-865-13

C Shaped Swivel Fixation Ring
 With Ø 13 mm With
 10 mm Opening and 9 Blunt Teeth


OP-865-16

C Shaped Swivel Fixation Ring
 With Ø 16 mm With
 10 mm Opening and 12 Blunt Teeth

Fine-Thornton Swivel Fixation Ring

Overall Length 12 cm (4 3/4")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Ophtalmologie Instruments
 Oftalmologia Strumenti



OP-867-10

Thornton Globe Fixation Ring
 Ø 16 mm Globe Fixation Ring
 With 12 Blunt Teeth
 Overall Length 10 cm (4")



OP-869-11

Hofmann-Thornton Globe Fixation Ring
 Ø 16 mm Globe Fixation Ring With
 Swivel Handle And 12 Blunt Teeth
 Overall Length 11 cm (4 1/4")



OP-871-12

Stevens Tenotomy Hook
 Curved 6 mm Wide Hook
 Overall Length 12.5 cm (5")




OP-872-12

Manson Double Strabismus Hook
 Angled 3.0 mm Wide Fork Shaped
 Tips With 7.0 mm Spread
 Overall Length 12.5 cm (5")

Ophthalmology Instruments
Ophthalmologie Instrumente
Oftalmología Instrumentos
Ophtalmologie Instruments
Oftalmologia Strumenti

OP

OPHTHALMIC SURGERY INSTRUMENTS

**OP-873-12**

Helveston Teaser Hook
Angled 6 mm Wide Blunt Tip Hook
Overall Length 12.5 cm (5")

OP-875-13

Gass Retinal Detachment Hook
Flattened 13 mm Wide Graefe
Type Hook With 1.5 mm Oval Eye
13 cm (5")

Helveston Finder Hook

Overall Length 12.5 cm (5")

OP-877-01

Small 8 mm Wide Hook
With Extra Delicate Cone
Shaped Tip

OP-877-02

Large 12 mm Wide Hook With
Extra Delicate Cone Shaped Tip

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Ophtalmologie Instruments
 Oftalmologia Strumenti



Helveston Muscle Hook
 Blunt Spatula Shaped Tip
 Overall Length 13 cm (5")



- OP-879-01**
 8 mm Wide
OP-879-02
 10 mm Wide
OP-879-03
 12 mm Wide



OP-881-01
Jameson Muscle Hook
 Small, 6 mm Wide Hook
 With 2 mm Bulbous Tip
 Overall Length 13 cm (5")



OP-881-02
Jameson Muscle Hook
 Large, 8 mm Wide Hook
 With 2 mm Bulbous Tip
 Overall Length 13.5 cm (5 1/4")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Optalmologie Instruments
 Oftalmologia Strumenti

OP

OPHTHALMIC SURGERY INSTRUMENTS



OP-883-13
Green Strabismus Hook
 Flat 8 mm Wide Hook
 Overall Length 13.5 cm (5 1/4")




OP-885-13
Chavasse Strabismus Hook
 12 mm Wide W Shaped Hook
 Overall Length 13 cm (5")



OP-887-11
Chayet Corneal Marker
 For LASIK Surgery
 Overall Length 11.5 cm (4 1/2")



OP-891-11
Kellen Capsulorrhesis Marker
 Ø 6 mm Marking Ring
 With 2.5 mm Extensions
 Overall Length 11.5 cm (4 1/2")



OP-900-00
Maloney Intra-Operative Keratometer
 For Qualitative Measurement of Astigmatism
 Cone Shaped With Maximum Ø 31.0 mm
 and 21.5 cm Height, Titanium
 Overall Length 11.5 cm (4 1/2")


TITANIUM


Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Ophtalmologie Instruments
 Oftalmologia Strumenti


OP-901-11
Mendez Degree Gauge

Ø 12.0 mm I.D. and 18.0 mm O.D.
 Measures 0-180° With 10° Increments
 Overall Length 11.5 cm (4 ½")


OP-905-11
Grene Visual Axis Marker

Ø 0.5 mm With Concave Tip
 Overall Length 11.5 cm (4 ½")


OP-907-11
McDonald Optic Zone Marker

Double Ended, With Cross Hairs
 Ø 6.5 mm and 7.0 mm
 Overall Length 11 cm (4 ½")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Optalmologie Instruments
 Oftalmologici Strumenti

OP

OPHTHALMIC SURGERY INSTRUMENTS



OP-911-01 Ø 3.0 mm	OP-911-02 Ø 3.25 mm	OP-911-03 Ø 3.5 mm	OP-911-04 Ø 3.75 mm	OP-911-05 Ø 4.0 mm	OP-911-06 Ø 4.25 mm	OP-911-07 Ø 4.5 mm	OP-911-08 Ø 4.75 mm
OP-911-09 Ø 5.0 mm	OP-911-10 Ø 6.0 mm	OP-911-11 Ø 7.0 mm	OP-911-12 Ø 8.0 mm	OP-911-13 Ø 9.0 mm			

Bores Optic Zone Marker

Without Cross Hairs

Overall Length 10.5 cm (4 1/4")



OP-913-01 Ø 3.0 mm	OP-913-02 Ø 3.25 mm	OP-913-03 Ø 3.5 mm	OP-913-04 Ø 3.75 mm	OP-913-05 Ø 4.0 mm	OP-913-06 Ø 4.25 mm	OP-913-07 Ø 4.5 mm	OP-913-08 Ø 4.75 mm
OP-913-09 Ø 5.0 mm	OP-913-10 Ø 5.5 mm	OP-913-11 Ø 6.0 mm	OP-913-12 Ø 7.0 mm	OP-913-13 Ø 8.0 mm	OP-914-14 Ø 9.0 mm		

Hoffer Optic Zone Marker

With Cross Hairs

Overall Length 10.5 cm (4 1/4")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Ophtalmologie Instruments
 Oftalmologia Strumenti


OP-915-08

With 8 Radial Blades

OP-915-12

With 12 Radial Blades

OP-915-16

With 16 Radial Blades

Anis Radial Marker


Overall Length 10 cm (4")


OP-917-10
Anis Corneal Suture Placement Marker
 Ø 8 mm Center Opening and 8 Radial Blades
 Overall Length 10 cm (4")


OP-919-09
Fechtner Trabeculectomy Marker
 With U-Shaped 2.75mm Marking
 Pattern and Side Landmarks
 Overall Length 9 cm (3 1/2")

OP-930-12
Helveston Scleral Ruler
 Vaulted 15 mm Long Blade
 With Notches in 5 mm Increments
 Overall Length 12 cm (4 3/4")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Optalmologie Instruments
 Oftalmologia Strumenti

OP
OPHTHALMIC SURGERY INSTRUMENTS


OP-932-01
Stahl Corneal Caliper
 Measures 10.0 mm to 13.0 mm
 With 0.5 mm Increments



OP-932-02
Stahl Corneal Caliper
 Measures 10.25 mm to 13.25 mm
 With 0.5 mm Increments



OP-934-30
Braunstein Caliper
 Double Ended Fixed Caliper
 For 3.0 and 3.5 mm Marking



OP-934-35
Braunstein Caliper
 Double Ended Fixed Caliper
 For 3.5 and 4.0 mm Marking



OP-935-01
Schioetz Tonometer
 Inclined Scale
 With 5.5, 7.5 & 10 Grams Weights



OP-941-01
Harms Trabeculotomy Probe
 Right, Two 10 mm Long Tips
 With 3 mm Spread
 Overall Length 5 cm (2")



OP-941-02
Harms Trabeculotomy Probe
 Left, Two 10 mm Long Tips
 With 3 mm Spread
 Overall Length 5 cm (2")

Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Ophtalmologie Instruments
 Oftalmologia Strumenti




OP-944-08
Infant Lacrimal Dilator
 \varnothing 0.4 mm Very Fine 21 mm Long Blunt Tip
 Overall Length 8 cm (3")

OP-948-01
Wilder Lacrimal Dilator
 Fig. 1, Fine, Tapered
 \varnothing 0.32 mm / 32 mm Long Point
 Overall Length 10.5 cm (4 1/4")



OP-946-10
Muldoon Lacrimal Dilator
 \varnothing 0.35 mm Cylindrical Tip, With 27 mm Long Taper
 Overall Length 10 cm (4")



OP-948-02
Wilder Lacrimal Dilator
 Fig. 2, Medium, Tapered
 \varnothing 0.75 mm / 23 mm Long Point
 Overall Length 10.5 cm (4 1/4")

OP-948-03
Wilder Lacrimal Dilator
 Fig. 3, Long, Tapered
 \varnothing 1.1 mm / 19 mm Long Point
 Overall Length 10.5 cm (4 1/4")



Ophthalmology Instruments
 Ophthalmologie Instrumente
 Oftalmología Instrumentos
 Optalmologie Instruments
 Oftalmologia Strumenti

OP

OP-950-14
Castroviejo Lacrimal Dilator

Double Ended, One Fine Tapered (\varnothing 0.32 mm / 32 mm Long) and
 One Medium Tapered (\varnothing 0.75 / 23 mm Long) Points


OP-953-14
Worst Pigtail Probe

Double Ended, 8.0 mm Curved Blunt Probe With Eye Holes
 Overall Length 16 cm (6 1/4")


OP-956-14
Crawford Retrieval Hook

Designed to Pick The Crawford Probe
 Overall Length 14.5 cm (5 3/4")


OP-960-00
Crawford Lacrimal Intubation Set

Including Silicone Tube and Two Olive Tip Probes \varnothing 0.6 mm
 Probes Length 11 cm (4 1/4") Tube Length 30.5 cm (12")

Scalpel Handles, Scalpel Blades
 Skalpellgriffe, Skalpellklingen
 Mangos de bisturíes, Hojas de bisturíes
 Manches pour bistouris, Lames pour bistouris
 Manici per bisturi, Lame per bisturi



SL-010-03
Bard-Parker Scalpel Handle # 3
 Solid
 Overall Length 12 cm (4 3/4")



Scalpel Handle
 For Micro Scalpel Blades

SL-012-03
Bard-Parker Scalpel Handle # 3
 Solid, CM Graduation
 Overall Length 12 cm (4 3/4")




SL-100-11
Scalpel Blade # 11
 Pack of 100
 Fig. 11



SL-100-15
Scalpel Blade # 15
 Pack of 100
 Fig. 15




- SL-206-10**
 10 cm (4")
- SL-206-13**
 13 cm (5")
- SL-206-15**
 15 cm (6")



SL-210-61
Scalpel Handle # 3
 Fig. 61
 Pack of 25



SL-210-62 Fig. 62 **SL-210-63** Fig. 63 **SL-210-64** Fig. 64 **SL-210-65** Fig. 65



SL-210-66 Fig. 66 **SL-210-67** Fig. 67 **SL-210-68** Fig. 68 **SL-210-69** Fig. 69

Blade Holders / Blade Breakers
 Klingenthaler / Kingenbrecher
 Porta-hojas / Rompe-hojas de afeitar
 Pincettes pour saisir et casser les lames de rasoir
 Pinze per afferrare e rompere le lame di rasoio

OP

OPHTHALMIC SURGERY INSTRUMENTS



SL-900-09
Straight



SL-901-09
Angled



SL-910-11
Barraquer Blade Holder/Breaker
Overall Length 11 cm (4 1/4")



SL-916-13
Castroviejo Blade Holder/Breaker
Round Handle
Overall Length 13 cm (5")

Troutman-Chris Blade Holder/Breaker
Round Handle
Overall Length 9 cm (3 1/2")

Dissecting Scissors, Gum Scissors, Iris Scissors
 Präparierscheren, Zahnfleischscheren, Irisscheren
 Tijeras de disección, Tijeras para encías, Tijeras para iridectomía
 Ciseaux à dissection, Ciseaux à gencives, Ciseaux à iridectomie
 Forbici per dissezione, Forbici per gengivi, Forbici per iridectomia



Ragnell (Kilner) Dissecting Scissor

SS-243-12
 12 cm (4 3/4")

SS-243-15
 15 cm (6")




Goldman-Fox Gum Scissor
 Overall Length 13 cm (5")



SS-400-13
 Straight


SUPERedge™



SS-400-13 SC
 Straight



SS-401-13
 Curved



SS-401-13 SC
 Curved



Iris Scissor

SS-430-09
 Straight
 9 cm (3 1/2")

SS-430-10
 Straight
 10 cm (4")

SS-431-09
 Curved
 9 cm (3 1/2")

SS-431-10
 Curved
 10 cm (4")

Iris Scissors
 Irischeren
 Tijeras para iridectomía
 Ciseaux à iridectomie
 Forbici per iridectomia

OP

OPHTHALMIC SURGERY INSTRUMENTS



Bonn Iris Scissor
 Overall Length 9 cm (3 ½")

SS-434-09
 Straight

SS-435-09
 Curved




Iris Scissor
 Overall Length 11.5 cm (4 ½")

SS-440-11
 Straight

SS-441-11
 Curved

SUPEREdge™



SS-440-11 SC
 Straight



SS-441-11 SC
 Curved



Dissecting Scissors
 Präparierscheren
 Tijeras de disección
 Ciseaux à dissection
 Forbici per dissezione

OP

OPHTHALMIC SURGERY INSTRUMENTS


SUPEREdge™


SS-452-10
 Straight
 Blunt/Blunt
 10 cm (4")

SS-452-11
 Straight
 Blunt/Blunt
 11 cm (4 ½")

SS-452-12
 Straight
 Blunt/Blunt
 12.5 cm (5")

SS-452-11 SC
 Straight
 Blunt/Blunt
 11 cm (4 ½")



SS-453-10
 Curved
 Blunt/Blunt
 10 cm (4")

SS-453-11
 Curved
 Blunt/Blunt
 11 cm (4 ½")

SS-453-11
 Curved
 Blunt/Blunt
 12.5 cm (5")

Stevens Tenotomy Scissor



Dissecting Scissors
Präparierscheren
Tijeras de disección
Ciseaux à dissection
Forbici per dissezione

OP

OPHTHALMIC SURGERY INSTRUMENTS

**SS-461-11****Littler Dissecting Scissor**

Curved, Tip With Eye For Suture
Overall Length 11.5 cm (4 1/2")


**SS-466-11**

Straight

Kilner Delicate Dissecting Scissor

Overall Length 11.5 cm (4 1/2")

Enucleation Scissors
Enukleationsschere
Tijeras de enucleación
Ciseaux d'énucléation
Forbici per enucleazione




Enucleation Scissor
Overall Length 11.5 cm (4 1/2")

Ribbon Ophthalmic Scissors
 Band-Augenschere
 Tijera oftálmica de cinta
 Ciseaux ophtalmiques à ruban
 Forbici oftalmica a nastro

OP

OPHTHALMIC SURGERY INSTRUMENTS



SS-484-09
Salyer Ribbon Scissor
 Straight, Flat Shanks
 Overall Length 9.5 cm (3 ¾")



Iris Scissor (Ribbon Type)
 Overall Length 9.5 cm (3 ¾")
SS-486-09
 Straight, Sharp, Flat Shanks



Iris Scissor (Ribbon Type)
 Overall Length 10 cm (4")
SS-488-10
 Straight, Blunt, Flat Shanks

SS-487-09
 Curved, Sharp, Flat Shanks

SS-489-10
 Curved, Blunt, Flat Shanks




Strabismus Scissor (Ribbon Type)
 Overall Length 10 cm (4")

SS-492-10
 Straight, Flat Shanks

SS-493-10
 Curved, Flat Shanks

Ribbon Ophthalmic Scissors
 Band-Augenschere
 Tijera oftálmica de cinta
 Ciseaux ophtalmiques à ruban
 Forbici oftalmica a nastro



Stevens Tenotomy Scissor (Ribbon Type)
 Overall Length 10 cm (4")

SS-496-10
 Straight, Flat Shanks, Sharp Tips

SS-497-10
 Curved, Flat Shanks, Sharp Tips

SS-498-10
 Straight, Flat Shanks, Blunt Tips

SS-499-10
 Curved, Flat Shanks, Blunt Tips



Gradle Scissor (Ribbon Type)
 Overall Length 10 cm (4")

SS-500-10
 Straight, Flat Shanks

SS-501-10
 Curved, Flat Shanks




SS-507-11
Kaye Facelift Scissor (Ribbon Type)

Curved, Flat Shanks, Toothed
 Overall Length 11 cm (4 ½")


Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

OP

OPHTHALMIC SURGERY INSTRUMENTS

**SS-626-01****Williamson-Noble Micro Scissor**

Straight, Sharp/Sharp 27 mm Long Blades With
Max. 14 mm Cutting Edge
Overall Length 11.5 cm (4 ¾")

**SS-630-12**
Straight**SS-631-12**
Curved**Noyes Iris Scissor**


Sharp/Sharp, 23 mm Long Blades With
Max. 14 mm Cutting Edge
Overall Length 12.5 cm (5")



Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici




SS-632-12
Straight



SS-633-12
Curved



SS-634-12
Straight



SS-635-12
Curved

Noyes Iris Scissor

Sharp/Blunt
23 mm Long Blades With
Max. 14 mm Cutting Edge
Overall Length 12.5 cm (5")

Noyes Iris Scissor

Blunt/Blunt
23 mm Long Blades With
Max. 14 mm Cutting Edge
Overall Length 12.5 cm (5")

Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

OP
OPHTHALMIC SURGERY INSTRUMENTS


Micro Scissor
13 mm Long Blades
With Max. 10 mm Cutting Edge
Overall Length 12 cm (4 $\frac{3}{4}$ ")

SS-638-10
Osher IOL Cutting Scissor
Straight
Micro Serrated Cutting Edge With
Blunt Tip 10 mm Long Blades With
Max. 5 mm Cutting Edge
Overall Length 10 cm (4")

SS-639-10
Osher IOL Cutting Scissor
Angled
Micro Serrated Cutting Edge With
Blunt Tip 10 mm Long Blades With
Max. 5 mm Cutting Edge
Overall Length 10 cm (4")

Micro Scissors
 Mikroscheren
 Micro-Tijeras
 Micro-Ciseaux
 Micro-Forbici



Wecker Iris Scissor

8 mm Long Blades With Max. 5 mm Cutting Edge
 Overall Length 11 cm (4 1/2")

Micro Scissors

Mikroscheren

Micro-Tijeras


Micro-Ciseaux

Micro-Forbici




SS-641-11

Angled, Sharp/Sharp



SS-643-11

Angled, Sharp/Blunt



SS-645-11

Angled, Blunt/Blunt



Barraquer Iris Scissor

9 mm Long Blades Max. 6 mm Cutting Edge
 Overall Length 5.5 cm (2 1/4")

SS-647-05

Angled, 5 mm

SS-647-07

Angled, 7 mm

SS-647-10

Angled, 10 mm

Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

OP
OPHTHALMIC SURGERY INSTRUMENTS


SS-649-01
Curved Right

SS-649-02
Curved Left



Castroviejo Corneoscleral Scissor
8 mm Long Blades With
Max. 4 mm Cutting Edge
Overall Length 10.5 cm (4")

SS-651-01
Curved Right

SS-651-02
Curved Left



Castroviejo Corneoscleral Scissor
12 mm Long Blades With
Max. 7 mm Cutting Edge
Overall Length 11 cm (4")

SS-653-01
Curved Right

SS-653-02
Curved Left

Micro Scissors
 Mikroscheren
 Micro-Tijeras
 Micro-Ciseaux
 Micro-Forbici



Troutman-Castroviejo Corneoscleral Scissor
 With Lock
 8 mm Long Blades With
 Max. 4 mm Cutting Edge
 Overall Length 10.5 cm (4")

SS-655-01

Curved Right

SS-655-02

Curved Left



Troutman-Castroviejo Corneoscleral Scissor
 With Lock
 12 mm Long Blades With
 Max. 7 mm Cutting Edge
 Overall Length 11 cm (4 1/2")

SS-657-01

Curved Right

SS-657-02

Curved Left

Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

OP

**Troutman-Katzin Corneal Transplant Scissor**

Strongly Curved, With Lock
8 mm Long Blades With
Max. 4 mm Cutting Edge
Overall Length 10.5 cm (4")

SS-659-01

Curved Right

SS-659-02

Curved Left

**Katzin Corneal Transplant Scissor**

Strongly Curved
13 mm Long Blades With
Max. 7 mm Cutting Edge
Overall Length 11 cm (4 1/2")


SS-661-01

Curved Right

SS-661-02

Curved Left

Micro Scissors
 Mikroscheren
 Micro-Tijeras
 Micro-Ciseaux
 Micro-Forbici



Castroviejo Universal Corneal Scissor

Curved, Blunt Tips
 Overall Length 11 cm (4 1/2")

SS-663-11


10 mm Long Blades With
 Max. 5 mm Cutting Edge

SS-665-11

12 mm Long Blades With
 Max. 6 mm Cutting Edge

SS-667-11

14 mm Long Blades With
 Max. 8 mm Cutting Edge



SS-669-11

Castroviejo Universal Corneal Scissor

Strongly Curved, Blunt Tips 15 mm Long Blades With
 Max. 9 mm Cutting Edge
 Overall Length 11.5 cm (4 1/2")

Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

OP

OPHTHALMIC SURGERY INSTRUMENTS

**SS-671-11**

Castroviejo Keratoplasty Scissor
Angled to Side, Blunt Tips
14 mm Long Blades With
Max. 9 mm Cutting Edge
Overall Length 11.5 cm (4 1/2")

SS-673-11

McPherson-Westcott Conjunctival Scissor
Curved, Blunt Tips
10 mm Long Blades With
Max. 5 mm Cutting Edge
Overall Length 11 cm (4 1/2")

SS-675-11

Westcott Tenotomy Scissor
Curved, Blunt Tips
14 mm Long Blades With
Max. 8 mm Cutting Edge
Overall Length 11 cm (4 1/2")

Micro Scissors
 Mikroscheren
 Micro-Tijeras
 Micro-Ciseaux
 Micro-Forbici



Westcott Tenotomy Scissor

Blunt Tips
 20 mm Long Blades With
 Max. 10 mm Cutting Edge
 Overall Length 11 cm (4 1/2")

SS-677-01

Curved Right

SS-677-02

Curved Left



SS-679-11

Tibolt Punctal Scissor
 One Conically Shaped Blade With
 Micro Grooves & Textured Surface
 15 mm Long Blades With
 Max. 8 mm Cutting Edge
 Overall Length 11.5 cm (4 1/2")

SS-681-11

McPherson-Westcott Stitch Scissor
 Curved, Very Sharp Pointed Tips, Small Blades
 10 mm Long Blades With
 Max. 5 mm Cutting Edge
 Overall Length 11 cm (4 1/2")



Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

OP

OPHTHALMIC SURGERY INSTRUMENTS

**SS-683-11**

Westcott Stitch Scissor
Curved, Very Sharp Pointed Tips
20 mm Long Blades With
Max. 10 mm Cutting Edge
Overall Length 11.5 cm (4 1/2")

SS-685-11

Jaffe Stitch Scissor
Curved, Very Sharp Pointed Tips
13 mm Long Blades With
Max. 7 mm Cutting Edge
Overall Length 11 cm (4 1/2")

SS-688-11

Iris Micro Scissor
Straight, Sharp Tips
10 mm Long Blades With
Max. 5 mm Cutting Edge
Overall Length 11 cm (4 1/2")

SS-689-11

Iris Micro Scissor
Curved, Sharp Tips
10 mm Long Blades With
Max. 5 mm Cutting Edge
Overall Length 11 cm (4 1/2")

Micro Scissors
 Mikroscheren
 Micro-Tijeras
 Micro-Ciseaux
 Micro-Forbici


SS-690-11

Fine Micro Scissor
 Straight, Sharp Tips
 17 mm Long Blades With
 Max. 9 mm Cutting Edge
 Overall Length 11.5 cm (4 1/2")


SS-691-11

Fine Micro Scissor
 Curved, Sharp Tips
 17 mm Long Blades With
 Max. 9 mm Cutting Edge
 Overall Length 11.5 cm (4 1/2")


SS-694-08

Vannas Capsulotomy Scissor
 Straight, Sharp Tips
 5 mm Long Blades With
 Max. 3 mm Cutting Edge
 Overall Length 8 cm (3 1/4")


SS-695-08

Vannas Capsulotomy Scissor
 Curved, Sharp Tips
 5 mm Long Blades With
 Max. 3 mm Cutting Edge
 Overall Length 8 cm (3 1/4")

Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

OP
OPHTHALMIC SURGERY INSTRUMENTS

SS-696-08
Vannas Capsulotomy Scissor

Angled Forward, Sharp Tips 5 mm Long Blades With
Max. 3 mm Cutting Edge
Overall Length 8 cm (3 1/4")


SS-697-08
Vannas Capsulotomy Scissor

Angled to Side, Sharp Tips 5 mm Long Blades With
Max. 3 mm Cutting Edge
Overall Length 8 cm (3 1/4")



SS-698-08
Gills-Vannas Capsulotomy Scissor

Straight, Sharp Tips 7 mm Long Blades With
Max. 3 mm Cutting Edge
Overall Length 8 cm (3 1/4")


SS-699-08
Gills-Vannas Capsulotomy Scissor

Curved, Sharp Tips 7 mm Long Blades With
Max. 3 mm Cutting Edge
Overall Length 8 cm (3 1/4")

Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici



SS-701-08
Gills-Vannas Capsulotomy Scissor
Angled Forward, Sharp Tips 7 mm
Long Blades With Max. 3 mm
Cutting Edge
Overall Length 8 cm (3 1/4")



SS-703-08
McPherson-Vannas Capsulotomy Scissor
Curved, Sharp Tips 8 mm
Long Blades With
Max. 3 mm Cutting Edge
Overall Length 8 cm (3 1/4")




SS-705-08
Holland DALK Scissor
Angled Blades With Blunt Tips
5 mm Long Blades With
Max. 3 mm Cutting Edge
Overall Length 8 cm (3 1/4")

Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

OP

OPHTHALMIC SURGERY INSTRUMENTS

**SS-706-10**

Gills-Vannas Micro Scissor
Straight, Extra Thin Blades
Sharp Tips
7 mm Long Blades With
Max. 3 mm Cutting Edge
Overall Length 10.5 cm (4")

SS-707-10

Gills-Vannas Micro Scissor
Curved, Extra Thin Blades
Sharp Tips
7 mm Long Blades With
Max. 3 mm Cutting Edge
Overall Length 10.5 cm (4")

SS-709-10

Gills-Vannas Micro Scissor
Angled Forward, Extra Thin Blades
Sharp Tips
7 mm Long Blades With
Max. 3 mm Cutting Edge
Overall Length 10.5 cm (4")

Micro Scissors
 Mikroscheren
 Micro-Tijeras
 Micro-Ciseaux
 Micro-Forbici



SS-710-11
Stern-Gills Micro Scissor
 Straight, Extra Thin Blades
 Sharp Tips
 10 mm Long Blades With
 Max. 4 mm Cutting Edge
 Overall Length 11 cm (4 1/2")



SS-711-11
Stern-Gills Micro Scissor
 Curved, Extra Thin Blades
 Sharp Tips
 10 mm Long Blades With
 Max. 4 mm Cutting Edge
 Overall Length 11 cm (4 1/2")




SS-713-11
Stern-Gills Micro Scissor
 Angled Forward, Extra Thin Blades
 Sharp Tips
 10 mm Long Blades With
 Max. 4 mm Cutting Edge
 Overall Length 11 cm (4 1/2")



SS-715-11
Gills Capsulotomy Scissor
 Angled Blades With Curved Sharp Tips
 10 mm Long Blades With
 Max. 4 mm Cutting Edge
 Overall Length 10.5 cm (4 1/2")

Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

OP
OPHTHALMIC SURGERY INSTRUMENTS


SS-720-20
Straight, Horizontal Opening 20 G



SS-721-20
Curved, Horizontal Opening 20 G
SS-721-23
Curved, Horizontal Opening 23 G



SS-721-30
Angled 30°, Horizontal Opening
20 G



SS-721-45
Angled 45°, Horizontal Opening
20 G



SS-723-23
Vertical Straight, 23 G



SS-723-60
Angled 60°, Horizontal Opening, 20 G



SS-723-80
Angled 80°, Horizontal Opening, 20 G

Vitrectomy Scissor
Overall Length 14 cm (5 1/2")

Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

**SS-725-23****Alio MICS Scissor**

Curved, 23 G

Shaft With Horizontal Micro Blades for Cutting Delicate Membranes

Overall Length 14 cm (5 1/2")

Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

OP
OPHTHALMIC SURGERY INSTRUMENTS

SS-726-12
Vannas Micro Scissor

Straight, 10 mm Long Blades With
Max. 4 mm Cutting Edge
Overall Length 12 cm (4 ¾")


SS-727-12
Vannas Micro Scissor

Curved, 10 mm Long Blades With
Max. 4 mm Cutting Edge
Overall Length 12 cm (4 ¾")


SS-729-12
Vannas Micro Scissor

Angled on Flat, 10 mm Long Blades With
Max. 4 mm Cutting Edge
Overall Length 12 cm (4 ¾")



UltraCut™ TC Stevens Tenotomy Scissor
Overall Length 10.5 cm (4 ¼")


SS-972-10

Straight
Blunt/Blunt


SS-973-10

Curved
Blunt/Blunt


SS-974-10

Straight
Sharp/Sharp


SS-975-10

Straight
Sharp/Sharp



Dissecting Forceps
Anatomische Pinzetten
Pinzas de disección
Pince à dissection
Pinze per dissezione



Adson Dressing Forcep
Serrated Tip

FR-050-12
12 cm (4 3/4")

FR-050-15
15 cm (6")

FR-050-18
18 cm (7")



Mini-Adson Dressing Forcep
Overall Length 12 cm (4 3/4")

FR-052-12
Cross Serrated Tip

Dissecting Forceps
Anatomische Pinzetten
Pinzas de disección
Pince à dissection
Pinze per dissezione

OP

OPHTHALMIC SURGERY INSTRUMENTS



FR-054-12
12 cm (4 3/4")

FR-054-15
15 cm (6")



FR-250-12
12 cm (4 3/4")

FR-250-15
15 cm (6")

FR-250-18
18 cm (7")



Micro-Adson Dressing Forcep
Cross Serrated Tip

Adson Dissecting Forcep
1 x 2 Teeth




Dissecting Forceps
Anatomische Pinzetten
Pinzas de disección
Pince à dissection
Pinze per dissezione



FR-252-12
Straight



FR-253-12
Curved



FR-256-12
12 cm (4 ¾")



FR-256-15
15 cm (6")

Mini-Adson Dissecting Forcep
1 x 2 Teeth, Cross Serrated
Overall Length 12 cm (4 ¾")

Micro-Adson Dissecting Forcep
1 x 2 Teeth

Cilia Forceps
Zillenzange
Pinzas para cilios
Pince à cils
Pinze Cilia

OP

OPHTHALMIC SURGERY INSTRUMENTS



Beer Cilia Forcep
Overall Length 9 cm (3 1/2")



FR-630-09
Smooth Jaws




Barraquer Cilia Forcep
Overall Length 10 cm (4")



FR-632-10
Smooth Jaws



Micro Forceps
Mikropinzen
Micro Pinzas
Micro-Pincers
Micro-Pinze



Vitrectomy Forcep
14 cm (5 1/2")



FR-680-20

Straight, Serrated Tip 20 G

FR-680-23

Straight, Serrated Tip 23 G



FR-682-20

Straight, End Grasping Jaws
Horizontal Opening Tip 20 G

FR-682-23

Straight, End Grasping Jaws
Horizontal Opening Tip 23 G



FR-685-20

Smooth Side, Gripping Jaws
Vertical Opening Tip 20 G



FR-687-20

Membrane Peeling, Smooth
45° Angled Jaws, 20 G



FR-689-23

Asymmetrical Jaws 23 G



FR-690-23

ILM 23 G

Micro Forceps
Mikropinzen
Micro Pinzas
Micro-Pincers
Micro-Pinze

OP

OPHTHALMIC SURGERY INSTRUMENTS

**FR-692-20**

ICL Lens Cartridge Loading Forcep
Straight, 20 G
Jaws Feature a Series of Concave
and Convex Surfaces

ICL® is a registered trademark of
STAAR Surgical Company

FR-695-23

Alio MICS Capsulorrhesis Forcep
23 G, Curved Shaft With Short Micro-Grasping
Tip Designed For Use Through a Very Small Paracentesis.



Eye Surgery Forceps
 Augenchirurgie Pinzetten
 Ojo Pinzas Cirugía
 Pince de chirurgie oculaire
 Eye Pinze Chirurgia

OP

OPHTHALMIC SURGERY INSTRUMENTS




FR-710-07
Botvin Iris Forcep
 0.5 mm, Serrated
 Overall Length 7.5 cm (3 1/4")



FR-712-07
Botvin Iris Forcep
 0.3 mm, 1 x 2 Teeth
 Overall Length 7.5 cm (3 1/4")



FR-714-07
Gills-Colibri Corneal Forcep
 Very Fine Pointed Tip, With Tying Platform
 Overall Length 7.5 cm (3 1/4")



FR-716-07
Colibri Corneal Forcep
 Very Delicate, 1 x 2 Teeth, 0.12 mm
 Overall Length 7.5 cm (3 1/4")



FR-718-07
Harms-Colibri Corneal Forcep
 Very Delicate, 1 x 2 Oblique Teeth
 0.12 mm With Tying Platform
 Overall Length 7.5 cm (3 1/4")



FR-718-07 T
Harms-Colibri Corneal Forcep
 Very Delicate, 1 x 2 Oblique Teeth
 0.12 mm With Tying Platform
 Overall Length 7.5 cm (3 1/4")


TITANIUM

Eye Surgery Forceps
 Augenchirurgie Pinzetten
 Ojo Pinzas Cirugía
 Pince de chirurgie oculaire
 Eye Pinze Chirurgia



FR-720-07
Colibri Corneal Forcep
 Standard 1 x 2 Teeth
 With Tying Platform, 0.4 mm
 Overall Length 7.5 cm (3 1/4")



FR-721-07
Pierce-Colibri Corneal Forcep
 1 x 1 Teeth, 0.12 mm
 Overall Length 7.5 cm (3 1/4")



FR-722-08
Hoskin Tissue Forcep
 Straight, 1 x 1 Teeth, 0.3 mm
 Overall Length 8 cm (3 1/4")




FR-723-08
Hoskin Tissue Forcep
 Curved, 1 x 1 Teeth, 0.3 mm
 Overall Length 8 cm (3 1/4")

Eye Surgery Forceps
 Augenchirurgie Pinzetten
 Ojo Pinzas Cirugía
 Pince de chirurgie oculaire
 Eye Pinze Chirurgia

OP

OPHTHALMIC SURGERY INSTRUMENTS


FR-724-10

Fechtner Conjunctival Forcep
 Very Delicate Ø 0.8 mm
 Ring Shaped Tips
 With 6 mm Long Tying Platform
 Overall Length 10.5 cm (4")

FR-726-11

Castroviejo-Colibri Corneal Forcep
 Very Delicate 1 x 2 Teeth, 0.12 mm
 Overall Length 11 cm (4 1/4")

FR-728-11

Harms-Colibri Corneal Forcep
 Very Delicate 1 x 2 Teeth, 0.12 mm
 With Tying Platform
 Overall Length 11 cm (4 1/4")



Eye Surgery Forceps
Augenchirurgie Pinzetten
Ojo Pinzas Cirugía
Pince de chirurgie oculaire
Eye Pinze Chirurgia



FR-730-08

Bonn Corneal Forcep

Very Delicate 1 x 2 Teeth 0.12 mm
Overall Length 8.5 cm (4")



FR-734-08

Bonn Corneal Forcep

Very Delicate, Serrated
Overall Length 8.5 cm (3 1/4")



FR-732-10

Bonn Corneal Forcep

Very Delicate 1 x 2 Teeth
With Tying Platform, 0.12 mm
Overall Length 10 cm (4")



1:2



FR-732-10 T

Bonn Corneal Forcep

Very Delicate 1 x 2 Teeth
With Tying Platform, 0.12 mm, Titanium
Overall Length 10 cm (4")




1:2



Eye Surgery Forceps
 Augenchirurgie Pinzetten
 Ojo Pinzas Cirugía
 Pince de chirurgie oculaire
 Eye Pinze Chirurgia

OP

OPHTHALMIC SURGERY INSTRUMENTS


FR-736-10
McPherson Suture Tying Forcep

Straight, 1 x 2 Teeth
 With Tying Platform, 0.4 mm
 Overall Length 10.5 cm (4")

FR-737-10
McPherson Suture Tying Forcep

Angled, 1 x 2 Teeth
 With Tying Platform, 0.4 mm
 Overall Length 10 cm (4")

FR-740-10
Kelman-McPherson Tissue Forcep

Angled, 1 x 2 Teeth
 With Tying Platform, 0.3 mm
 Overall Length 10 cm (4")

Eye Surgery Forceps
 Augenchirurgie Pinzetten
 Ojo Pinzas Cirugía
 Pince de chirurgie oculaire
 Eye Pinze Chirurgia



FR-742-10
Hoffer-McPherson Tissue Forcep
 Very Delicate, 1 x 2 Teeth
 11 mm, 0.12 mm
 Overall Length 10 cm (4")



1:2



Castroviejo Suture Tying Forcep
 1 x 2 Oblique Teeth With Tying Platform
 Overall Length 11 cm (4 1/4")

FR-750-01 0.12 mm

FR-750-03 0.3 mm

FR-750-05 0.5 mm

FR-750-09 0.9 mm




1:2



FR-750-01 T
Castroviejo Suture Tying Forcep
 1 x 2 Oblique Teeth, With Tying Platform
 0.12 mm, Titanium
 Overall Length 11 cm (4 1/4")


TITANIUM

Eye Surgery Forceps
 Augenchirurgie Pinzetten
 Ojo Pinzas Cirugía
 Pince de chirurgie oculaire
 Eye Pinze Chirurgia

OP
OPHTHALMIC SURGERY INSTRUMENTS

Castroviejo Suture Tying Forcep

Angled, 1 x 2 Oblique Teeth
 With Tying Platform
 Overall Length 11 cm (4 1/4")

FR-751-01 0.12 mm

FR-751-03 0.3 mm

FR-751-05 0.5 mm

FR-751-09 0.9 mm

FR-751-01 T
Castroviejo Suture Tying Forcep

Angled, 1 x 2 Oblique Teeth With Tying Platform
 0.12 mm, Titanium
 Overall Length 11 cm (4 1/4")

FR-752-05
Stern-Castroviejo Suture Tying Forcep

1 x 2 Teeth, 0.5 mm
 Overall Length 11 cm (4 1/4")





Eye Surgery Forceps
 Augenchirurgie Pinzetten
 Ojo Pinzas Cirugía
 Pince de chirurgie oculaire
 Eye Pinze Chirurgia

OP

OPHTHALMIC SURGERY INSTRUMENTS


FR-762-09

Bishop-Harmon Dressing Forcep
 0.8 mm, Cross Serrated Jaws
 Overall Length 8.5 cm (3 1/4")


FR-764-09

St. Martin's Suture Tying Forcep
 0.3 mm, 1 x 2 Oblique Teeth
 With 5.0 mm Tying Platform
 Overall Length 9 cm (3 1/2")


FR-765-09

Paufique Suture Tying Forcep
 0.5 mm, 1 x 2 Oblique Teeth
 With 5.0 mm Tying Platform
 Overall Length 9 cm (3 1/2")


FR-766-09

Sauer Fixation Forcep
 0.3 mm, 1 x 2 Oblique Teeth
 Overall Length 9 cm (3 1/2")

Eye Surgery Forceps
 Augen chirurgie Pinzetten
 Ojo Pinzas Cirugía
 Pince de chirurgie oculaire
 Eye Pinze Chirurgia



FR-768-09

Manhattan Eye & Ear Fixation Forcep
 1 x 2 Teeth, 0.3 mm, With Tying Platform
 Overall Length 9 cm (3 1/2")



FR-770-11

Elschnig Fixation Forcep
 0.9 mm, 1 x 2 Oblique Teeth
 Overall Length 11 cm (4 1/4")



FR-772-01

Graefe Fixation Forcep
 4.5 mm Wide Jaws
 Without Lock
 Overall Length 11 cm (4 1/4")



FR-772-02

Graefe Fixation Forcep
 4.5 mm Wide Jaws
 With Lock
 Overall Length 11 cm (4 1/4")

Eye Surgery Forceps
 Augenchirurgie Pinzetten
 Ojo Pinzas Cirugía
 Pince de chirurgie oculaire
 Eye Pinze Chirurgia

OP

OPHTHALMIC SURGERY INSTRUMENTS



FR-774-10
Green Fixation Forcep
 9.0 mm Wide Jaws
 Overall Length 10.5 cm (4")



FR-774-11
Troutman Superior Rectus Forcep
 Angled, 0.5 mm 1 x 2 Oblique Teeth
 Overall Length 11 cm (4 1/4")



FR-776-01
 0.6 mm
 1 x 2 Teeth



FR-776-02
 0.6 mm
 2 x 3 Teeth




Lester Fixation Forcep
 Overall Length 9.5 cm (3 3/4")

Eye Surgery Forceps
 Augenchirurgie Pinzetten
 Ojo Pinzas Cirugía
 Pince de chirurgie oculaire
 Eye Pinze Chirurgia




FR-777-01
Moody Fixation Forcep
 Curved Right
 Overall Length 11 cm (4 1/4")



FR-777-02
Moody Fixation Forcep
 Curved Left
 Overall Length 11 cm (4 1/4")



FR-780-10
McPherson Suture Tying Forcep
 Straight, With Tying Platform, 5.0 mm
 Overall Length 10 cm (4")



FR-780-10 T
McPherson Suture Tying Forcep
 Straight, With Tying Platform, 5.0 mm
 Titanium
 Overall Length 10 cm (4")



Eye Surgery Forceps
Augenchirurgie Pinzetten
Ojo Pinzas Cirugía
Pince de chirurgie oculaire
Eye Pinze Chirurgia

OP

OPHTHALMIC SURGERY INSTRUMENTS

**FR-781-10****McPherson Suture Tying Forcep**

Angled, With Tying Platform, 5 mm
Overall Length 10 cm (4")

**FR-781-10 T****McPherson Suture Tying Forcep**

Angled, With Tying Platform, 5 mm Titanium
Overall Length 10 cm (4")



Eye Surgery Forceps
 Augenchirurgie Pinzetten
 Ojo Pinzas Cirugía
 Pince de chirurgie oculaire
 Eye Pinze Chirurgia



FR-784-10
Kelman-McPherson Suture Tying Forcep
 Angled, With Tying Platform, 7.5 mm
 Overall Length 10 cm (4")



FR-784-10 T
Kelman-McPherson Suture Tying Forcep
 Angled, With Tying Platform, 7.5 mm Titanium
 Overall Length 10 cm (4")



FR-790-10
Kurstin Flap Stretching Forcep
 Disc Shaped, Serrated Jaws
 Overall Length 10.5 cm (4")



FR-793-10
Buratto LASIK Flap Forcep
 Curved, Disc Shaped, Serrated Jaws
 Overall Length 10.5 cm (4")

Eye Surgery Forceps
 Augenchirurgie Pinzetten
 Ojo Pinzas Cirugía
 Pince de chirurgie oculaire
 Eye Pinze Chirurgia

OP

OPHTHALMIC SURGERY INSTRUMENTS



FR-795-10
Mendez Multipurpose LASIK Forcep
 Angled, Spatulated Jaws & Tip
 Overall Length 10 cm (4")



FR-797-10
Kraft Capsule Tag Forcep
 Angled, 11 mm, 1.5 mm
 Overall Length 10.5 cm (4")



FR-801-10
Utrata Capsulorrhesis Forcep
 Very Delicate Triangular Grasping Tips Extremely Thin, 11 mm
 Over Length 10.5 cm (4")



FR-801-10 T
Utrata Capsulorrhesis Forcep
 Very Delicate Triangular Grasping Tips Extremely Thin, 11 mm, Titanium
 Overall Length 10.5 cm (4")


TITANIUM




Eye Surgery Forceps
Augenchirurgie Pinzetten
Ojo Pinzas Cirugía
Pince de chirurgie oculaire
Eye Pinze Chirurgia


FR-803-10

Jaffe Capsulorhexis Forcep
Very Delicate Round Grasping Tip
Extremely Thin, 11 mm
Overall Length 10.5 cm (4")


FR-805-10

O'Gawa-Utrata Capsulorhexis Forcep
Very Delicate, Round Grasping Tips
Extremely Thin
Overall Length 11.5 cm (4 1/2")


FR-807-11

Inamura Capsulorhexis Forcep
Delicate Tapered 12 mm Long Shanks
Overall Length 11.5 cm (4 1/2")


FR-809-10

Market Capsulorhexis Forcep
With Very Delicate Triangle Grasping Tip & Very
Delicate Vaulted Shanks, Length 11 mm
Overall Length 11.5 cm (4 1/2")

Eye Surgery Forceps
 Augenchirurgie Pinzetten
 Ojo Pinzas Cirugía
 Pince de chirurgie oculaire
 Eye Pinze Chirurgia

OP

OPHTHALMIC SURGERY INSTRUMENTS


FR-810-10

Jaffe Suture Tying Forcep
 Straight, Very Delicate
 Smooth Jaws, 6 mm
 Overall Length 10.5 cm (4")


FR-811-10

Jaffe Suture Tying Forcep
 Curved, Very Delicate
 Smooth Jaws 6 mm
 Overall Length 10.5 cm (4")


FR-816-10

Tennant Suture Tying Forcep
 Straight, Round Handle With
 Guide Pin Extra Delicate
 Smooth Jaws, 6 mm
 Overall Length 10.5 cm (4")


FR-816-10 T


Tennant Suture Tying Forcep
 Straight, Round Handle With
 Guide Pin Extra Delicate Smooth
 Jaws, 6 mm
 Overall Length 10.5 cm (4")



Eye Surgery Forceps
 Augenchirurgie Pinzetten
 Ojo Pinzas Cirugía
 Pince de chirurgie oculaire
 Eye Pinze Chirurgia



FR-817-10
Tennant Suture Tying Forcep
 Curved, Round Handle With Guide Pin
 Extra Delicate Smooth Jaws, 6 mm
 Overall Length 10.5 cm (4")



FR-817-10 T
Tennant Suture Tying Forcep
 Curved, Round Handle With Guide Pin
 Extra Delicate Smooth Jaws, 6 mm
 Titanium
 Overall Length 10.5 cm (4")



FR-820-10
Harms Suture Tying Forcep
 Straight
 Overall Length 10.5 cm (4")



FR-821-10
Harms Suture Tying Forcep
 Curved
 Overall Length 10.5 cm (4")



Eye Surgery Forceps
 Augenchirurgie Pinzetten
 Ojo Pinzas Cirugía
 Pince de chirurgie oculaire
 Eye Pinze Chirurgia

OP


OPHTHALMIC SURGERY INSTRUMENTS




FR-824-10
Helveston Suture Tying Forcep
 Straight
 Overall Length 10.5 cm (4")



FR-825-10
Fine Suture Tying Forcep
 Curved
 Overall Length 10.5 cm (4")



FR-827-08
Gerl Suture Tying Forcep
 Smooth Jaws, Curved Shanks, 6 mm
 Overall Length 8.5 cm (3 1/4")




FR-827-08 T
Gerl Suture Tying Forcep
 Smooth Jaws, Curved Shanks, 6 mm
 Titanium
 Overall Length 8.5 cm (3 1/4")



Eye Surgery Forceps
 Augen chirurgie Pinzetten
 Ojo Pinzas Cirugía
 Pince de chirurgie oculaire
 Eye Pinze Chirurgia




FR-830-10
Gradle Cilia Forcep
 Smooth Jaws
 Overall Length 10 cm (4")




Bonaccolto Conjunctival Forcep
 Overall Length 10 cm (4")

FR-834-12
 Longitudinal Serrated Jaws
 With Cross Serrated Tips, 1.2 mm

FR-834-17
 Longitudinal Serrated Jaws
 With Cross Serrated Tips, 1.7 mm



FR-836-10
Lister Conjunctival Forcep
 1 x 2 Teeth
 Overall Length 10 cm (4")



FR-838-10
Moorfield Conjunctival Forcep
 Overall Length 10 cm (4")

Eye Surgery Forceps
 Augenchirurgie Pinzetten
 Ojo Pinzas Cirugía
 Pince de chirurgie oculaire
 Eye Pinze Chirurgia

OP

OPHTHALMIC SURGERY INSTRUMENTS



Nugent Utility Forcep
 Overall Length 10 cm (4")

FR-840-10
 Smooth Jaws

FR-842-10
 Serrated Jaws

FR-844-11
Lister Utility Forcep
 Serrated Jaws, 1 x 2 Teeth
 Overall Length 11 cm (4 1/4")

FR-847-08
Scleral Plug Removing Forcep
 Overall Length 8.5 cm - (3 1/2")

FR-846-10
Gold Punctal Plug Forcep
 With Longitudinal Groove
 Pointed Tip
 Overall Length 11 cm (4 1/4")

Eye Surgery Forceps
 Augenchirurgie Pinzetten
 Ojo Pinzas Cirugía
 Pince de chirurgie oculaire
 Eye Pinze Chirurgia



FR-848-11
Osher Iris Forcep
 Very Fine Textured Jaws
 Overall Length 11 cm (4 1/4")



FR-850-11
Foreign Body Forcep
 Straight
 Overall Length 8.5 cm (3 1/4")



FR-855-11
Jaffe Utility Forcep
 Curved, Very Delicate
 Smooth Jaws
 Overall Length 11 cm (4 1/4")



FR-857-11
Castroviejo Capsule Forcep
 Overall Length 11 cm (4 1/4")



FR-859-11
Arruga Capsule Forcep
 Curved, Very Delicate
 Overall Length 11 cm (4 1/4")

Eye Surgery Forceps
 Augenchirurgie Pinzetten
 Ojo Pinzas Cirugía
 Pince de chirurgie oculaire
 Eye Pinze Chirurgia

OP

OPHTHALMIC SURGERY INSTRUMENTS


FR-861-11

Akahoshi Phaco PreChopper
 Delicate Jaws With Sharp Tip
 4.5 mm
 Overall Length 11 cm (4 1/4")

FR-863-11

Ernest Nucleus Cracker
 Extra Delicate Cross Action
 Serrated Paddle Tip
 Overall Length 11 cm (4 1/4")

FR-867-11

Dodick Nucleus Cracker
 Cross Action, Longitudinally
 Serrated Paddle Shaped Jaws
 Overall Length 11 cm (4 1/4")

FR-869-10

Alfonso Nucleus Grasping Forcep
 With Two Short Rows of Delicate
 Interlocking Teeth
 Overall Length 10.5 cm (4")

Eye Surgery Forceps
 Augenchirurgie Pinzetten
 Ojo Pinzas Cirugía
 Pince de chirurgie oculaire
 Eye Pinze Chirurgia



FR-871-11
Rosenwasser Donor Lamella Inserting Forcep
 With Vaulted Shanks
 Overall Length 11 cm (4 1/4")



FR-875-11
Dodick IOL Folding Forcep
 For Short IOLs
 Overall Length 11 cm (4 1/4")



FR-877-10
Livernois IOL Folding Forcep
 Concave Jaw Surfaces, For Soft IOLs
 Overall Length 11 cm (4 1/4")

Eye Surgery Forceps
Augenchirurgie Pinzetten
Ojo Pinzas Cirugía
Pince de chirurgie oculaire
Eye Pinze Chirurgia

OP

OPHTHALMIC SURGERY INSTRUMENTS

**FR-879-10**

Steinert IOL Inserting Forcep
Conically Shaped Jaws
Overall Length 10 cm (4")

**FR-882-10**

Livernois IOL Holding Forcep
Straight, Highly Polished Round
Jaws For Holding Soft IOLs
Overall Length 10.5 cm (4")

**FR-883-10**

Livernois IOL Holding Forcep
Angled, Highly Polished Round
Jaws For Holding Soft IOLs
Overall Length 10.5 cm (4")

Eye Surgery Forceps
 Augenchirurgie Pinzetten
 Ojo Pinzas Cirugía
 Pince de chirurgie oculaire
 Eye Pinze Chirurgia



FR-885-10
Livernois-McDonald IOL Inserting Forcep
 Concave Jaws, 5.0 mm
 Overall Length 10.5 cm (4")



FR-887-10
Faulkner Lens Holding Forcep
 Round, Smooth Jaws
 Overall Length 10 cm (4")



FR-889-10
Buratto IOL Inserting Forcep
 Extra Delicate Highly Polished Convex Jaws
 For Soft IOLs
 Overall Length 10 cm (4")

Eye Surgery Forceps
 Augenchirurgie Pinzetten
 Ojo Pinzas Cirugía
 Pince de chirurgie oculaire
 Eye Pinze Chirurgia

OP

OPHTHALMIC SURGERY INSTRUMENTS


FR-891-12

Bechert Lens Holding Forcep
 Very Delicate Flat Jaws
 Overall Length 12.5 cm (5")

FR-893-12

Kratz Lens Holding Forcep
 Very Delicate Narrow Jaws
 Overall Length 12.5 cm (5")

FR-895-11


Blaydes Lens Holding Forcep
 Very Delicate Narrow Jaws
 Angled 45°
 Overall Length 11.5 cm (4 1/2")

FR-897-12

Shepard Lens Holding Forcep
 Gently Curved, Duckbill Shaped Jaws
 Overall Length 12.5 cm (5")



Eye Surgery Forceps
Augenchirurgie Pinzetten
Ojo Pinzas Cirugía
Pince de chirurgie oculaire
Eye Pinze Chirurgia



FR-899-12

Clayman Lens Holding Forcep
Gently Curved, Blades With
Delicately Angled Jaws
Overall Length 12.5 cm (5")



FR-900-10

Bangerter Muscle Clamp
With 10 Delicate Teeth
and Push Button
Overall Length 10 cm (4")



FR-901-01

Jameson Muscle Forcep
Right, With Slide Lock,
4 Teeth, Child Size
Overall Length 9.5 cm (3 ¾")




FR-901-02

Jameson Muscle Forcep
Left, With Slide Lock,
4 Teeth, Child Size
Overall Length 9.5 cm (3 ¾")

Eye Surgery Forceps
 Augenchirurgie Pinzetten
 Ojo Pinzas Cirugía
 Pince de chirurgie oculaire
 Eye Pinze Chirurgia

OP

OPHTHALMIC SURGERY INSTRUMENTS


FR-903-01

Jameson Muscle Forcep
 Right, With Slide Lock,
 6 Teeth, Adult Size
 Overall Length 10 cm (4")

FR-903-02

Jameson Muscle Forcep
 Left, With Slide Lock,
 6 Teeth, Adult Size
 Overall Length 10 cm (4")

FR-905-20

Berke Ptosis Forcep
 Longitudinally Grooved
 Jaws, With Side Lock, 20 mm
 Overall Length 11 cm (4 1/4")

FR-905-27

Berke Ptosis Forcep
 Longitudinally Grooved
 Jaws, With Side Lock, 27 mm
 Overall Length 11 cm (4 1/4")




Eye Surgery Forceps
 Augen chirurgie Pinzetten
 Ojo Pinzas Cirugía
 Pince de chirurgie oculaire
 Eye Pinze Chirurgia



FR-906-10
Putterman-Muellers Muscle Conjunctival Resection Ptosis Clamp
 With Slide Lock
 Overall Length 10 cm (4")



FR-907-11
Watzke Sleeve Spreading Forcep
 Jaw Length 10 mm
 Overall Length 11.5 cm (4 1/2")




FR-908-10
Erhardt Lid Forcep
 Overall Length 10.5 cm (4")

Eye Surgery Forceps
 Augenchirurgie Pinzetten
 Ojo Pinzas Cirugía
 Pince de chirurgie oculaire
 Eye Pinze Chirurgia

OP

OPHTHALMIC SURGERY INSTRUMENTS



FR-910-09
Downes Lid Clamp
 Overall Length 9.5 cm (3 ¾")



FR-912-08
Lambert Chalazion Forcep
 Round Jaws, Ø 8.0 mm
 Overall Length 9 cm (3 ½")



FR-912-12
Hunt Chalazion Forcep
 Round Jaws, Ø 12.0 mm
 Overall Length 9.5 cm (3 ¾")



FR-914-09
Baird Chalazion Forcep
 Small, Oval Jaws, 11 x 8 mm
 Overall Length 9.5 cm (3 ¾")

Eye Surgery Forceps
 Augen chirurgie Pinzetten
 Ojo Pinzas Cirugía
 Pince de chirurgie oculaire
 Eye Pinze Chirurgia


FR-916-09
Francis Chalazion Forcep

Medium, Oval Jaws, 14 x 11 mm
 Overall Length 9.5 cm (3 3/4")


FR-918-09
Desmarres Chalazion Forcep

Overall Length 9 cm (3 1/2")

FR-919-09

Large, Oval Jaws, 25 x 13 mm


FR-920-09
Wies Chalazion Forcep

Grooved Upper Jaws
 Overall Length 9.5 cm (3 3/4")


FR-921-09
Lordan Chalazion Forcep

Triangular Jaws
 Overall Length 9.5 cm (3 3/4")

Eye Surgery Forceps
 Augenchirurgie Pinzetten
 Ojo Pinzas Cirugía
 Pince de chirurgie oculaire
 Eye Pinze Chirurgia

OP

OPHTHALMIC SURGERY INSTRUMENTS


FR-922-01

Ayer Chalazion Forcep
 Round Jaws, Ø 12.0 mm
 Overall Length 9.5 cm (3 ¾")


FR-923-01

Desmarres Entropian Forcep
 Right
 Overall Length 10 cm (4")


FR-923-02

Desmarres Entropian Forcep
 Left
 Overall Length 10 cm (4")


FR-925-10

Collin Meibomian Gland Expression Forcep
 Ø 10 mm
 Overall Length 10 cm (4")

Bulldog Clamps, Haemostatic Forceps
 Bulldogklemmen, Arterienklemmen
 Clamps Bulldog, Pinzas hemostáticas
 Clamps Bulldog, Pinces hémostatiques
 Pinze Bulldog, Pinze emostatische



AF-062-03
Straight

AF-063-03
Curved

Dieffenbach Bulldog Clamp
3.5 cm (1 1/2")



AF-140-10

Straight
10 cm (4")

AF-141-10

Curved
9.5 cm (3 3/4")



AF-160-12
Straight
12.5 cm (4 3/4")

AF-161-12
Curved
12 cm (4 3/4")

Halsted-Mosquito Haemostatic Forcep



AF-170-10
Straight
10 cm (4")

AF-171-10
Curved
9.5 cm (3 3/4")

AF-172-10
Straight
1 x 2 Teeth
10 cm (4")

AF-173-10
Curved
1 x 2 Teeth
10 cm (4")

Micro-Mosquito Haemostatic Forcep

Towel Clamps
Tuchklemmen
Pinzas de campo
Pince à champs
Pinze fissatelle

OP

OPHTHALMIC SURGERY INSTRUMENTS

**AF-761-05****Jones Towel Clamp**

Sharp

Overall Length 5.5 cm (2 1/2")

**AF-773-08****Backhaus Towel Clamp**

Overall Length 9 cm (3 1/2")

AF-761-09

Sharp

Overall Length 8.5 cm (3 1/4")

Sphenoid Bone Punches

Keilbeinstanzen

Pinzas sacabocadas para esfenoides

Pinces emporte-pièces pour sinus sphénoidal

Pinze taglienti per seno sfenoidale


Citelli Mastoid Rongeur
 Shaft Length 17 cm (6 3/4")

□
OT-751-02

Bite Size 2.0 mm

□
OT-751-03

Bite Size 3.0 mm


Kerrison Sphenoid Punch
□
RH-305-03

Up Cutting

Bite Size Ø 3.0 mm x 3.0 mm

Shaft Length 14 cm (5 1/2")

□
RH-305-04

Up Cutting

Bite Size Ø 4.0 mm x 4.0 mm

Shaft Length 9 cm (3 1/2")

□
RH-305-05

Up Cutting

Bite Size Ø 5.0 mm x 5.0 mm

Shaft Length 9 cm (3 1/2")

□
RH-305-06

Up Cutting

Bite Size Ø 6.0 mm x 6.0 mm

Shaft Length 9 cm (3 1/2")

Septum Chisels, Septum Gouges, Nasal Specula
 Septummeißel, Septum-Hohlmeißel, Nasenspekula
 Cinceles para tabique nasal, Gubias para tabique nasal, Espéculos nasales
 Élévateurs pour le périoste, Gouges pour la cloison nasale, Spécums nasaux
 Elevatori per il perioste, Sgorbie per il setto nasale, Specoli nasali

OP

OPHTHALMIC SURGERY INSTRUMENTS



Hartmann-Halle Nasal Speculum
 Overall Length 15.5 cm (6")




RH-033-01

27 mm, Fig. 1



RH-033-02

31 mm, Fig. 2



RH-033-03

35 mm, Fig. 3



RH-475-18

Freer Septum Elevator

Double Ended, Sharp/Blunt, Blades Size 5.0 mm, 5.0 mm
 Overall Length 18.5 cm (7 1/4")



RH-600-16

Freer Septum Chisel

Straight, 4.0 mm

Overall Length 16.5 cm (6 1/2")



RH-601-16

Freer Septum Chisel

Curved, 4.0 mm

Overall Length 16.5 cm (6 1/2")



RH-604-16

West Septum Chisel

Overall Length 16 cm (6 1/4")



RH-605-16

West Septum Gouge

Overall Length 16 cm (6 1/4")


Micro Needle Holders
 Mikro-Nadelhalter
 Micro Porta-agujas
 Micro Porte-aiguilles
 Micro Portagli



Halsey Needle Holder
 Overall Length 13 cm (5")



SU-024-13
 Smooth Jaws **SU-026-13**
 Cross Serrated Jaws



SU-040-13
Brown Needle Holder
 Overall Length 13 cm (5")



McPherson Micro Needle Holder
 Very Delicate 0.5 x 9 mm Tips, Without Lock
 Overall Length 10 cm (4")



SU-200-10
 Straight **SU-201-10**
 Curved



McPherson Micro Needle Holder
 Very Delicate 0.5 x 9 mm Tips, With Lock
 Overall Length 10 cm (4")



SU-202-10
 Straight **SU-203-10**
 Curved

Micro Needle Holders
Mikro-Nadelhalter
Micro Porta-agujas
Micro Porte-aiguilles
Micro Portagli

OP

OPHTHALMIC SURGERY INSTRUMENTS



SU-210-14
Straight



SU-211-14
Curved



SU-212-14
Straight



SU-213-14
Curved



Castroviejo Micro Needle Holder

Delicate 1.0 x 12 mm Tips, Without Lock
Overall Length 14 cm (5 ½")

Castroviejo Micro Needle Holder

Delicate 1.0 x 12 mm Tips, With Lock
Overall Length 14 cm (5 ½")

Micro Needle Holders
 Mikro-Nadelhalter
 Micro Porta-agujas
 Micro Porte-aiguilles
 Micro Portagli



SU-220-14
Straight



SU-221-14
Curved



SU-222-14
Straight



SU-223-14
Curved

Castroviejo Micro Needle Holder

Standard 1.5 x 12 mm Tips, Without Lock
 Overall Length 14 cm (5 1/2")

Castroviejo Micro Needle Holder

Standard 1.5 x 12 mm Tips, With Lock
 Overall Length 14 cm (5 1/2")

Micro Needle Holders
 Mikro-Nadelhalter
 Micro Porta-agujas
 Micro Porte-aiguilles
 Micro Portagli

OP

OPHTHALMIC SURGERY INSTRUMENTS



SU-300-10
Straight



SU-301-10
Curved



SU-310-10
Straight
With Lock
10.5 cm (4")



SU-310-11
Straight
With Lock
11.5 cm (4 1/2")



SU-311-10
Curved
With Lock
10.5 cm (4")



SU-311-11
Curved
With Lock
11.5 cm (4 1/2")

Barraquer-Troutman Micro Needle Holder
 Very Delicate 0.6 x 8 mm Tips, Without Lock
 Round Handle
 Overall Length 10 cm (4")

Barraquer-Troutman Micro Needle Holder
 Very Delicate 0.6 x 8 mm Tips, Round Handle

Micro Needle Holders
 Mikro-Nadelhalter
 Micro Porta-agujas
 Micro Porte-aiguilles
 Micro Portagli



SU-313-11
Curved



Barraquer Micro Needle Holder
 Very Delicate 0.8 x 10 mm Tips
 Round Handle
 Overall Length 11.5 cm (4 1/2")

SU-315-11
Curved



Barraquer Micro Needle Holder
 Very Delicate 0.8 x 10 mm Tips
 Round Handle, With Lock
 Overall Length 11.5 cm (4 1/2")

SU-315-11 T
Curved



Barraquer Micro Needle Holder
 Very Delicate 0.8 x 10 mm Tips
 Round Handle, With Lock
 Overall Length 11.5 cm (4 1/2")



Micro Needle Holders
Mikro-Nadelhalter
Micro Porta-agujas
Micro Porte-aiguilles
Micro Portagli

OP

OPHTHALMIC SURGERY INSTRUMENTS



SU-320-13
Straight



SU-321-13
Curved



SU-322-13
Straight



SU-323-13
Curved

Barraquer Micro Needle Holder

Very Delicate 0.8 x 11 mm Tips, Without Lock
Round Handle
Overall Length 13 cm (5")


Barraquer Micro Needle Holder

Very Delicate 0.8 x 11 mm Tips, With Lock
Round Handle
Overall Length 13 cm (5")


Micro Needle Holders
 Mikro-Nadelhalter
 Micro Porta-agujas
 Micro Porte-aiguilles
 Micro Portagli




Barraquer Micro Needle Holder
 Very Delicate 0.8 x 12 mm Tips
 Round Handle
 Overall Length 14 cm (5 1/2")



SU-331-14
 Straight



Barraquer Micro Needle Holder
 Very Delicate 0.8 x 12 mm Tips
 Round Handle, Titanium
 Overall Length 14 cm (5 1/2")



SU-331-14 T
 Straight



Barraquer Micro Needle Holder
 Very Delicate 0.8 x 12 mm Tips
 Round Handle, With Lock
 Overall Length 14 cm (5 1/2")




SU-333-14
 Curved



Micro Needle Holders
 Mikro-Nadelhalter
 Micro Porta-agujas
 Micro Porte-aiguilles
 Micro Portagli

OP

OPHTHALMIC SURGERY INSTRUMENTS



Barraquer Micro Needle Holder
 Very Delicate 0.8 x 12 mm Tips
 Round Handle, With Lock, Titanium
 Overall Length 14 cm (5 1/2")



Troutman Micro Needle Holder
 Very Delicate 0.6 x 10 mm Tips Strongly Curved
 Round Handle
 Overall Length 12 cm (4 3/4")



Micro Needle Holders
 Mikro-Nadelhalter
 Micro Porta-agujas
 Micro Porte-aiguilles
 Micro Portagli



SU-500-13
Curved

UltraGrip™ TC Barraquer Micro Needle Holder
 1.0 x 8 mm Tungsten Carbide Tip, Smooth Jaws
 Without Lock
 Overall Length 13 cm (5")

ULTRA GRIP™



SU-502-13
Curved

UltraGrip™ TC Barraquer Micro Needle Holder
 1.0 x 8 mm Tungsten Carbide Tip, Smooth Jaws
 Without Lock
 Overall Length 13 cm (5")

ULTRA GRIP™



Micro Needle Holders
 Mikro-Nadelhalter
 Micro Porta-agujas
 Micro Porte-aiguilles
 Micro Portagli

OP

OPHTHALMIC SURGERY INSTRUMENTS



UltraGrip™ TC Castroviejo Micro Needle Holder
 1.5 x 8 mm Tungsten Carbide Tip, With Lock
 Overall Length 14.5 cm (5 $\frac{3}{4}$ ")



UltraGrip™ TC Castroviejo Micro Needle Holder
 1.5 x 8 mm Tungsten Carbide Tip, With Lock
 Overall Length 14.5 cm (5 $\frac{3}{4}$ ")

ULTRA GRIP X

ULTRA GRIP X


Micro Needle Holders
 Mikro-Nadelhalter
 Micro Porta-agujas
 Micro Porte-aiguilles
 Micro Portagli



SU-512-14
Straight



SU-513-14
Curved

UltraGrip™ TC Castroviejo Micro Needle Holder
 1.5 x 8 mm Tungsten Carbide Tip, Smooth Jaws,
 With Lock
 Overall Length 14.5 cm (5 ¾")

UltraGrip™ TC Castroviejo Micro Needle Holder
 1.5 x 8 mm Tungsten Carbide Tip, Smooth Jaws,
 With Lock
 Overall Length 14.5 cm (5 ¾")

ULTRA GRIP™

ULTRA GRIP™



Calipers
Schieblehren
Calibres
Pieds à coulisse
Calibro



DI-850-20
Castroviejo Caliper
Straight
Measuring Range 20 mm
Overall Length 14.5 cm (5 ¾")



DI-880-15
Ruler
Graduated in mm & inches
Measuring Range 150 mm
Overall Length 10 cm (4")



Graefe Nerve Hook
Blunt
Overall Length 14 cm (5 ½")



Knapp Fine Wound Retractor
Overall Length 13.5 cm (5 ¼")



Joseph Mucus Hook
Overall Length 16.5 cm (6 ½")



Desmarres Lid Retractor
Overall Length 16 cm (6 ¼")



Bone Rongeurs
Hohlmeißelzangen
Pinzas gubias
Pinces-gouges
Pinze-sgorbie



OR-470-19

Cottle Bone Mallet

Plain, Convex, Ø 30.0 mm,
300 Grams
Overall Length 19 cm (7 1/2")


Friedman Bone Rongeur

Overall Length 13.5 cm (5 1/4")



OR-933-01



OR-934-02



OR-937-16

Mead Bone Rongeur

Overall Length 16.5 cm (6 1/2")




PR-200-00

Bowman Lacrimal Probes Complete
Set of 6 Probes

Sterilizing Cases
Sterilisierbehälter Koffer
Los casos de esterilización
Stérilisation cas
Sterilizzazione Custodie i

OP

OPHTHALMIC SURGERY INSTRUMENTS

**ST-110-01****Storing Case Plastic**


Single Level With Perforated
Silicone Mat Size 160 x 70 x 25

ST-110-02**Storing Case Plastic**

Single Level With Perforated
Silicone Mat Size 200 x 110 x 25

ST-110-03**Storing Case Plastic**

Single Level With Perforated
Silicone Mat Size 265 x 160 x 25

**ST-112-26****Storing Case Plastic**

Double Level With Perforated
Silicone Mat Size 265 x 160 x 50

Ophthalmic Surgery Sets



SURTEX
INSTRUMENTS

Crafted to Excellence

Basic Eye Surgery Set

Codes	Description	Quantity
SL-010-03	Bard-Parker Scalpel Handle No. 3 Solid 12 cm - 4 3/4"	1
SL-030-07	Scalpel Handle No. 7 Solid 16 cm - 6 1/4"	1
SL-206-13	Scalpel Handle For Micro Scalpel Blades 13 cm - 5"	1
SS-430-10	Iris Scissor Straight 10 cm - 4"	1
SS-431-10	Iris Scissor Curved 10 cm - 4"	1
AF-170-12	Micro-Mosquito Haemostatic Forcep Straight 12 cm - 4 3/4"	4
AF-171-12	Micro-Mosquito Haemostatic Forcep Curved 11.5 cm - 4 1/2"	6
AF-773-08	Backhaus Towel Clamp 9 cm - 3 1/2"	4
SS-423-09	Very Delicate Operating Scissor Curved 9 cm - 3 1/2"	1
OP-010-15	Barraquer Wire Speculum Standard Wire Blade Size 15 mm	1
OP-083-15	Lieberman Eye Speculum Adjustable Blade Size 15 mm	1
OP-091-14	Castroviejo Eye Speculum Adjustable Mechanism Blade Size 14 mm	1
OP-175-15	Desmarres Lid Retractor Thin Solid Blades - Size 2 13 cm - 5" Blade Width 15 mm	1
RT-048-02	Graefe Nerve Hook Fig. 2 14 cm - 5 1/2"	1
OP-745-50	Castroviejo Cyclodialysis Spatula 13 cm - 5" Width 0.5 mm	1
DI-850-20	Castroviejo Caliper Straight 8 cm - 3" Measuring Range 20 mm	1
FR-700-01	Graefe Eye Dressing Forcep 10 cm - 4"	1
FR-700-02	Graefe Eye Dressing Forcep 10 cm - 4"	1
FR-702-01	Graefe Dissecting Forcep 1 x 2 Teeth 10 cm - 4"	1
FR-702-02	Graefe Dissecting Forcep 1 x 2 Teeth 10 cm - 4"	1
FR-772-01	Graefe Fixation Forcep 11 cm - 4 1/4"	1
FR-750-03	Castroviejo Suture Tying Forcep 1 x 2 Teeth with Tying Platform 11 cm - 4 1/4" Tip Size 0.3 mm	1
FR-859-11	Arruga Capsule Forcep 11 cm - 4 1/4"	1
FR-918-09	Desmarres Chalazion Forcep Large - Oval Jaws 9 cm - 3 1/2" Jaw Size 20 x 13 mm	1
SS-453-11	Stevens Tenotomy Scissor Curved - Blunt/Blunt 11 cm - 4 1/2"	1
SS-653-01	Castroviejo Corneoscleral Scissor Right - Medium Blades 11 cm - 4 1/2"	1
SS-653-02	Castroviejo Corneoscleral Scissor Left - Medium Blades 11 cm - 4 1/2"	1
SS-630-12	Noyes Iris Scissor Straight - Sharp/Sharp 12.5 cm - 5"	1
SS-631-12	Noyes Iris Scissor Curved - Sharp/Sharp 12.5 cm - 5"	1
SS-683-11	Westcott Stitch Scissor Curved - Very Sharp Pointed Tips - Standard Blades 11.5 cm - 4 1/2"	1
SU-222-14	Castroviejo Micro Needle Holder Straight - Standard - With Lock 14 cm - 5 1/2"	1
OP-180-02	Scleral Shield Large (Pair) Blade Width 22 x 30 mm	2
ST-112-26	Storing Case Double Level 265 x 160 x 50 mm Plastic	1

Planned Extra-Capsular Cataract Extraction

Codes	Description	Quantity
OP-471-19	Bishop-Harmon Anterior Chamber Irrigator Gauge 19	1
AF-761-05	Jones Towel Clamp Sharp 5.5 cm - 2 1/4"	2
AF-140-10	Hartmann Haemostatic Forcep Straight 10 cm - 4"	1
AF-062-03	Dieffenbach Bulldog Clamp Straight 3.5 cm - 1 1/2"	1
OP-083-15	Lieberman Eye Speculum Adjustable Blade Size 15 mm	1
OP-010-15	Barraquer Wire Speculum Standard Wire Blade Size 15 mm	1
SL-206-13	Scalpel Handle For Micro Scalpel Blades 13 cm - 5"	2
SL-210-69 A	Micro Scalpel Blade No. 69A Pack of 25	1
OP-232-15	Disposable Stab Knife Straight 15 Degree Cutting Edge Box of 6	1
OP-320-01	Diamond Step Knife With Adjustable Cutting Depth of 0/0.3/0.4/0.5/0.6/0.8/1.0/2.5 mm	1
OP-865-13	Single 45° Cutting Edge 12 cm - 4 3/4" Width 1.0 mm Titanium Handle	1
SS-651-01	Fine-Thronton Swivel Fixation Ring With Swivel Handle 12 cm - 4 3/4" Tip Ø 13 mm	1
SS-651-02	Castroviejo Corneoscleral Scissor Right - Small Blades 10.5 cm - 4"	1
SS-675-11	Castroviejo Corneoscleral Scissor Left - Small Blades 10.5 cm - 4"	1
SS-683-11	Westcott Tenotomy Scissor Curved - Blunt Tips - Medium Blades 11 cm - 4 1/2"	1
SS-695-08	Westcott Stitch Scissor Curved - Very Sharp Pointed Tips - Standard Blades 11.5 cm - 4 1/2"	1
SS-713-11	Vannas Capsulotomy Scissor Curved - Sharp Tips 8 cm - 3 1/4 Blade Size 5 mm	1
DI-850-20	Stern-Gills Micro Scissor Angled Forward - Sharp Tips - Extra Thin 11 cm - 4 1/2" Blade Size 10 mm	1
FR-834-12	Castroviejo Caliper Straight 8 cm - 3" Measuring Range 20 mm	1
FR-718-07	Bonaccolto Conjunctival Forcep Longitudinal Serrated Jaws with Cross Serrations at Tips 10.5 cm - 4"	1
FR-732-10	Tip Width 1.2 mm Stainless Steel	1
FR-760-08	Harms-Colibri Corneal Forcep Very Delicate 1 x 2 Teeth with Tying Platform 7.5 cm - 3 1/4" Tip Size 0.12 mm	1
FR-780-10	Bonn Corneal Forcep Very Delicate 1 x 2 Teeth with Tying Platform 10 cm - 4" Tip Size 0.12 mm	1
FR-781-10	Bishop-Harmon Iris Forcep Standard 1 x 2 Teeth 8.5 cm - 3 1/4"	1
FR-784-10	McPherson Suture Tying Forcep Straight - With Tying Platform 10 cm - 4" Jaws Length 5 mm	1
FR-801-10	McPherson Suture Tying Forcep Angled - With Tying Platform 10 cm - 4" Jaws Length 5 mm	1
OP-677-11	Kelman-McPherson Suture Tying Forcep Angled - With Tying Platform 10 cm - 4" Jaws Length 7.5 mm	1
OP-826-04	Utrata Capsulorrhesis Forcep Very Delicate Triangular Grasping Tips - Extremely Thin 10.5 cm - 4" Shanks	1
OP-827-02	Length 11 mm	1
OP-635-02	Nagahara Nucleus Chopper And Spatula 11 cm - 4 1/4"	1
SU-313-11	Sinskey II Lens Manipulating Hook Angled 11.5 cm - 4 1/2" Tip Diameter 0.2 mm	1
SU-213-14	Kuglen Iris Hook and Lens Manipulator Push Pull Modle - Angled 11.5 cm - 4 1/2"	1
OP-750-02	Barraquer Iris Spatula Round - Smooth 12 cm - 4 3/4" Width 0.5 mm	2
OP-765-10	Barraquer Micro Needle Holder Curved - Delicate - Round Handle 11.5 cm - 4 1/2"	1
OP-775-02	Castroviejo Micro Needle Holder Curved - Delicate - With Lock 14 cm - 5 1/2"	1
OP-869-10	Lewis Lens Loop Large 14 cm - 5 1/2"	1
OP-542-25	Alfonso Nucleus Trisector For Use With Ref:- OP-775-02 10.5 cm - 4 1/4" Width 1.5 mm	1
OP-542-40	Alfonso Cutting Platform 12.5 cm - 5" Width 2 mm Stainless Steel	1
OP-519-30	Alfonso Nucleas Grasping Forcep With Two Short Rows of Delicate Interlocking Teeth 10.5 cm - 4"	1
OP-421-27	Alfonso Nucleus Irrigating Spoon 2.5 mm Wide	1
OP-517-26	Alfonso Nucleus Irrigating Spoon 4.0 mm Wide	1
OP-467-23	Gimbel Cortex Aspirating Cannula End Opening Gauge 30	1
OP-451-23	Air Injection Cannula Angled at 5 mm Gauge 27	1
ST-112-26	McIntyre Anterior Chamber Cannula Angled Gauge 26	1
OP-451-23	Jansen Capsul Polisher Cannula 23 Gauge	1
OP-451-23	Simcoe Irrigation Aspirating Cannula Original Modle Gauge 23	1
OP-451-23	Storing Case Double Level 265 x 160 x 50 mm Plastic	1

Corneal Graft Surgery Set

Codes	Description	Quantity
OP-471-19	Bishop-Harmon Anterior Chamber Irrigator Gauge 19	1
AF-761-05	Jones Towel Clamp Sharp 5.5 cm - 2 1/4"	1
OP-100-00	Flieringa Scleral Fixation Ring Set of 9 Ref:- OP-100-14 to OP-100-22	1
OP-111-16	McNeill-Goldman Scleral Fixation Ring & Blepharostat Medium, I.D Blade Size 16 mm	1
OP-010-15	Barraquer Wire Speculum Standard Wire Blade Size 15 mm	1
OP-087-15	Lieberman Eye Speculum Adjustable Blade Size 15 mm	1
AF-140-10	Hartmann Haemostatic Forcep Straight 10 cm - 4"	2
AF-062-03	Dieffenbach Bulldog Clamp Straight 3.5 cm - 1 1/2"	2
SS-659-01	Troutman-Katzin Corneal Transplant Scissor Right - Strongly Curved - Small Blades - With Lock 10.5 cm - 4"	1
SS-659-02	Troutman-Katzin Corneal Transplant Scissor Left - Strongly Curved - Small Blades - With Lock 10.5 cm - 4"	1
SS-675-11	Westcott Tenotomy Scissor Curved - Blunt Tips - Medium Blades 11 cm - 4 1/2"	1
SS-685-11	Jaffe Stitch Scissor Curved - Very Sharp Pointed Tips - Medium Blades 11 cm - 4 1/2"	1
OP-685-11	Rosen Phaco Splitter Angled 60° - Wedge Shaped 11.5 cm - 4 1/2"	1
OP-400-06	Corneal Trehine Blade Diameter 7.5 mm	1
OP-400-08	Corneal Trehine Blade Diameter 8.0 mm	1
OP-401-06	Corneal Trehine Handle Diameter 7.5 mm	1
OP-401-08	Corneal Trehine Handle Diameter 8.0 mm	1
FR-720-07	Colibri Corneal Forcep Standard 1 x 2 Teeth with Tying Platform 7.5 cm - 3 1/4" Tip Size 0.4 mm	1
FR-750-01	Castroviejo Suture Tying Forcep 1 x 2 Teeth with Tying Platform 11 cm - 4 1/4" Tip Size 0.12 mm	1
FR-750-03	Castroviejo Suture Tying Forcep 1 x 2 Teeth with Tying Platform 11 cm - 4 1/4" Tip Size 0.3 mm	1
FR-718-07	Harms-Colibri Corneal Forcep Very Delicate 1 x 2 Teeth with Tying Platform 7.5 cm - 3 1/4" Tip Size 0.12 mm	1
FR-816-10	Tennant Suture Tying Forcep Straight - Round Handle with Guide Pin - Extra Delicate Smooth Jaws 10.5 cm - 4" Jaws Length 6 mm	1
FR-817-10	Tennant Suture Tying Forcep Curved - Round Handle with Guide Pin - Extra Delicate Smooth Jaws 10.5 cm - 4" Jaws Length 6 mm	1
OP-793-14	Paton Spatula & Spoon For Corneal Button Transfer 14 cm - 5 1/2"	1
SU-323-13	Barraquer Micro Needle Holder Curved - Very Delicate - Round Handle - With Lock 13 cm - 5"	1
SU-213-14	Castroviejo Micro Needle Holder Curved - Delicate - With Lock 14 cm - 5 1/2"	1
OP-421-30	Air Injection Cannula Angled at 5 mm Gauge 30	1
OP-517-26	McIntyre Anterior Chamber Cannula Angled Gauge 26	1
ST-112-26	Storing Case Double Level 265 x 160 x 50 mm Plastic	1

Phacoemulsification Cataract Extraction Set

Codes	Description	Quantity
OP-471-19	Bishop-Harmon Anterior Chamber Irrigator Gauge 19	1
OP-083-15	Lieberman Eye Speculum Adjustable Blade Size 15 mm	2
OP-030-15	Barraquer Wire Speculum Open Blades Blade Size 15 mm	1
OP-316-03	Diamond Knife Tri-Facet Cutting Edge 12.5 cm - 5" Width 1.0 mm Titanium Handle	1
OP-341-28	Diamond Phaco Knife Angled - Spear Shaped 13 cm - 5" Blade Width 2.8 mm Titanium Handle	1
OP-230-02	Disposable Crescent Knife Angled Box of Six	1
OP-232-15	Disposable Stab Knife Straight 15 Degree Cutting Edge Box of 6	1
SL-206-13	Scalpel Handle For Micro Scalpel Blades 13 cm - 5"	1
AF-140-10	Hartmann Haemostatic Forcep Straight 10 cm - 4"	2
AF-761-05	Jones Towel Clamp Sharp 5.5 cm - 2 1/4"	2
OP-659-02	Chang Pahco Chopper Left Hand 13 cm - 5" Tip Length 0.9 mm	1
OP-532-27	Chang Hydrodissection Cannula Flat Tip With Bevelled Opening 27 Gauge	1
OP-683-11	Rosen Phaco Splitter Wedge Shaped 11.5 cm - 4 1/2"	1
OP-827-02	Kuglen Iris Hook and Lens Manipulator Push Pull Modle - Angled 11.5 cm - 4 1/2"	1
OP-689-02	Jaffe-Bechert Nucleus Rotator Blunt Forked Tip - Angled At 9 mm 11.5 cm - 4 1/2"	1
OP-729-12	Lester IOL Manipulator Angled 12 cm - 4 3/4" Tip Diameter 0.25 mm	1
SS-675-11	Westcott Tenotomy Scissor Curved - Blunt Tips - Medium Blades 11 cm - 4 1/2"	1
DI-850-20	Castroviejo Caliper Straight 8 cm - 3" Measuring Range 20 mm	1
SS-713-11	Stern-Gills Micro Scissor Angled Forward - Sharp Tips - Extra Thin 11 cm - 4 1/2" Blade Size 10 mm	1
FR-716-07	Colibri Corneal Forcep Very Delicate 1 x 2 Teeth 7.5 cm - 3 1/4" Tip Size 0.12 mm	1
FR-750-01	Castroviejo Suture Tying Forcep 1 x 2 Teeth with Tying Platform 11 cm - 4 1/4" Tip Size 0.12 mm	1
FR-724-10	Fechtner Conjunctival Forcep Very Delicate Ring Shaped Tips with Tying Platform 10.5 cm - 4"	1
FR-816-10	Tennant Suture Tying Forcep Straight - Round Handle with Guide Pin - Extra Delicate Smooth Jaws 10.5 cm - 4"	1
	Jaws Length 6 mm	
FR-817-10	Tennant Suture Tying Forcep Curved - Round Handle with Guide Pin - Extra Delicate Smooth Jaws 10.5 cm - 4"	1
	Jaws Length 6 mm	
FR-801-10	Utrata Capsulorrhesis Forcep Very Delicate Triangular Grasping Tips - Extremely Thin 10.5 cm - 4" Shanks Length 11 mm	1
FR-784-10	Kelman-McPherson Suture Tying Forcep Angled - With Tying Platform 10 cm - 4" Jaws Length 7.5 mm	1
FR-885-10	Livernois-McDonald IOL Inserting Forcep Concave Jaws 10.5 cm - 4" Jaws Length 5 mm	1
FR-877-10	Livernois IOL Folding Forcep Concave Jaw Surfaces - For Soft IOLs 10.5 cm - 4"	1
SU-315-11	Barraquer Micro Needle Holder Curved - Delicate - Round Handle - With Lock 11.5 cm - 4 1/2"	1
FR-863-11	Ernest Nucleus Cracker Extra Delicate Cross Action Serrated Paddle Tips 11 cm - 4 1/4"	2
OP-451-23	Simcoe Irrigation Aspirating Cannula Original Modle Gauge 23	1
OP-475-25	Hydrodissection Cannula Angled at 8 mm - Flat Tip Gauge 25	1
OP-421-27	Air Injection Cannula Angled at 5 mm Gauge 27	1
OP-467-23	Jansen Capsul Polisher Cannula 23 Gauge	1
OP-441-19	Stevens Sub Tenon's Anesthesia Cannula Gently Curved - Flat Tip Gauge 19	1
OP-865-13	Fine-Thronton Swivel Fixation Ring With Swivel Handle 12 cm - 4 3/4" Tip Diameter 13 mm	1
OP-465-23	McIntyre Cystotome Gauge 23	1
SS-695-08	Vannas Capsulotomy Scissor Curved - Sharp Tips 8 cm - 3 1/4 Blade Size 5 mm	1
ST-112-26	Storing Case Double Level 265 x 160 x 50 mm Plastic	1

Phacofracture Cataract Extraction Set

Codes	Description	Quantity
OP-085-15	Lieberman Eye Speculum Adjustable Blade Size 15 mm	1
SS-440-11	Iris Scissor Straight 11.5 cm - 4 1/2"	1
SS-675-11	Westcott Tenotomy Scissor Curved - Blunt Tips - Medium Blades 11 cm - 4 1/2"	1
FR-750-01	Castroviejo Suture Tying Forcep 1 x 2 Teeth with Tying Platform 11 cm - 4 1/4" Tip Size 0.12 mm	1
SL-206-13	Scalpel Handle For Micro Scalpel Blades 13 cm - 5"	1
SL-210-69	A Micro Scalpel Blade No. 69A Pack of 25	1
OP-232-15	Disposable Stab Knife Straight 15 Degree Cutting Edge Box of 6	1
DI-850-20	Castroviejo Caliper Straight 8 cm - 3" Measuring Range 20 mm	1
FR-801-10	Ultrata Capsulorrhesis Forcep Very Delicate Triangular Grasping Tips - Extremely Thin 10.5 cm - 4" Shanks Length 11 mm	1
OP-827-02	Kuglen Iris Hook and Lens Manipulator Push Pull Modle - Angled 11.5 cm - 4 1/2"	1
OP-775-02	Alfonso Cutting Platform 12.5 cm - 5" Width 2 mm	1
OP-775-03	Alfonso Cutting Platform 12.5 cm - 5" Width 3 mm	1
OP-765-10	Alfonso Nucleus Trisector For Use With Ref:- OP-775-02 10.5 cm - 4 1/4" Width 1.5 mm	1
OP-767-10	Bucher-Kansas Nucleus Trisector For Use With Ref:- OP-775-03 10.5 cm - 4 1/4" Width 2.5 mm	1
OP-771-11	Kansas Nucleus Bisector 11.5 cm - 4 1/2" Stainless Steel	1
FR-869-10	Alfonso Nucleus Grasping Forcep With Two Short Rows of Delicate Interlocking Teeth 10.5 cm - 4"	1
OP-542-25	Alfonso Nucleus Irrigating Spoon 2.5 mm Wide	1
OP-542-40	Alfonso Nucleus Irrigating Spoon 4.0 mm Wide	1
OP-421-27	Air Injection Cannula Angled at 5 mm Gauge 27	1
OP-525-21	Simcoe Cortex Extractor Smooth Closed Tip - 21 mm Long Gauge 21	1
FR-875-11	Dodick IOL Folding Forcep For Soft IOLs 11 cm - 4 1/4"	1
FR-889-10	Buratto Soft IOL Inserting Forcep Extra Delicate Highly Polished Convex Jaws 10 cm - 4"	1
ST-112-26	Storing Case Double Level 265 x 160 x 50 mm Plastic	1

Glaucoma (Optic Nerve Repair) Surgery Set

Codes	Description	Quantity
OP-471-19	Bishop-Harmon Anterior Chamber Irrigator Gauge 19	1
AF-761-05	Jones Towel Clamp Sharp 5.5 cm - 2 1/4"	1
AF-140-10	Hartmann Haemostatic Forcep Straight 10 cm - 4"	1
OP-064-15	Murdoch Eye Speculum Blade Size 15 mm	1
OP-010-15	Barraquer Wire Speculum Standard Wire Blade Size 15 mm	1
OP-083-15	Lieberman Eye Speculum Adjustable Blade Size 15 mm	1
SS-440-11	Iris Scissor Straight 11.5 cm - 4 1/2"	1
SS-675-11	Westcott Tenotomy Scissor Curved - Blunt Tips - Medium Blades 11 cm - 4 1/2"	1
SS-647-07	Barraquer Iris Scissor Angled 5.5 cm - 2 1/4" Blade Size 7 mm	1
SS-707-10	Gills-Vannas Micro Scissor Curved - Sharp Tips - Extra Thin 10.5 cm - 4" Blade Size 7 mm	1
SS-681-11	McPherson-Westcott Stitch Scissor Curved - Very Sharp Pointed Tips - Small Blades 11 cm - 4 1/2"	1
FR-724-10	Fechtner Conjunctival Forcep Very Delicate Ring Shaped Tips with Tying Platform 10.5 cm - 4"	2
FR-716-07	Colibri Corneal Forcep Very Delicate 1 x 2 Teeth 7.5 cm - 3 1/4" Tip Size 0.12 mm	1
FR-760-08	Bishop-Harmon Iris Forcep Standard 1 x 2 Teeth 8.5 cm - 3 1/4"	1
OP-919-09	Fechtner Trabeclectomy Marker 9 cm - 3 1/2"	1
DI-850-20	Castroviejo Caliper Straight 8 cm - 3" Measuring Range 20 mm	1
SL-206-13	Scalpel Handle For Micro Scalpel Blades 13 cm - 5"	2
OP-230-02	Disposable Crescent Knife Angled Box of Six	1
SL-210-69 A	Micro Scalpel Blade No. 69A Pack of 25	1
OP-232-15	Disposable Stab Knife Straight 15 Degree Cutting Edge Box of 6	1
OP-395-00	Luntz-Dodick Trabeclectomy Punch Complete With Punch Ref:- OP-455-01 Diameter - Deep Bite 1.0 mm - 0.5 mm	1
OP-397-00	Kelly Descemet's Membrane Punch Diameter - Deep Bite 1.0 mm - 0.5 mm	1
OP-587-13	Blumenthal Conjunctiva Dissector With Blunt Disc Shaped Tip 13 cm - 5"	1
SU-321-13	Barraquer Micro Needle Holder Curved - Very Delicate - Round Handle 13 cm - 5"	1

Glaucoma (Optic Nerve Repair) Surgery Set

Codes	Description	Quantity
OP-421-27	Air Injection Cannula Angled at 5 mm Gauge 27	1
OP-941-01	Harms Trabeculotomy Probe Right 5 cm - 2"	1
OP-941-02	Harms Trabeculotomy Probe Left 5 cm - 2"	1
ST-112-26	Storing Case Double Level 265 x 160 x 50 mm Plastic	1

LASIK Set

Codes	Description	Quantity
OP-083-15	Lieberman Eye Speculum Adjustable Blade Size 15 mm	1
OP-083-15 A	Lieberman Eye Speculum With Aspirating Ports Adjustable Blade Size 15 mm	1
OP-591-12	Buratto Lasik Flap Protector 12.5 cm - 5"	1
OP-597-11	Mannis-Buratto Lasik Flap Protector Blades Angled 45° 11.5 cm - 4 1/2"	1
OP-713-14	Hersh Flap Spatula 14 cm - 5 1/2"	1
OP-717-12	Fukasaaku Lasik Spatula 12 cm - 4 3/4"	1
OP-720-13	Katzen Lasik Flap Unzipper 13.5 cm - 5 1/4"	1
OP-863-11	Gimbel Stabilization Ring With Swivel Handle 11.5 cm - 4 1/2" Tip Diameter 13 mm	1
OP-887-11	Chayet Corneal Marker for LASIK Surgery 11.5 cm - 4 1/2"	1
FR-793-10	Buratto LASIK Flap Forcep Curved - Disc Shaped Serrated Jaws 10.5 cm - 4"	1
OP-507-25	Banaji Lasik Cannula With 6 Ports Gauge 25	1
OP-495-26	Spatulated Lasik Cannula Flattened Spatulated Tip - End Opening Gauge 26	1
OP-509-23	Buratto Lasik Cannula With 3 Ports - Slightly Flattened Gauge 23	1
ST-110-03	Storing Case Large 265 x 160 x 25 mm Plastic	1

LASEK Surgery Set

Codes	Description	Quantity
OP-083-15 A	Lieberman Eye Speculum With Aspirating Ports Adjustable Blade Size 15 mm	1
OP-385-08	Camellin Lasek Trehpine & Alcohol Retaining Well 12.5 cm - 5" Diameter 8 mm - 8.5 mm	1
OP-385-09	Camellin Lasek Trehpine & Alcohol Retaining Well 12.5 cm - 5" Diameter 9 mm - 9.5 mm	1
OP-619-13	Sloane Lasek Micro Hoe & Epi Peeler 13.5 cm - 5 1/4"	1
OP-627-12	Sloane Lasek Flap Reppositor Gently Curved Spatulated Tip 12 cm - 4 3/4"	1
OP-439-25	Lasek Alcohol Cannula Delicate Olive Shaped Tip Gauge 25	1
ST-110-03	Storing Case Small 165 x 70 x 25 mm Plastic	1

Chalazion Surgery Set

Codes	Description	Quantity
OP-192-11	Jaeger Lid Plate 11 cm - 4 1/4"	1
OP-550-01	Meyerhoefer Chalazion Curette Size 1 12.5 cm - 5" Diameter 1.75 mm	1
OP-550-02	Meyerhoefer Chalazion Curette Size 2 12.5 cm - 5" Diameter 2.25 mm	1
OP-550-03	Meyerhoefer Chalazion Curette Size 3 12.5 cm - 5" Diameter 2.5 mm	1
OP-550-04	Meyerhoefer Chalazion Curette Size 4 12.5 cm - 5" Diameter 3.5 mm	1
SS-441-11	Iris Scissor Curved 11.5 cm - 4 1/2"	1
FR-768-09	Manhattan Eye and Ear Fixation Forcep 1 x 2 Teeth with Tying Platform 9 cm - 3 1/2" Tip Size 0.3 mm	1
FR-912-12	Hunt Chalazion Forcep Round Jaws 9.5 cm - 3 3/4" Jaw Diameter 12 mm Ø	1
FR-916-09	Francis Chalazion Forcep Medium - Oval Jaws 9.5 cm - 3 3/4" Jaw Size 14 x 11 mm	1
SU-211-14	Castroviejo Micro Needle Holder Curved - Delicate 14 cm - 5 1/2"	1
SL-206-13	Scalpel Handle For Micro Scalpel Blades 13 cm - 5"	1
SL-210-67	Micro Scalpel Blade No. 67 Pack of 25	1
ST-110-03	Storing Case Large 265 x 160 x 25 mm Plastic	1

Strabismus / Squint Surgery Set

Codes	Description	Quantity
OP-087-10	Lieberman Eye Speculum Adjustable - Solid Blades Blade Size 10 mm	1
OP-087-15	Lieberman Eye Speculum Adjustable Blade Size 15 mm	1
AF-141-10	Hartmann Haemostatic Forcep Curved 9.5 cm - 3 3/4"	2
AF-062-03	Dieffenbach Bulldog Clamp Straight 3.5 cm - 1 1/2"	2
OP-211-07	Helveston Tissue Retractor Thin Curved Blade 12.5 cm - 5" Blade Width 7 mm	1
OP-211-09	Helveston Tissue Retractor Thin Curved Blade 12.5 cm - 5" Blade Width 9 mm	1
OP-211-11	Helveston Tissue Retractor Thin Curved Blade 6.5 cm - 2 1/2" Blade Width 11 mm	2
OP-873-12	Helveston Teaser Hook Extra Delicate - Blunt 12.5 cm - 5"	2
OP-877-01	Helveston Finder Hook Large 12.5 cm - 5"	1
OP-877-02	Helveston Finder Hook Small 12 cm - 4 3/4"	1
OP-879-01	Helveston Muscle Hook Fig. 1 13 cm - 5" Tip Length 8 mm	1
OP-879-02	Helveston Muscle Hook Fig. 2 13 cm - 5" Tip Length 10 mm	1
OP-879-03	Helveston Muscle Hook Fig. 3 13 cm - 5" Tip Length 12 mm	1
OP-930-12	Helveston Scleral Ruler 15 mm Long in Notches in 5 mm increments 12 cm - 4 3/4"	1
DI-850-20	Castroviejo Caliper Straight 8 cm - 3" Measuring Range 20 mm	1
SS-675-11	Westcott Tenotomy Scissor Curved - Blunt Tips - Medium Blades 11 cm - 4 1/2"	1
SS-683-11	Westcott Stitch Scissor Curved - Very Sharp Pointed Tips - Standard Blades 11.5 cm - 4 1/2"	1
FR-732-10	Bonn Corneal Forcep Very Delicate 1 x 2 Teeth with Tying Platform 10 cm - 4" Tip Size 0.12 mm	1
FR-777-01	Moody Fixation Forcep Curved Right 11 cm - 4 1/4"	1
FR-777-02	Moody Fixation Forcep Curved Left 11 cm - 4 1/4"	1
FR-776-02	Lester Fixation Forcep Angled - 2 x 3 Teeth 9.5 cm - 3 3/4"	3
FR-824-10	Helveston Suture Tying Forcep Straight 10.5 cm - 4"	1
OP-881-02	Jameson Muscle Hook Large 13.5 cm - 5 1/4"	1
SU-222-14	Castroviejo Micro Needle Holder Straight - Standard - With Lock 14 cm - 5 1/2"	1
SU-323-13	Barraquer Micro Needle Holder Curved - Very Delicate - Round Handle - With Lock 13 cm - 5"	1
OP-471-19	Bishop-Harmon Anterior Chamber Irrigator Gauge 19	1
SL-206-13	Scalpel Handle For Micro Scalpel Blades 13 cm - 5"	1
SL-210-69	Micro Scalpel Blade No. 69 Pack of 25	1
ST-112-26	Storing Case Double Level 265 x 160 x 50 mm Plastic	1

Lacrimal (Dacryocystorhinostomy) Surgery Set

Codes	Description	Quantity
OP-471-19	Bishop-Harmon Anterior Chamber Irrigator Gauge 19	1
AF-773-08	Backhaus Towel Clamp 9 cm - 3 1/2"	2
AF-160-12	Halsted-Mosquito Haemostatic Forcep Straight 12.5 cm - 5"	2
AF-161-12	Halsted-Mosquito Haemostatic Forcep Curved 12 cm - 4 3/4"	2
RH-305-03	Kerrison Sphenoid Punch Up Cutting 14 cm - 5 1/2" Bite size 3.0 x 3.0 mm	1
RH-033-03	Hartmann-Halle Nasal Speculum Fig. 3 15.5 cm - 6" Blade Length 35 mm	1
OR-937-16	Mead Bone Rongeur 16.5 cm - 6 1/2"	1
OR-470-19	Cottle Bone Mallet Plain - Convex 19 cm - 7 1/2" Head Ø 30.0 mm - Weight 300 Grams	1
SL-206-13	Scalpel Handle For Micro Scalpel Blades 13 cm - 5"	1
SL-210-65	Micro Scalpel Blade No. 65 Pack of 25	1
OP-220-00	Stevenson Lacrimal Sac Retractor 3 x 3 Blunt Prongs - Adjustable	1
OP-410-23	Lacrimal Cannula Straight Gauge 23	1
OP-411-23	Lacrimal Cannula Curved Gauge 23	1
OP-950-14	Castroviejo Lacrimal Dilator Double Ended Fine and Meedium Taper 14 cm - 5 1/2"	1
RH-475-18	Freer Septum Elevator Double Ended - Sharp/Blunt 18.5 cm - 7 1/4" Blade Size 1 - Blade Size 2 5.0 mm - 5.0 mm	1
RH-604-16	West Septum Chisel 16 cm - 6 1/4"	1
RH-601-16	Freer Septum Chisel Curved 16 cm - 6 1/4" Blade Width 4.0 mm	1
OP-954-14	Worst Pigtail Probe With Hook 16 cm - 6 1/4"	1
OP-956-14	Crawford Retrieval Hook 14.5 cm - 5 3/4"	1
OP-960-00	Crawford Lacrimal Intubation Set Including 30.5 cm - 12" Silicone Tube and Two Olive Tip Probes Ø 0.6 mm	1
PR-200-00	Bowman Probes Complete Set of 7	1
FR-846-10	Gold Punctal Plug Forcep With Longitudinal Groove - Pointed Tips 11 cm - 4 1/4"	1
FR-700-01	Graefe Eye Dressing Forcep 10 cm - 4"	1
FR-766-09	Sauer Fixation Forcep 1 x 2 Teeth 9 cm - 3 1/2" Tip Size 0.3 mm	1
FR-834-12	Bonaccolto Conjunctival Forcep Longitudinal Serrated Jaws with Cross Serrations at Tips 10.5 cm - 4" Tip Width 1.2 mm	1
FR-750-03	Castroviejo Suture Tying Forcep 1 x 2 Teeth with Tying Platform 11 cm - 4 1/4" Tip Size 0.3 mm	1
RT-068-13	Knapp Fine Wound Retractor 4 Blunt Prongs 13.5 cm - 5 1/4"	1
SS-453-11	Stevens Tenotomy Scissor Curved - Blunt/Blunt 11 cm - 4 1/2"	1
SS-440-11	Iris Scissor Straight 11.5 cm - 4 1/2"	1
SU-213-14	Castroviejo Micro Needle Holder Curved - Delicate - With Lock 14 cm - 5 1/2"	1
SS-683-11	Westcott Stitch Scissor Curved - Very Sharp Pointed Tips - Standard Blades 11.5 cm - 4 1/2"	1
SS-679-11	Tibolt Punctal Scissor One Conically Shaped Blade with Micro Grooves and Textured Surface 11.5 cm - 4 1/2"	1
ST-112-26	Storing Case Double Level 265 x 160 x 50 mm Plastic	1

Enucleation / Evisceration Set

Codes	Description	Quantity
OP-078-00	Lester-Burch Eye Speculum Adjustable	1
AF-161-12	Halsted-Mosquito Haemostatic Forcep Curved 12 cm - 4 3/4"	1
OP-802-15	Wells Enucleation Spoon 15 cm - 6"	1
SS-469-11	Enucleation Scissor Curved 11.5 cm - 4 1/2"	1
OP-800-01	Bunge Evisceration Spoon Small 13.5 cm - 5 1/4"	1
OP-800-02	Bunge Evisceration Spoon Small 14 cm - 5 1/2"	1
RT-048-02	Graefe Nerve Hook Fig. 2 14 cm - 5 1/2"	1
SS-452-11	Stevens Tenotomy Scissor Straight - Blunt/Blunt 11 cm - 4 1/2"	1
FR-702-01	Graefe Dissecting Forcep 1 x 2 Teeth 10 cm - 4"	1
AF-161-12	Halsted-Mosquito Haemostatic Forcep Curved 12 cm - 4 3/4"	1
ST-110-03	Storing Case Large 265 x 160 x 25 mm Plastic	1

Pterygium Excision Set

Codes	Description	Quantity
SL-206-13	Scalpel Handle For Micro Scalpel Blades 13 cm - 5"	1
SL-210-64	Micro Scalpel Blade No. 64 Pack of 25	1
OP-076-16	Lancaster Eye Speculum Adjustable Blade Size 16 mm	1
FR-780-10	McPherson Suture Tying Forcep Straight - With Tying Platform 10 cm - 4" Jaws Length 5 mm	1
FR-750-01	Castroviejo Suture Tying Forcep 1 x 2 Teeth with Tying Platform 11 cm - 4 1/4" Tip Size 0.12 mm	1
OP-260-12	Gill Corneal Knife Curved Cutting Edge - Straight 12.5 cm - 5"	1
SU-211-14	Castroviejo Micro Needle Holder Curved - Delicate 14 cm - 5 1/2"	1
SS-440-11	Iris Scissor Straight 11.5 cm - 4 1/2"	1
SS-430-10	Iris Scissor Straight 10 cm - 4"	1
FR-834-12	Bonaccolto Conjunctival Forcep Longitudinal Serrated Jaws with Cross Serrations at Tips 10.5 cm - 4" Tip Width 1.2 mm	1
FR-776-02	Lester Fixation Forcep Angled - 2 x 3 Teeth 9.5 cm - 3 3/4"	1
ST-110-03	Storing Case Large 265 x 160 x 25 mm Plastic	1

Scleral Buckling (Retinal Detachment) Set

Codes	Description	Quantity
OP-010-15	Barraquer Wire Speculum Standard Wire Blade Size 15 mm	1
DI-850-20	Castroviejo Caliper Straight 8 cm - 3" Measuring Range 20 mm	1
OP-160-14	Schochet Scleral Depressor Double End With Pocket Clip 14 cm - 5 1/2"	1
FR-760-08	Bishop-Harmon Iris Forcep Standard 1 x 2 Teeth 8.5 cm - 3 1/4"	1
FR-750-01	Castroviejo Suture Tying Forcep 1 x 2 Teeth with Tying Platform 11 cm - 4 1/4" Tip Size 0.12 mm	1
FR-700-01	Graefe Eye Dressing Forcep 10 cm - 4"	1
FR-781-10	McPherson Suture Tying Forcep Angled - With Tying Platform 10 cm - 4" Jaws Length 5 mm	1
FR-716-07	Colibri Corneal Forcep Very Delicate 1 x 2 Teeth 7.5 cm - 3 1/4" Tip Size 0.12 mm	1
OP-881-01	Jameson Muscle Hook Small 13 cm - 5"	1
SU-333-14	Barraquer Micro Needle Holder Curved - Very Delicate - Round Handle - With Lock 14 cm - 5 1/2"	1
SU-222-14	Castroviejo Micro Needle Holder Straight - Standard - With Lock 14 cm - 5 1/2"	1
DI-880-15	Ruler Graduated in mm and inches 16.5 cm - 6 1/2" Measuring Range 150 mm	1
SS-441-11	Iris Scissor Curved 11.5 cm - 4 1/2"	1
SS-453-11	Stevens Tenotomy Scissor Curved - Blunt/Blunt 11 cm - 4 1/2"	1
SS-683-11	Westcott Stitch Scissor Curved - Very Sharp Pointed Tips - Standard Blades 11.5 cm - 4 1/2"	1
AF-062-03	Dieffenbach Bulldog Clamp Straight 3.5 cm - 1 1/2"	2
OP-750-01	Lewis Lens Loop Small 13 cm - 5"	1
FR-750-03	Castroviejo Suture Tying Forcep 1 x 2 Teeth with Tying Platform 11 cm - 4 1/4" Tip Size 0.3 mm	1
FR-750-05	Castroviejo Suture Tying Forcep 1 x 2 Teeth with Tying Platform 11 cm - 4 1/4" Tip Size 0.5 mm	1
AF-160-12	Halsted-Mosquito Haemostatic Forcep Straight 12.5 cm - 5"	2
FR-903-01	Jameson Muscle Forcep Right - With Slide Lock 6 Teeth - Adult Size 10 cm - 4"	1
FR-903-02	Jameson Muscle Forcep Left - With Slide Lock 6 Teeth - Adult Size 10 cm - 4"	1
FR-784-10	Kelman-McPherson Suture Tying Forcep Angled - With Tying Platform 10 cm - 4" Jaws Length 7.5 mm	1
FR-702-01	Graefe Dissecting Forcep 1 x 2 Teeth 10 cm - 4"	1
FR-840-10	Nugent Utility Forcep Smooth Jaws 10.5 cm - 4"	1
AF-755-09	Schaedel Towel Clamp Sharp 9 cm - 3 1/2"	2
OP-875-13	Gass Retinal Detachment Hook With Oval Hole 13 cm - 5"	1
OP-883-13	Green Muscle Hook 13 cm - 5"	1
OP-849-01	Twist Hook Right 11 cm - 4 1/4"	1
SU-333-14	Barraquer Micro Needle Holder Curved - Very Delicate - Round Handle - With Lock 14 cm - 5 1/2"	1
OP-171-15	Schepens Orbital Retractor 15 cm - 6"	1
SS-643-11	Wecker Iris Scissor Angled - Sharp/Blunt 11 cm - 4 1/2"	1
SS-440-11	Iris Scissor Straight 11.5 cm - 4 1/2"	1
OP-230-02	Disposable Crescent Knife Angled Box of 6	1
OP-232-15	Disposable Stab Knife Straight 15 Degree Cutting Edge Box of 6	1
ST-112-26	Storing Case Double Level 265 x 160 x 50 mm Plastic	1

Blepharoplasty Set (Eyelid Surgery Set)

Codes	Description	Quantity
DI-880-15	Ruler Graduated in mm and inches 16.5 cm - 6 1/2" Measuring Range 150 mm	1
SL-012-03	Scalpel Handle No. 3 Solid, cm-Graduation 12 cm - 4 3/4"	1
SS-440-11 SC	SuperEdge™ Iris Scissor Straight 11.5 cm - 4 1/2"	1
SS-441-11 SC	SuperEdge™ Iris Scissor Curved 11.5 cm - 4 1/2"	1
AF-160-12	Halsted-Mosquito Haemostatic Forcep Straight 12.5 cm - 5"	2
AF-161-12	Halsted-Mosquito Haemostatic Forcep Curved 12 cm - 4 3/4"	2
SU-222-14	Castroviejo Micro Needle Holder Straight - Standard - With Lock 14 cm - 5 1/2"	2
SS-243-12	Kilner Dissecting Scissor Curved 12 cm - 4 3/4"	1
RT-074-02	Joseph Mucosa Hook 16.5 cm - 6 1/2" Width 2.0 mm	2
AF-773-08	Backhaus Towel Clamp 9 cm - 3 1/2"	4
FR-250-12	Adson Dissecting Forcep 1 x 2 Teeth 12 cm - 4 3/4"	1
RT-120-14	Desmarres Lid Retractor 16 cm - 6 1/4" Blade Width 14.0 mm	1
DI-850-20	Castroviejo Caliper Straight 8 cm - 3" Measuring Range 20 mm	1
FR-760-08	Bishop-Harmon Iris Forcep Standard 1 x 2 Teeth 8.5 cm - 3 1/4"	1
ST-110-02	Storing Case Medium 200 x 110 x 25 mm Plastic	1

Oculoplasty Set (Eyelid Surgery Set)

Codes	Description	Quantity
OP-078-00	Lester-Burch Eye Speculum Adjustable	1
OP-171-15	Schepens Orbital Retractor 15 cm - 6"	1
OP-175-15	Desmarres Lid Retractor Thin Solid Blades - Size 2 13 cm - 5" Blade Width 15 mm	1
OP-180-02	Scleral Shield Large (Pair) Blade Width 22 x 30 mm	1
OP-192-11	Jaeger Lid Plate 11 cm - 4 1/4"	1
OP-197-14	Wright Fascia Needle Oval Hole 14 cm - 5 1/2" Width - Length 1 x 6 mm	1
RT-068-13	Knapp Fine Wound Retractor 4 Blunt Prongs 13.5 cm - 5 1/4"	1
OP-550-03	Meyerhoefer Chalazion Curette Size 3 12.5 cm - 5" Diameter 2.5 mm	1
OP-800-02	Bunge Evisceration Spoon Small 14 cm - 5 1/2"	1
OP-845-02	Fixation Hook Large (Sharp) 13 cm - 5"	1
RT-048-02	Graefe Nerve Hook Fig. 2 14 cm - 5 1/2"	1
DI-880-15	Ruler Graduated in mm and inches 16.5 cm - 6 1/2" Measuring Range 150 mm	1
SS-452-11	Stevens Tenotomy Scissor Straight - Blunt/Blunt 11 cm - 4 1/2"	1
SS-469-11	Enucleation Scissor Curved 11.5 cm - 4 1/2"	1
FR-702-01	Graefe Dissecting Forcep 1 x 2 Teeth 10 cm - 4"	1
FR-906-10	Putterman Mullers Muscle Conjunctival Resection Ptosis Clamp - With Side Lock 10.5 cm - 4"	1
FR-916-09	Francis Chalazion Forcep Medium - Oval Jaws 9.5 cm - 3 3/4" Jaw Size 14 x 11 mm	1
FR-910-09	Downes Lid Clamp 9.5 cm - 3 3/4"	1
AF-140-10	Hartmann Haemostatic Forcep Straight 10 cm - 4"	1
SU-222-14	Castroviejo Micro Needle Holder Straight - Standard - With Lock 14 cm - 5 1/2"	1
RH-475-18	Freer Septum Elevator Double Ended - Sharp/Blunt 18.5 cm - 7 1/4" Blade Size 1 - Blade Size 2 5.0 mm - 5.0 mm	1
RH-604-16	West Septum Chisel 16 cm - 6 1/4"	1
OP-948-02	Wilder Lacrimal Dilator Medium Taper - Size 2 10.5 cm - 4 1/4"	1
PR-204-34	Bowman Probe Fig. 3/4 13 cm - 5" Tip Diameter 1.3 mm - 1.4 mm	1
OP-954-14	Worst Pigtail Probe With Hook 16 cm - 6 1/4"	1
SL-012-03	Scalpel Handle No. 3 Solid, cm-Graduation 12 cm - 4 3/4"	1
SL-100-11	Scalpel Blade No. 11 Pack of 100	1
SL-100-13	Scalpel Blade No. 15 Pack of 100	1
ST-112-26	Storing Case Double Level 265 x 160 x 50 mm Plastic	1

23 G Vitreoretinal Set

Codes	Description	Quantity
OP-091-16	Castroviejo Eye Speculum Adjustable Mechanism Blade Size 16 mm	1
AF-062-03	Dieffenbach Bulldog Clamp Straight 3.5 cm - 1 1/2"	4
AF-160-12	Halsted-Mosquito Haemostatic Forcep Straight 12.5 cm - 5"	2
SS-440-11	Iris Scissor Straight 11.5 cm - 4 1/2"	1
OP-471-19	Bishop-Harmon Anterior Chamber Irrigator Gauge 19	1
SS-452-11	Stevens Tenotomy Scissor Straight - Blunt/Blunt 11 cm - 4 1/2"	1
SS-453-11	Stevens Tenotomy Scissor Curved - Blunt/Blunt 11 cm - 4 1/2"	1
SS-675-11	Westcott Tenotomy Scissor Curved - Blunt Tips - Medium Blades 11 cm - 4 1/2"	1
SS-683-11	Westcott Stitch Scissor Curved - Very Sharp Pointed Tips - Standard Blades 11.5 cm - 4 1/2"	1
SL-206-13	Scalpel Handle For Micro Scalpel Blades 13 cm - 5"	1
SL-210-64	Micro Scalpel Blade No. 64 Pack of 25	1
OP-234-19	Disposable MVR Knife Straight 19 Ga Box of 6	1
OP-171-15	Schepens Orbital Retractor 15 cm - 6"	1
RT-048-02	Graefe Nerve Hook Fig. 2 14 cm - 5 1/2"	1
OP-875-13	Gass Retinal Detachment Hook With Oval Hole 13 cm - 5"	1
OP-160-14	Schocket Scleral Depressor Double End With Pocket Clip 14 cm - 5 1/2"	1
DI-850-20	Castroviejo Caliper Straight 8 cm - 3" Measuring Range 20 mm	1
OP-151-01	Gass Scleral Marker Tip Diameter 1.5 mm	1
FR-907-11	Watzke Sleeve Spreading Forcep 11.5 cm - 4 1/2" Jaws Length 10 mm	1
FR-834-12	Bonaccolto Conjunctival Forcep Longitudinal Serrated Jaws with Cross Serrations at Tips 10.5 cm - 4" Tip Width 1.2 mm	1
FR-840-10	Nugent Utility Forcep Smooth Jaws 10.5 cm - 4"	1
FR-750-01	Castroviejo Suture Tying Forcep 1 x 2 Teeth with Tying Platform 11 cm - 4 1/4" Tip Size 0.12 mm	1
FR-750-03	Castroviejo Suture Tying Forcep 1 x 2 Teeth with Tying Platform 11 cm - 4 1/4" Tip Size 0.3 mm	1
FR-760-08	Bishop-Harmon Iris Forcep Standard 1 x 2 Teeth 8.5 cm - 3 1/4"	2
SU-213-14	Castroviejo Micro Needle Holder Curved - Delicate - With Lock 14 cm - 5 1/2"	1
SU-212-14	Castroviejo Micro Needle Holder Straight - Delicate - With Lock 14 cm - 5 1/2"	1
FR-780-10	McPherson Suture Tying Forcep Straight - With Tying Platform 10 cm - 4" Jaws Length 5 mm	1
OP-211-07	Helveston Tissue Retractor Thin Curved Blade 12.5 cm - 5" Blade Width 7 mm	1
OP-211-09	Helveston Tissue Retractor Thin Curved Blade 12.5 cm - 5" Blade Width 9 mm	1
OP-429-22	Viscoelastic Aspirating Cannula Angled at 10 mm Gauge 22	1
OP-427-19	Randolph Cyclodialysis Cannula Flat Tip Gauge 19	1
OP-934-35	Braunstein Caliper 3.5 / 4.0 mm	1
SS-721-23	Virectomy Scissor Curved - Horizontal Opening 23 G	1
FR-680-23	Virectomy Forcep Straight Horizontal Opening - Serrated Tips 23 G	1
FR-682-23	Virectomy Forcep Straight End Gripping Jaws Horizontal Opening Tips 23 G	1
FR-689-23	Virectomy Forcep Straight Asymmetrical Jaws 23 G	1
FR-690-23	ILM Virectomy Forcep Straight 23 G	1
ST-112-26	Storing Case Double Level 265 x 160 x 50 mm Plastic	1

Micro Incision Cataract Surgery (MICS) Surgery Set

Codes	Description	Quantity
OP-357-18	MICS Diamond Knife Angled Trapezoid Cutting Edge 12 cm - 4 3/4" Blade Width 1.8 to 2.2 mm Titanium Handle	1
OP-735-12	Alio-Kelman IOL Hook 12 cm - 4 3/4" Tip Diameter 0.2 mm	1
SS-725-23	Alio MICS Virectomy Scissor Straight Horizontal Opening With Curved Shaft 23 G	1
FR-695-23	Alio MICS Virectomy Capsulorrhesis Forcep Angled Shaft 23 G	1
FR-807-11	Inamura Capsulorrhesis Forcep Delicate Tapered Shanks 11.5 cm - 4 1/2" Shanks Length 12 mm	1
FR-863-11	Ernest Nucleus Cracker Extra Delicate Cross Action Serrated Paddle Tips 11 cm - 4 1/4"	1
OP-549-01 T	Irrigation Handpiece For Bimanual Technique 2 x 0.5 mm (21G) Side Ports 10.5 cm - 4" Titanium	1
OP-549-02 T	Aspiration Handpiece For Bimanual Technique 2 x 0.5 mm (21G) Side Ports 10.5 cm - 4" Titanium	1
ST-110-02	Storage Case Medium 200 x 110 x 25 mm Plastic	1

Deep Anterior Lamellar Keratoplasty (DALK) Set

Codes	Description	Quantity
OP-083-15	Lieberman Eye Speculum Adjustable Blade Size 15 mm	1
OP-629-01	DALK Corneal Dissector 12 mm Long Flat Vaulted Spatula With Vaulted Tip 12.5 cm - 5"	1
OP-629-02	Rosenwasser DALK Trisector 10 mm Long Triangular Shaft With Spatula Shaped Tip 12.5 cm - 5"	1
OP-629-03	Bonfadini DALK Spatula Vaulted 1.0 mm Wide 8 mm Spatula Shaped Blade With Central Groove 12.5 cm - 5"	1
SS-649-01	Holland DALK Scissor With Inner Spatulated Blade Curved Right 10.5 cm - 4"	1
SS-649-02	Holland DALK Scissor With Inner Spatulated Blade Curved Left 10.5 cm - 4"	1
SS-705-08	Holland DALK Scissor Angled Blade With Blunt Tips 8 cm - 3 1/4"	1
OP-417-27	DALK Air Injection Cannula Vaulted Shank With Flattened Tip Gauge 27	1
OP-419-27	Rosenwasser DALK Air Injection Cannula Flattened Cobra Shaped Tip Gauge 27	1
FR-816-10	Tennant Suture Tying Forcep Straight - Round Handle with Guide Pin - Extra Delicate Smooth Jaws 10.5 cm - 4" Jaws Length 6 mm	1
FR-817-10	Tennant Suture Tying Forcep Curved - Round Handle with Guide Pin - Extra Delicate Smooth Jaws 10.5 cm - 4" Jaws Length 6 mm	1
SU-313-11	Barraquer Micro Needle Holder Curved - Delicate - Round Handle 11.5 cm - 4 1/2"	1

General Terms and Conditions of Sales

1. Definitions

“Company” means Surtex Instruments Ltd (registered in the United Kingdom under number 7035724).

“Contract” means the Order accepted by the Company and incorporating these conditions.

“Customer” means the company, firm or organisation who’s Order for the Goods is accepted by the Company.

“Goods” means the goods which the Company agrees to supply in accordance with these conditions and which are identified on the Order. “Order” means the order for the purchase of the Goods submitted by the Customer.

“Goods” means the goods agreed to be provided to the customs by Surtex Instrument Ltd.

“Quotation” means the quotation by Surtex Instruments stating the price at which and the time period during which it will supply the goods to the customer.

1. Any reference to a law is a reference to it as in force for the time being taking account of any amendment, extension, application or re-enactment and includes any subordinate legislation for the time being in force made under it.

2. Where a condition provides for a list of items following the word “including” or “includes”, then such list is not to be interpreted as being an exhaustive list.

2. General

2.1 Headings are purely for the ease of reference and do not form part of effect the interpretation of Conditions.

2.2 In the event of a conflict the contractual documents shall be applied in the following order:

- a) the Quotation
- b) the Conditions

2.3 No other terms and conditions or warrantied are to be incorporated into the Contract unless agreed expressly in writing with Surtex Instruments. In particular, the Customer’s terms and conditions are expressly excluded whether provided in full, referred to in other documents and whether before or after the Contract is formed.

2.4 Surtex Instruments may update or amend these Conditions from time to time to comply with the law or to meet the changing business requirements without notice to the customer.

3. Ordering

3.1 Each Order (whether submitted by telephone, fax, post, email or electronic ordering system) shall be deemed to be an offer by the Customer to purchase the Goods pursuant to these conditions. The Company may, at its complete discretion, refuse to accept any Order. The Contract will be formed upon confirmation of acceptance of the Customer’s purchase order by Surtex Instruments.

3.2 The Customer shall be responsible to the Company for ensuring the accuracy of the terms of any Order (including the description of the Goods ordered).

Samples and Free of Charge Goods

3.3 For the avoidance of doubt, where no order is placed and the Company provides Samples, to familiarize the Customer with the goods, as part of a sampling program or where the Company provides Free of Charge Goods to the Customer as part of a clinical trial or compassionate use program these conditions shall, to the extent applicable, apply to such Samples and Free of Charge Goods.

4. Prices

4.1 The price is exclusive of VAT and all other applicable taxes and duties, which the Customer shall be additionally liable to pay the Surtex Instruments.

4.2 All prices quoted are subject to change up to the date of delivery, subject to prior written notice by Surtex Instruments.

4.3 Surtex Instruments reserves the right, by giving notice to the Customer at any time before delivery, to increase the price of the Goods to reflect any increase in the cost of Surtex Instruments which is due to any factor beyond the reasonable control of Surtex Instruments.

4.4 Errors on Surtex Instruments website prices and availability of products and services are subject to change without notice errors will be corrected where discovered, and Surtex Instruments reserves the right to revoke any stated offer and to correct any errors, inaccuracies or omissions including after an order has been submitted and whether or not the order has been confirmed and your credit card charged. Surtex Instruments shall have the right to refuse or cancel any orders placed for product listed at the incorrect price. Surtex Instruments shall have the right to refuse or cancel any such orders whether or not the order has been confirmed and your credit card charged. If your credit card has already been charged for the purchase and your order is cancelled, Surtex Instruments will issue a credit to your credit card account in the amount of the charge. Individual bank policies will dictate when this amount is credited to your account. If you are not fully satisfied with your purchase, you may return it in accordance with Surtex Instruments’ Return Policy.

5. Delivery

5.1 Delivery of the Goods shall be made to the delivery address specified as stated on the Pro Forma Invoice or communicated separately whichever is greater after the date of the Surtex Instruments acceptance of the Order.

5.2 Surtex Instruments reserves the right to withhold delivery of any Goods if any sum due under any contract with the Customer is overdue or if, in the reasonable opinion of Surtex Instruments, the financial standing of the Customer has been impaired for any other reason.

5.3 Time for delivery shall not be of the essence and, Surtex Instruments shall not be liable for any loss or damage whatever due to any failure by the Surtex Instruments to deliver the Goods (or any of them) promptly or at all. The customer will have no right damages or to cancel the order for any failure to meet the delivery times stated.

5.4 The delivery of the Goods may happen by separate instalments. Each separate instalment shall be invoiced and paid for in accordance with the provisions of the Contract. Each instalment shall be a separate Contract and no cancellation or termination of any one Contract relating to an instalment shall entitle the Customer to repudiate or cancel any other Contractor instalment.

5.5 The Customer shall take delivery of the Goods promptly whenever they are tendered for delivery (which shall include signing the relevant delivery note) If, for any reason, the Customer fails to take such delivery of the Goods or otherwise causes or requests a delay in delivery:

a) The Goods shall be deemed to have been delivered;

b) Risk in the Goods shall pass to the Customer, and

c) Surtex Instruments shall be entitled to store or arrange for the storage of the Goods until actual delivery and charge the Customer for the reasonable costs (including insurance) of such storage and any additional charges incurred in delivering the Goods at a later date.

5.6 The Customer shall ensure that adequate and safe facilities and procedures exist at the Customer’s premises for the delivery of the Goods.

5.7 The Customer shall immediately notify Surtex Instruments Customer Support Centre if the Goods are not received within 10 days of the specified delivery date.

5.8 The customer outside of the UK shall be solely responsible for obtaining all import authorisations.

6. Title and Risk

6.1 The Goods shall be at the Customer’s risk at the time of delivery unless the parties have agreed in writing that the Goods are to be collected from the Surtex Instruments premises, in which case, the Goods shall be at the Customer’s risk when Surtex Instruments notifies the Customer that the Goods are available for collection.

1. Property in the Goods shall not pass from Surtex Instruments until:

a) The Customer has paid the price in full; and

b) No other sums whatsoever are due from the Customer to Surtex Instruments.

2. Until property in the Goods passes to the Customer in accordance with the condition 6.1, the Customer shall:
 - a) Hold the Goods on a fiduciary basis as bailed for Surtex Instruments
 - b) Store the Goods separately from all other goods in its possession and mark them in such a way that they are clearly identified as Surtex Instrument's property; and
 - c) Maintain the Goods in satisfactory condition and keep them insured on Surtex Instrument's behalf for their full price.
3. The Customer may sell at full market value or use the Goods only in the ordinary course of the Customer's business. Any such sale or dealing shall be a sale or use of Surtex Instrument's property on the Customer's own behalf and the Customer shall deal as principal when making such sales or dealings. Until property in the Goods passes from Surtex Instruments, the proceeds of sale or otherwise of the Goods, to the full extent of the price of the Goods, shall be held in trust for Surtex Instruments and shall not be mixed with other money or paid into any overdrawn bank account and shall be at all material times identified as Surtex Instruments money.
4. Surtex Instruments shall be entitled to recover the price notwithstanding that property in any of the Goods has not passed from Surtex Instruments.
5. Until such time as the property in the Goods passes from Surtex Instruments, if the Customer: a) Enters into liquidation or has a winding-up order made against it or has an administrator, receiver or manager appointed in respect of its assets; or b) Fails to make payment in accordance with these conditions, the Customer shall upon request deliver up Surtex Instruments such of the Goods as have not ceased to be in existence or resold. If the Customer fails to do so, Surtex Instruments may enter upon any premises owned. Occupied or controlled by the Customer where the Goods are situated and repossess the Goods.
6. The Customer shall not pledge or in any way charge by way of security for any indebtedness, any of the Goods which are the property Surtex Instruments. If the Customer does so, all sums whatever owing by the Customer to Surtex Instruments shall forthwith become due and payable.
7. On termination of the Contract, howsoever caused, Surtex Instruments rights contained in this condition 6 shall remain in effect.

7. Terms of Payment

- 7.1 Payment is due to Surtex Instruments within 30 days of invoice (the "Due Date"). In the case of payment other than by cash, payment shall only be deemed to have been made on receipt of cleared funds.
- 7.2 Each consignment or partial delivery shall be invoiced and paid separately.
- 7.3 No disputes arising under the contract no delays in the delivery period beyond the control of Surtex Instruments shall interfere with prompt payment in full by Customer. The Customer shall not in any circumstances or for any reason whatsoever be entitled to make any deduction or withhold any sum from the price by way of set-off or otherwise.
- 7.4 If the Customer fails to make the payment on the due date, then Surtex Instruments without prejudice to other rights and remedies shall be entitled to charge the Customer interest (both before and after judgment) on the amount unpaid at the rate of 4% per annum above the basic rate of HSBC plc during the delay and suspend all further deliveries without notice.
- 7.5 All payments payable to the Company under the Contract shall become due immediately on its termination, despite any other provision.

8. Warranties

- 8.1 Surtex Instruments will use reasonable endeavours to ensure that the Goods
 - a) Correspond to their description and any applicable specification;
 - b) Be of satisfactory quality and fit for purpose held out by Surtex Instruments or made to known to the customer by Surtex Instruments
 - c) Where applicable, be free from defects in design, materials and workmanship and remain so for a period of 12 months after delivery; and
 - d) Comply with all applicable statutory and regulatory requirements relating to the manufacture, labelling, packaging, storage, handling and delivery of Goods
- 8.3 In case of Goods not manufactured by Surtex Instruments
 - a) Surtex Instruments gives no assurance or warranty whatsoever that the sale or use of the Goods will not infringe the patent design right registered, design copyright or other industrial property rights of other person firm or Company and
 - b) The warranty given to the Customer shall be limited to the warranty (if any) which Surtex Instruments receives from the manufacturer or supplier of the Goods to Surtex Instruments.

9. Missing and Defective Goods

- 9.1 The Customer must inspect the Goods immediately on arrival at their premises.
- 9.2 In case of missing goods or defects apparent on inspection:
 - a) a written complaint must be made to Surtex Instruments by the customer within seven days of receipt of the Goods specifying the shortage and/or defect and
 - b) Surtex Instruments shall be given an opportunity to inspect the Goods and investigate any complaint before any use is made of the Goods.
- 9.3 In the case of defects not apparent on inspection:
 - a) the Customer shall send a written complaint to Surtex Instruments as soon as reasonably practicable after the defect is noticed and ensure no use is made of or alterations made to the Goods thereafter;
 - b) the Customer shall allow Surtex Instruments 30 days to inspect the Goods and shall take such steps as are necessary to enable Surtex Instruments to do so including returning of the Goods to its head office in the United Kingdom if requested.
- 9.4 In any potential case of a product being defective pursuant to Part 1 of the Consumer Protection Act 1987, the customer must send the complaint within 12 months of the date of delivery by Surtex Instruments or in the case of items not manufactured by Surtex Instruments, within the guarantee period specified by the manufacturer of such item.
- 9.5 Any complaints made outside of the periods specified in clause 9 will not be considered.
- 9.6 All allegations of total non-delivery of any consignment of the Goods must be made by notice in writing by the customer to Surtex Instruments within 10 days of the specified delivery date. Failure by the Customer to give such notice shall render the Customer liable for any temporary or permanent loss of the Goods and all additional costs and expenses incurred by Surtex Instruments in relation to such loss.

10. Liability

- 10.1 Save where Goods manufactured by Surtex Instruments are held to be defective as defined in Part 1 of the Consumer Protection Act 1987, Surtex Instruments sole obligation in the event of any missing or defective Goods, shall be to repair or replace the Goods.
- 10.2 Surtex Instruments shall not be liable in respect of consequential loss, loss of profits, damage to the property, loss of goodwill, loss of business opportunity, wasted expenditure, cost of mitigation arising out of or occasioned by any fault or defect in the Goods supplied or any other losses howsoever caused.
- 10.3 The maximum total liability under or arising out of in connection with this Agreement will not exceed the total of the Contract value.
- 10.4 Nothing in this Agreement shall exclude or limit either party's liability for death or personal injury resulting from negligence or their fraudulent acts, omissions or statements or for infringement of third party's intellectual property rights.
- 10.5 For the purpose of the Health and Safety at Work Act 1974, Surtex Instruments relies upon the testing examination and research carried out by the suppliers to it of propriety equipment and component parts of Goods.
- 10.6 This clause 10 shall survive the termination of the contract.

11. Return of Goods

- 11.1 In cases where Goods are not defective and the customer requests to return the Goods, Surtex Instruments shall be under no contractual obligation to accept the return except in the event of any error on its parts as to the amount or type of Goods delivered.

- 11.2 Any requests to return must be made to Surtex Customer Support Centre within ten working days of the date of actual delivery.
- 11.3 The return of Goods to correct an error in delivery or ordering will only be considered if the Customer Care Team is notified immediately or in any event within 48 hours of delivery to the Customer. Goods must be returned to Surtex Instruments within 5 days of delivery to the customer to correct an error in delivery or ordering.
- 11.4 For the avoidance of doubt, unless defective, the returns of Goods as detailed in sections 11.2, 11.4, 11.5 and 11.6 must be unused, undamaged and in the original packaging. The Goods must be the same batch as originally supplied by Surtex Instruments and are not falsified/counterfeit. All goods are individually assessed by Surtex Instruments before being authorised to return.
- 11.5 Return of Goods for any reason shall require prior consent from Surtex Instruments which shall be sent solely at the discretion of Surtex Instruments. To obtain consent, the Customer Care Team must be notified. Consideration to requests for a return may be given if:
- a) of a bona fide nature;
 - b) made within 10 working days of the actual delivery; and
 - c) unless defective, the Goods must be unused, undamaged and in the original undamaged packaging. The Goods must be the same batch as originally supplied by Surtex Instruments and not falsified/counterfeit.
- 11.6 On the customer contacting Surtex Instruments to request a return, the Customer will be issued with a returns merchandise authorisation number by Surtex Instruments. Goods returned without the prior written approval of Surtex Instruments may be returned to the customer or disposed of at the absolute discretion of the company.
- 11.7 Goods returned will be subjected to the following charges:
- a) for palletised items there will be a standard charge of £75.00 per pallet for any return uplifted by Surtex Instruments.
 - b) for returns of small parcel there will be a minimum £15.00 charge per parcel and a 30% re-stocking charge.
 - c) the event of return resulting from an error on the part of Surtex Instruments in relation to the amount or type of Goods delivered, no charges under this clause 11.7 will apply.
- 11.8 In the event, that any Goods are subject to be called by Surtex Instruments for the purpose of recall, the customer shall provide to Surtex Instruments as much information as it is reasonably possible regarding the whereabouts of the Goods.

12. Cancellation Policy

- 12.1 In the event of cancellation before delivery and if returned after delivery where products are not defective, there will be a re-stocking charge of not less than 30% of the invoice value charged by Surtex Instruments.
- 12.2 Cancellation of an order will only be agreed by Surtex Instruments at its sole discretion, on condition that all costs and expenses incurred by Surtex Instruments up to the time of cancellation and loss of profits and other loss or damage resulting to Surtex Instruments of such cancellation will be reimbursed by the Customer to Surtex Instruments forthwith.

13. Specifications

Surtex Instruments reserve the right to alter or change the specification and dimension of the Goods supplied with reasonable limits having regard to the nature of the Goods. These do not materially affect the quality or fitness for purpose of Goods. Dimensions specified by Surtex Instruments are to be treated as an approximate only unless the Customer specifically states in writing that exact dimensions are required. Our website and catalogues display the product images as accurately as possible, However, Surtex Instruments cannot guarantee the finish/colour you see matches the actual product as the display of the finish and/or colour depends in part, upon your monitor and its settings.

14. Marking and Instructions

- 14.1 No name, mark, numbering, colouring, appearance or log on the goods will be obscured, removed or concealed by the Customer.
- 14.2 The Customer will not repackage or alter the presentation of Goods and it will not assist, cause, or enable any other party to do any of the said acts or deal in Goods or packaging which have been subject to such acts.
- 14.3 The Customer acknowledges that it is under a duty to pass on to its customers (where appropriate) all instructions, information and warnings supplied by Surtex Instruments with the Goods.

15. Confidential Information

- 15.1 A party shall treat the Contract as private and confidential and neither party shall publish or disclose any particulars without the prior written consent of the other with the exceptions of 15.2 below.
- 15.2 The obligations expressed in 15.1 above shall not apply to any information which:
- a) Is or subsequently comes into the public domain other than by breach of this clause;
 - b) The disclosure of which is required to be disclosed by law, any governmental or regulatory authority or by a court of competent jurisdiction;
 - c) Comes into the possession of one of the parties prior to its disclosure by the other, or which is acquired lawfully and in good faith from an independent third party.
- 15.3 On the conclusion or termination of the Contract, both parties shall destroy all copies of confidential information obtained from the other and the obligation relating to the confidentiality shall continue for a period of 10 years.
- 15.4 The Customer shall not, without the prior written consent of Surtex Instruments, advertise or make known to any third that Surtex Instruments provided Goods to the Customer.

16. Insolvency

- 16.1 Without limiting its other rights or remedies, Surtex Instruments shall be entitled to cancel the contract, in whole or in part, by notice in writing if:
- a) The Customer makes a voluntary arrangement with its creditors or (being an individual firm) becomes bankrupt or (being a company) becomes subject to an administration order or goes into liquidation (otherwise than for purpose of amalgamation or reconstruction);
 - b) An encumbrancer takes possession, or a receiver is appointed, of any of the property or assets of the Customer.
 - c) The Customer ceases, or threatens to cease, to carry on business; d) Surtex Instruments reasonably apprehends that any of the events mentioned above is about to occur in relation to the Customer and notifies the Customer accordingly.

17. Sub-Contractors

- 17.1 Surtex Instruments may assign the Contractor any part of it to any company, firm or organization.
- 17.2 The Customer shall not be entitled to assign the Contractor any part of it without the prior written consent of Surtex Instruments (such consent not to be unreasonably withheld or delayed).

18. Force Majeure

- 18.1 Surtex Instruments shall not be responsible for any failure to fulfil its obligations under these conditions and/or in relation to the Goods if such failure is caused by an event which is beyond its reasonable control including acts of God, war, armed conflict, acts of terrorism, riot, strikes, industrial action (including action involving Surtex Instruments employees), fire, flood, any natural disaster, any defaults by Surtex Instruments suppliers or subcontractors, difficulties obtaining workmen or materials, and breakdown of machinery.
- 18.2 If such event continues for more than one month, Surtex Instruments may terminate the Contract without liability.

19. Waiver

19.1 waiver of any right or remedy under the contract or law is only effective if given in writing and shall not be deemed a waiver of any subsequent breach or default. No failure or delay by a party to exercise any right or remedy provided under the Contract or by law shall constitute a waiver of that or any other right or remedy. No single or partial exercise of such right or remedy shall prevent or restrict the further exercise of that or any other right or remedy.

20. Severance

20.1 any provision or part-provision of the Contract is or becomes invalid, illegal or unenforceable, it shall be deemed modified to the minimum extent necessary to make it valid, legal and enforceable whilst retaining the risk and reward between the parties. If such modification is not possible, the relevant provision or part-provision under this clause shall not affect the validity and enforceability of the rest of the Contract.

21. Variation

21.1 as set out in these Conditions, no variation of the Contract, including the introduction for any additional terms and conditions, shall become effective unless it is agreed in writing and signed by Surtex Instruments.

22. General

22.1 rights and remedies of Surtex Instruments set out in these conditions shall be in addition and without prejudice to any other rights and remedies which may be available to Surtex Instruments at common law or under the statute.

23. Governing Law and Jurisdiction

23.1 The contract and any dispute or claim arising out of or in connection with it or its subject matter or formation (including non-contractual disputes or claims) shall be governed by and construed in accordance with the law of England and Wales.

23.2 Each party irrevocably agrees that the courts of England and Wales shall have exclusive jurisdiction to settle any dispute or claim arising out of or in connection with his agreement or its subject matter or formation (including non-contractual disputes or claims).

23.3 The schedules to the Uniform Law of International Sales Act 1967 shall not in any circumstances apply to the contract.

24. Entire Agreement

24.1 This contract constitutes the whole agreement and understanding between the parties in relation to its subject matter and supersedes any previous arrangement, understanding, representations, (including any information or data set out in any sales and marketing literature), or agreement between them relating to the subject matter of this Agreement.

24.2 Each party acknowledges that, in entering into this Agreement, it has not relied on any statement, assurance or warranty (whether made negligently or innocently) other than those expressly set out in this Agreement and has no rights and remedies in respect thereof.

25. Third-Party Rights

25.1 term of this agreement shall be enforceable under the Contracts (Rights of Third Parties) Act 1999 by any legal entity who is not a party to this Agreement.

26. Notices

26.1 communications between the parties about the Contract shall be in writing and delivered by pre-paid first-class post or by hand to the party's registered office, principal place of business or such other address as may at the relevant time have been notified pursuant to this provision to the party giving notice. In the absence of evidence of earlier receipt, any notices shall be deemed to be duly served:

- a) If sent by first class post, two working days after posting; or
- b) If delivered personally, when left at the relevant address.

ALPHABETICAL INDEX



SURTEX
INSTRUMENTS

Crafted to Excellence

Alphabetical Index

A	Pages		B	Pages			
Adson			Binder	58			
Dressing Forcep	123		Flap Separator & Elevator				
Dissecting Forcep	124		Bishop-Harmon	32			
Agarwal			Anterior Chamber Irrigator				
Karate Phaco Chopper	50		Iris Forcep 137				
Agricola,			Dressing Forcep	138			
Lacrimal Sac Retractor	10		Blaydes				
AK			Lens Holding Forcep	156			
Diamond Knife	19		Blumenthal				
Akahoshi			Conjunctiva Dissector	40			
Nucleus Sustainer	47		Bonaccutto				
Phaco PreChopper	152		Conjunctival Forcep	149			
Aker			Bonfadini				
Lens Pusher	62		DALK Spatula	43			
Alfonso			Bonn				
Eye Speculum	2		Scleral Knife	11			
Nucleus Irrigating Spoon	37		Micro Iris Hook	74			
Nucleus Trisector	66		Corneal Forcep	133			
Cutting Platform	68		Iris Scissor	96			
Nucleus Grasping Forcep	152		Botvin				
Alio			Iris Forcep	130			
LASIK Irrigation Cannula	33		Bowman				
MICS Scissor	121		Lacrimal Probe	179			
MICS Capsulorrhesis Forcep	128		Bracken				
Alio-Kelman			Irrigation Cannula	33			
IOL Hook	62		Braunstein				
Anis			Caliper	90			
Nucleus Manipulator	56		Brierley				
Radial Marker	89		LASIK Cannula	35			
Arruga			Nucleus Splitter	51			
Micro Burr	17		Brown				
Capsule Forcep	151		Needle Holder	167			
Atkinson			Bucher-Kansas				
Retrobulbar Needle	30		Nucleus Trisector	66			
Ayer			Bunge				
Chalazion Forcep	162		Evisceration Spoon	72			
B							
Backhaus			Buratto				
Towel Clamp	164		LASIK Cannula	35			
Baird			LASIK Flap Protector	40-41			
Chalazion Forcep	160		LASIK Flap Forcep	143			
Banaji			IOL Inserting Forcep	155			
LASIK Cannula	34						
Bangerter			C				
Eye Speculum	3		Camellin				
Muscle Clamp	157		LASEK Trepbine	25			
Bard-Parker			Carter				
Scalpel Handle	93		LASIK Cannula	34			
Barkan			Sphere Introducer	71			
Goniotomy Knife	15		Castroviejo				
Barraquer			Eye Speculum	5			
Eye Speculum	1-2		Synechia Spatula	63			
Iris Spatula	45		Cyclodialysis Spatula	92			
Iris Scissor	107		Corneoscleral Scissor	108			
Cilia Forcep	126		Universal Corneal Scissor	111			
Micro Needle Holder	171-175		Keratoplasty Scissor	112			
Blade Holder/Breaker	94		Suture Tying Forcep	135-136			
Barraquer-Trotman			Capsule Forcep	151			
Micro Needle Holder	170		Micro Needle Holder	168-169, 176-177			
Barrett			Caliper	178			
Nucleus Expressor	63		Blade Holder/Breaker	94			
Bechert			Castroviejo-Colibri				
Nucleus Rotator	54		Corneal Forcep	132			
Lens Pusher	62		Chang				
Lens Holding Forcep	156		Hydrodissection Cannula	36			
Beer			Phaco Chopper	49-50			
Cilia Forcep	126		Charleux				
Berke			Cortex Extractor	36			
Forcep			Chavasse				
Ptosis	158		Strabismus Hook	86			

Alphabetical Index

C	Pages		Pages
Chayet			
Corneal Marker	86		
Christensen			
Punctum Punch	26		
Citelli			
Mastoid Rongeur	165		
Clayman			
Nucleus Rotator	54		
Lens Holding Forcep	157		
Clayman-Jaffe			
Lens Manipulating Hook	75		
Colibri			
Corneal Forcep	130-131		
Collin			
Meibomian Gland Expression Forcep	162		
Corbin			
Sub Tenon's Anesthesia Cannula	30		
Cottle			
Bone Mallet	179		
Crawford			
Retrieval Hook	92		
Lacrimal Intubation Set	92		
Crescent			
Diamond Knife	21		
Culler			
Iris Spatula	44		
Lens Spoon	69		
D			
DALK			
Air Injection Cannula	28		
Corneal Dissector	43		
Dastoor			
Gouge & Lacrimal Sac Dissector	11		
Pupil and Iris Repositor	44		
Erisophake (Bell)	72		
De La O			
Aspirating Cannula	36		
Deitz			
ICL Manipulator	80		
Desmarres			
Lid Retractor	8, 178		
Chalazion Forcep	161		
Entropian Forcep	162		
Dix			
Foreign Spud & Gouge	17		
DLEK			
Corneal Dissector	12		
Dodick			
Nucleus Cracker	152		
IOL Folding Forcep	153		
Dodick-Kammann			
Nucleus Chopper	51		
Double-X			
Eye Speculum	5		
Downes			
Lid Clamp	160		
Drews			
Irrigation Vectis	36		
Drof			
Infusion Cannula	37		
Drysdale			
Nucleus Manipulator	55		
E			
Elliot			
Scleral Trephine	25		
Ellis			
Foreign Body Spud	16		
Elschnig			
Cyclodialysis Spatula	62		
Fixation Forcep	139		
Enucleation			
Scissor	101		
Erhardt			
Lid Forcep	159		
Ernest			
Nucleus Cracker	152		
Escaf			
Nucleus Manipulator	46		
Nucleus Sustainer	47		
Nucleus Divider	52		
F			
Faulkner			
Lens Holding Forcep	155		
Feaster			
Wire Speculum	1		
Fechtner			
Trabeculectomy Marker	89		
Conjunctival Forcep	132		
Fenzl			
Lens Manipulating Hook	79		
Fine			
Micro Scissor	115		
Suture Tying Forcep	148		
Fine-Thornton			
Fixation Ring	82		
Flieringa			
Fixation Ring	6		
Foster			
Enucleation Snare	71		
Francis			
Foreign Body Spud	16		
Chalazion Forcep	161		
Frankel			
Eye Speculum	2		
Friedman			
Bone Rongeur	179		
Fukasaku			
Hockey Knife	14		
Sub Tenon's Anesthesia Cannula	30		
LASIK Cannula	34		
Hydro-Chop Cannula	36		
Snapper Hook	52		
Snapper Hook & Small Pupil	53		
LASIK Spatula	58		
Fukasaku-Lieberman			
Nucleus Manipulator and Hook	48		

Alphabetical Index

G	Pages	H	Pages
Gass		Hersh	58
Scleral Marker	7	IntraLase Flap Spatula	
Corneoscleral Punch	26	Heslin	38
Retinal Detachment Hook	84	Infusion Cannula	
Gerl		Hirschmann	
Suture Tying Forcep	148	Lens Spatula	60
Giannetti		Iris Hook	76
Nucleus Capture Spoon	68	Hoffer	
Nucleus Capture Spatula	68	Optic Zone Marker	88
Nucleus Capture Hook	69	Hoffer-McPherson	
Push Pull Manipulator	80	Tissue Forcep	135
Gills		Hofmann-Thornton	
Corneal Knife	14	Fixation Ring	83
Irrigating - Aspirating Cannula	31	Holland	
Aspirating Cannula	36	DALK Scissor	108, 117
Capsulotomy Scissor	119	Hoskin	
Gills-Colibri		Tissue Forcep	131
Corneal Forcep	130	Hunkeler	
Gills-Vannas		IOL Manipulator	61
Capsulotomy Scissor	116-118	Hwang	
Gimbel		Phaco Chopper	5
Cortex Aspirating Cannula	35	I	
Nucleus Spatula	45	ICL	
Stabilization Ring	82	Lens Cartridge Loading Forcep	128
Girard		ILM	
Synechia Spatula	57	Vitrectomy Scissor	127
Corneoscleral Forcep	137	Inamura	
Gold		Capsulorrhesis Forcep	145
Punctal Plug Forcep	150	Iris	
Goldman-Fox		Scissor	95-97, 102
Gum Scissor	95	Micro Scissor	114
Gonnin		J	
Amsler Marker	7	Jaeger	
Gradle		Lid Plate	9
Ribbon Scissor	103	Keratome	11
Cilia Forcep	149	Jaffe	
Graefe		Lid Retractor	8
Cataract Knife	13	Lens Spatula	60
Dressing Forcep	129	Lens Manipulating Hook	75
Dissecting Forcep	129	Stitch Scissor	114
Fixation Forcep	139	Capsulorrhesis Forcep	145
Nerve Hook	178	Suture Tying Forcep	146
Graether		Utility Forcep	151
Collar Button	73	Jaffe-Bechert	
Grene		Nucleus Rotator	54
Visual Axis Marker	87	Jaffe-Maltzman	
Guimaraes		Lens Manipulating Hook	80
LASIK Spatula	59	Jameson	
H	Pages	Muscle Hook	85
Halsted-Mosquito		Muscle Forcep	157-158
Haemostatic Forcep	163	Jansen	
Harms-Colibri		Capsule Polisher Cannula	32
Corneal Forcep	130	Jones	
Harris-Sinskey		Towel Clamp	164
Lens Manipulating Hook	75	Joseph	
Hartmann		Mucus Hook	178
Haemostatic Forcep	163		
Hartmann-Halle			
Nasal Speculum	166		
Hebra			
Chalazion Curette	39		
Helveston			
Tissue Retractor	9		
Teaser Hook	84		
Finder Hook	84		
Muscle Hook	85		
Scleral Ruler	89		
Suture Tying Forcep	148		

Alphabetical Index

K	Pages	L	Pages
Kansas		Lewis	
Lens Loop	65	Lens Loop	64
Nucleus Trisector	67	Lieberman	
Nucleus Bisector	67	Eye Speculum	4-5
Nucleus Vectis	67	Nucleus Manipulator	47
Katzen		Nucleus Manipulator and Spatula	48
LASIK Flap Spatula	59	Lippman	
Katzin		Nucleus Spatula	56
Corneal Transplant Scissor	110	Lister	
Kaye		Conjunctival Forcep	149
Ribbon Facelift Scissor	103	Utility Forcep	150
Kellen		Livernois	
Capsulorrhesis Marker	86	IOL Folding Forcep	153
Kelly		IOL Holding Forcep	154
Descemet's Membrane Punch	26	Livernois-McDonald	
Kelman		IOL Inserting Forcep	155
Cystotome	31	Lordan	
Kelman-McPherson		Chalazion Forcep	161
Tissue Forcep	134	LRI	
Suture Tying Forcep	143	Diamond Knife	18
Kerrison		Lu-Mendez	
Sphenoid Punch	165	Fixation Ring and Guide	82
Kilner		Luntz-Dodick	
Dissecting Scissor	95, 100	Trabeculectomy Punch	
Kirby			
Lens Loop and Muscle Hook	65		
Knapp			
Wound Retractor	178		
Knolle			
Anterior Chamber Irrigation Cannula	33		
		Maddox	
Nucleus Spatula	46	LASIK Spatula	59
Knolle-Pearce		Maloney	
Irrigation Vectis	37	Intra-Operative Keratometer	86
Koch		Maltzman-Fenzl	
Nucleus Spatula	46	Lens Manipulating Hook	79
Kraff		Manhattan	
Capsule Tag Forcep	144	Eye & Ear Fixation Forcep	139
Kratz		Mannis-Buratto	
Cystotome	32	LASIK Flap Protector	41
		Manson	
Lens Holding Forcep	156	Strabismus Hook	83
Kratz-Barraquer		Masket	
Eye Speculum	1	Phaco Spatula	49
Kuglen		Capsulorrhesis Forcep	145
Iris Hook and Lens Manipulator	77-78	McDonald	
Kuhnt		Optic Zone Marker	87
Corneal Scarifier	12	McIntyre	
Kurstin		Lacrimal Cannula	28
Flap Stretching Forcep	143	Cystotome	32
		Anterior Chamber Cannula	35
		McNeill-Goldman	
		Scleral Fixation Ring	6
L		McPherson	
		Suture Tying Forcep	134, 141-142
Lacrimal		Micro Needle Holder	167
Cannula	28	McPherson-Vannas	
Lambert		Capsulotomy Scissor	117
Chalazion Forcep	160	McPherson-Westcott	
Lancaster		Conjunctival Scissor	112
Eye Speculum	4	Stitch Scissor	113
LASEK		Mead	
Alcohol Cannula	30	Bone Rongeur	179
LASIK		Mellinger	
Flap Manipulator	60	Eye Speculum	3
Lester		Mendez	
IOL Manipulator	61	Degree Gauge	87
Fixation Forcep	140	Multipurpose LASIK Forcep	144
Lester-Burch		Meyerhoefer	
Eye Speculum	4	Chalazion Curette	39
Lewicky		Micro	
Lens Manipulating Hook	75	Viscocanalostomy Knife	15
		Scissor	106

Alphabetical Index

Pages	
M	
Micro-Adson	
Dressing Forcep	124
Dissecting Forcep	125
Micro-Mosquito	
Haemostatic Forcep	163
MICS	
Diamond Knife	24
Minardi	
Phaco Chopper	49
Mini-Adson	
Dressing Forcep	123
Dissecting Scissor	125
Moody	
Fixation Forcep	141
Moorfield	
Conjunctival Forcep	149
Moran	
LASIK Enhancement Spatula	57
Morrison	
Lens Loop and Probe	65
Muldoon	
Lacrimal Dilator	91
Murdoch	
Eye Speculum	3
N	
Nagahara	
Nucleus Lens Spatula & Chopper	52
Nucleus Chopper	52
Neuhann	
Nucleus Divider	52
New Orleans	
Lens Loop	64
Nightingale	
Capsule Polishing Curette	40
Noyes	
Iris Scissor	104-105
Nugent	
Utility Forcep	150
O	
Obstbaum	
Synechia Spatula	57
Osher	
Hook	78
Iris Forcep	151
IOL Cutting Scissor	106
O'Gawa-Utrata	
Capsulorrhesis Forcep	145
P	
Paton	
Spatula & Spoon	69
Paufique	
Graft Knife	14
Suture Tying Forcep	138
Pierce-Colibri	
Corneal Forcep	131
Pisacano	
Nucleus Rotator	55
Puttermann-Muellers	
Muscle Conjunctival Resection Ptosis Clamp	159
R	
Ragnell	
Dissecting Scissor	95
Rainin	
Air Injection Cannula	29
Randolph	
Cyclodialysis Cannula	29
Rentsch	
Boat Hook	78
Ribbon	
Scissor	102-103
Rosen	
Phaco Splitter	53
Rosenwasser	
DALK Cannula	28
DALK Corneal Dissector	43
Lamella Donor Shovel	70
Donor Lamella Inserting Forcep	153
Rowen	
LASIK Flap Protector/Spatula	41
Phaco Spatula 53	
S	
Salyer	
Ribbon Scissor	102
Sauer	
Fixation Forcep	138
Saunders	
Eye Speculum	5
Sauter	
Orbital Spatula	8
Hydrodissection Cannula	33
Schepens	
Scleral Depressor	7
Schoiotz	
Tonometer	90
Schocket	
Scleral Depressor	7
Scleral	
Shield 8	
Plug Removing Forcep	150
Seeley	
Hydrodissection Cannula	32
Sheets	
Irrigation Vectis	36
Shepard	
Iris Hook	74
Lens Holding Forcep	156
Shepherd	
Capsule Polishing Curette	40
Siepser	
Endocapsular Controller	55
Silverman	
Fixation Hook	81
Simcoe	
Irrigating - Aspirating Cannula	31
Cortex Extractor	35
Irrigation Lens Nucleus Loop	37
Simcoe (Reversed)	
Irrigating - Aspirating Cannula	31
Sinskey	
Lens Manipulating Hook	75-76
Sinskey (Reversed)	
Lens Manipulating Hook	74
Sinskey II	
Lens Manipulating Hook	77
Sloane	
LASIK Micro Hoe	42
Epi Peeler	42
LASIK Micro Hoe & Epi Peeler	42

Alphabetical Index

S	Pages	
LASIK Epi-Dissector	43	
LASIK Flap Repositor	43	
St. Martin's		
Suture Tying Forcep	138	
Stahl		
Trabeculotomy Probe	90	
Steinert		
IOL Inserting Forcep	154	
Stern-Castroviejo		
Suture Tying Forcep	136	
Stern-Gills		
Micro Scissor	119	
Stevens		
Sub Tenon's Anesthesia Cannula	30	
Tenotomy Hook	83	
Tenotomy Scissor	97-98, 103, 122	
Stevenson		
Lacrimal Sac Retractor	10	
Strabismus		
Scissor	99, 102	
T		
Tanna-Lieberman		
Eye Speculum	4	
Tennant		
Nucleus Rotator	56	
Suture Tying Forcep	146-147	
Thornton		
Fixation Ring	83	
Tibolt		
Punctal Scissor	113	
Tooke		
Corneal Knife	12	
Troutman		
Olive Tip Cannula	29	
Superior Rectus Forcep	140	
Micro Needle Holder	174	
Troutman-Castroviejo		
Corneoscleral Scissor	109	
Troutman-Katzin		
Corneal Transplant Scissor	110	
Troutman-Chris		
Blade Holder/Breaker	94	
U		
Utrata		
Capsulorrhesis Forcep	144	
V		
Vannas		
Capsulotomy Scissor	115-116	
Micro Scissor	122	
Vitrectomy		
Scissor	120	
Forcep	127	
Viscoelastic		
Aspirating Cannula	29	

W	Pages	
Walter		
Cataract Needle	13	
Walton		
Conjunctiva Retractor	9	
Foreign Body Gouge	16	
Watzke		
Sleeve Spreading Forcep	159	
Wecker		
Iris Scissor	107	
Wells		
Encucleation Spoon	72	
Westcott		
Tenotomy Scissor	112-113	
Stitch Scissor	114	
Wheeler		
Spatula	45	
Wies		
Chalazion Forcep	161	
Williams		
Eye Speculum	3	
Williamson-Noble		
Micro Scissor	104	
Wills		
Cautery	9	
Wright		
Fascia Needle	9	
Z		
Ziegler		
Cataract Knife	13	

Numerical Index

Cat. No.	Page	Cat. No.	Page	Cat. No.	Page	Cat. No.	Page
OP-010-10	1	OP-211-11	9	OP-359-22	24	OP-463-22	32
OP-010-15	1	OP-220-00	10	OP-370-10	25	OP-465-23	32
OP-012-15	1	OP-222-00	10	OP-370-15	25	OP-467-23	32
OP-020-10	1	OP-226-00	10	OP-370-20	25	OP-467-25	32
OP-020-15	1	OP-228-01	11	OP-370-90	25	OP-471-19	32
OP-022-20	1	OP-228-02	11	OP-381-08	25	OP-467-25	32
OP-030-13	1	OP-239-40	11	OP-381-09	25	OP-477-26	32
OP-030-15	1	OP-239-70	11	OP-381-10	25	OP-481-23	33
OP-032-15	1	OP-240-01	11	OP-385-08	25	OP-481-27	33
OP-040-10	2	OP-240-02	11	OP-385-09	25	OP-481-30	33
OP-042-00	2	OP-240-03	11	OP-385-10	25	OP-485-27	33
OP-046-00	2	OP-241-01	11	OP-391-00	26	OP-487-19	33
OP-050-09	2	OP-241-02	11	OP-395-00	26	OP-489-23	33
OP-060-08	2	OP-241-03	11	OP-397-00	26	OP-491-23	33
OP-060-11	2	OP-250-11	12	OP-399-00	26	OP-493-21	33
OP-060-15	2	OP-251-01	12	OP-400-01	27	OP-495-26	34
OP-064-12	3	OP-251-02	12	OP-400-02	27	OP-501-23	34
OP-064-15	3	OP-253-12	12	OP-400-03	27	OP-503-25	34
OP-070-15	3	OP-255-12	12	OP-400-04	27	OP-507-25	34
OP-072-02	3	OP-256-01	13	OP-400-05	27	OP-509-23	35
OP-074-11	3	OP-256-02	13	OP-400-06	27	OP-509-25	35
OP-074-15	3	OP-256-03	13	OP-400-07	27	OP-511-25	35
OP-076-16	4	OP-256-04	13	OP-400-08	27	OP-516-23	35
OP-078-00	4	OP-256-05	13	OP-400-09	27	OP-516-26	35
OP-081-15	4	OP-257-05	13	OP-400-10	27	OP-516-30	35
OP-083-10	4	OP-257-06	13	OP-400-11	27	OP-517-23	35
OP-083-15	4	OP-257-07	13	OP-400-12	27	OP-517-26	35
OP-083-15 A	4	OP-258-12	13	OP-400-13	27	OP-517-30	35
OP-085-15	4	OP-259-12	13	OP-400-14	27	OP-519-30	35
OP-087-10	5	OP-260-12	14	OP-401-01	27	OP-521-23	35
OP-087-15	5	OP-261-12	14	OP-401-02	27	OP-523-23	35
OP-089-15	5	OP-263-12	14	OP-401-03	27	OP-525-21	35
OP-089-15 A	5	OP-265-12	14	OP-401-04	27	OP-527-23	36
OP-091-14	5	OP-266-11	15	OP-401-05	27	OP-527-25	36
OP-091-16	5	OP-268-01	15	OP-401-06	27	OP-529-05	36
OP-095-04	5	OP-271-12	16	OP-401-07	27	OP-529-07	36
OP-100-00	6	OP-273-12	16	OP-401-08	27	OP-529-08	36
OP-100-14	6	OP-275-11	16	OP-401-09	27	OP-529-10	36
OP-100-15	6	OP-276-12	16	OP-401-10	27	OP-531-23	36
OP-100-16	6	OP-280-18	17	OP-401-11	27	OP-532-27	36
OP-100-17	6	OP-283-11	17	OP-401-12	27	OP-533-23	36
OP-100-18	6	OP-285-00	17	OP-401-13	27	OP-535-20	36
OP-100-19	6	OP-291-11	17	OP-410-23	28	OP-537-23	37
OP-100-20	6	OP-300-12	18	OP-410-27	28	OP-539-25	37
OP-100-21	6	OP-302-12	18	OP-411-23	28	OP-541-25	37
OP-100-22	6	OP-304-12	18	OP-411-27	28	OP-542-25	37
OP-111-15	6	OP-305-11	19	OP-412-23	28	OP-542-40	37
OP-111-16	6	OP-311-01	19	OP-413-23	28	OP-543-22	37
OP-111-17	6	OP-311-02	19	OP-415-23	28	OP-545-23	38
OP-151-01	7	OP-311-03	19	OP-417-27	28	OP-545-25	38
OP-153-01	7	OP-316-01	20	OP-419-27	28	OP-547-22	38
OP-153-02	7	OP-316-02	20	OP-421-27	28	OP-549-01 T	38
OP-153-03	7	OP-316-03	20	OP-421-30	28	OP-549-02 T	38
OP-160-14	7	OP-316-04	20	OP-423-27	28	OP-550-00	39
OP-163-00	7	OP-316-05	20	OP-423-30	28	OP-550-01	39
OP-171-15	7	OP-320-01	21	OP-427-19	29	OP-550-02	39
OP-172-01	8	OP-320-02	21	OP-427-21	29	OP-550-03	39
OP-172-02	8	OP-320-03	21	OP-429-19	29	OP-550-04	39
OP-172-03	8	OP-331-25	21	OP-429-22	29	OP-552-10	39
OP-173-13	8	OP-340-28	22	OP-433-25	29	OP-552-15	39
OP-173-15	8	OP-340-30	22	OP-437-27	29	OP-552-20	39
OP-175-11	8	OP-340-32	22	OP-437-30	29	OP-552-25	39
OP-175-17	8	OP-341-28	22	OP-439-25	30	OP-552-30	39
OP-175-15	8	OP-341-30	22	OP-441-19	30	OP-552-35	39
OP-175-17	8	OP-341-32	22	OP-443-23	30	OP-571-12	40
OP-180-01	8	OP-351-25	23	OP-445-23	30	OP-573-12	40
OP-180-02	8	OP-351-28	23	OP-447-24	30	OP-587-13	40
OP-192-11	9	OP-351-30	23	OP-451-23	31	OP-591-12	40
OP-197-14	9	OP-351-32	23	OP-453-23	31	OP-594-11	41
OP-199-14	9	OP-355-27	23	OP-455-23	31	OP-597-11	41
OP-201-12	9	OP-357-18	24	OP-457-23	31	OP-599-13	41
OP-211-07	9	OP-357-20	24	OP-459-23	31	OP-611-12	42
OP-211-09	9	OP-359-18	24	OP-462-22	31	OP-615-11	42

Numerical Index

Cat. No.	Page	Cat. No.	Page	Cat. No.	Page	Cat. No.	Page
OP-619-13	42	OP-735-12	62	OP-865-13	82	SL-100-11	93
OP-623-11	43	OP-738-12	62	OP-865-16	82	SL-100-15	93
OP-627-12	43	OP-739-11	62	OP-867-10	83	SL-206-10	93
OP-629-01	43	OP-741-13	62	OP-869-11	83	SL-206-13	93
OP-629-02	43	OP-743-13	63	OP-871-12	83	SL-206-15	93
OP-629-03	43	OP-745-25	63	OP-872-12	83	SL-210-61	93
OP-629-14	44	OP-745-50	63	OP-873-12	84	SL-210-62	93
OP-630-01	44	OP-749-12	63	OP-875-13	84	SL-210-63	93
OP-630-02	44	OP-750-01	64	OP-877-01	84	SL-210-64	93
OP-632-13	45	OP-750-02	64	OP-877-02	84	SL-210-65	93
OP-635-01	45	OP-752-13	64	OP-879-01	85	SL-210-66	93
OP-635-02	45	OP-754-13	64	OP-879-02	85	SL-210-67	93
OP-637-12	45	OP-757-15	65	OP-879-03	85	SL-210-68	93
OP-639-12	46	OP-758-14	65	OP-881-01	85	SL-210-69	93
OP-640-12	46	OP-760-13	65	OP-881-02	85	SL-900-09	94
OP-643-01	46	OP-761-13	66	OP-883-13	86	SL-901-09	94
OP-643-02	46	OP-765-10	66	OP-885-13	86	SL-910-11	94
OP-645-11	47	OP-767-10	66	OP-887-11	86	SL-916-13	94
OP-647-12	47	OP-769-11	67	OP-891-11	86	SS-243-12	95
OP-649-01	47	OP-771-11	67	OP-900-00	86	SS-243-15	95
OP-649-02	47	OP-773-12	67	OP-901-11	87	SS-400-13	95
OP-651-01	48	OP-775-02	68	OP-905-11	87	SS-400-13 SC	95
OP-651-02	48	OP-775-03	68	OP-907-11	87	SS-401-13	95
OP-653-01	48	OP-779-12	68	OP-911-01	88	SS-401-13 SC	95
OP-653-02	48	OP-781-12	68	OP-911-02	88	SS-430-09	95
OP-657-13	49	OP-783-12	69	OP-911-03	88	SS-430-10	95
OP-659-01	49	OP-790-13	69	OP-911-04	88	SS-431-09	95
OP-659-02	49	OP-793-14	69	OP-911-05	88	SS-431-10	95
OP-661-11	49	OP-795-04	70	OP-911-06	88	SS-434-09	96
OP-663-13	50	OP-795-05	70	OP-911-07	88	SS-435-09	96
OP-665-11	50	OP-796-00	71	OP-911-08	88	SS-440-11	96
OP-667-11	50	OP-798-00	71	OP-911-09	88	SS-441-11	96
OP-669-12	51	OP-798-90	71	OP-911-10	88	SS-440-11 SC	96
OP-670-15	51	OP-799-13	71	OP-911-11	88	SS-441-11 SC	96
OP-670-20	51	OP-800-01	72	OP-911-12	88	SS-446-12	97
OP-671-15	51	OP-800-02	72	OP-911-13	88	SS-447-12	97
OP-671-20	51	OP-802-15	72	OP-913-01	88	SS-450-10	97
OP-673-12	52	OP-803-00	72	OP-913-02	88	SS-450-11	97
OP-675-12	52	OP-804-12	73	OP-913-03	88	SS-451-10	97
OP-677-11	52	OP-805-12	73	OP-913-04	88	SS-451-11	97
OP-678-11	52	OP-811-11	73	OP-913-05	88	SS-450-11 SC	97
OP-679-11	52	OP-813-11	74	OP-913-06	88	SS-451-11 SC	97
OP-681-11	53	OP-815-12	74	OP-913-07	88	SS-452-10	98
OP-683-11	53	OP-817-12	74	OP-913-08	88	SS-452-11	98
OP-685-11	53	OP-819-01	75	OP-913-09	88	SS-452-12	98
OP-687-13	53	OP-819-02	75	OP-913-10	88	SS-452-11 SC	98
OP-689-01	54	OP-821-01	75	OP-913-11	88	SS-453-10	98
OP-689-02	54	OP-821-02	75	OP-913-12	88	SS-453-11	98
OP-691-11	54	OP-823-12	75	OP-913-13	88	SS-453-12	98
OP-693-11	55	OP-825-11	76	OP-915-08	89	SS-456-11	99
OP-695-11	55	OP-826-01	76	OP-915-12	89	SS-457-11	99
OP-697-11	55	OP-826-02	76	OP-915-16	89	SS-456-11 SC	99
OP-699-11	55	OP-826-03	77	OP-917-10	89	SS-457-11 SC	99
OP-700-12	56	OP-826-04	77	OP-919-09	89	SS-461-11	100
OP-701-12	56	OP-827-01	77	OP-930-12	89	SS-466-11	100
OP-703-12	56	OP-827-02	77	OP-932-01	90	SS-469-11	101
OP-705-11	56	OP-829-01	78	OP-932-02	90	SS-484-09	102
OP-707-11	57	OP-831-12	78	OP-934-30	90	SS-486-09	102
OP-709-12	57	OP-833-12	78	OP-934-35	90	SS-487-09	102
OP-711-13	57	OP-835-01	79	OP-935-01	90	SS-488-10	102
OP-713-14	58	OP-835-02	79	OP-941-01	90	SS-489-10	102
OP-715-13	58	OP-837-01	79	OP-941-02	90	SS-492-10	102
OP-717-12	58	OP-837-02	79	OP-944-08	91	SS-493-10	102
OP-719-12	59	OP-839-02	80	OP-946-10	91	SS-496-10	103
OP-720-13	59	OP-841-01	80	OP-948-01	91	SS-497-10	103
OP-723-12	59	OP-843-02	80	OP-948-02	91	SS-498-10	103
OP-725-11	60	OP-845-01	81	OP-948-03	91	SS-499-10	103
OP-726-11	60	OP-845-02	81	OP-950-14	92	SS-500-10	103
OP-727-12	60	OP-849-01	81	OP-953-14	92	SS-501-10	103
OP-728-12	61	OP-849-02	81	OP-956-14	92	SS-507-11	103
OP-729-12	61	OP-851-09	81	OP-960-00	92	SS-626-01	104
OP-731-11	61	OP-861-11	82	SL-010-03	93	SS-630-12	104
OP-733-12	61	OP-863-11	82	SL-012-03	93	SS-631-12	104

Numerical Index

Cat. No.	Page	Cat. No.	Page	Cat. No.	Page	Cat. No.	Page
SS-632-12	105	SS-972-10	122	FR-765-09	138	FR-900-10	157
SS-633-12	105	SS-973-10	122	FR-766-09	138	FR-901-01	157
SS-634-12	105	SS-974-10	122	FR-768-09	139	FR-901-02	157
SS-635-12	105	SS-975-10	122	FR-770-11	139	FR-903-01	158
SS-636-12	106	FR-050-12	123	FR-772-01	139	FR-903-02	158
SS-637-12	106	FR-050-15	123	FR-772-02	139	FR-905-20	158
SS-638-10	106	FR-050-18	123	FR-774-10	140	FR-905-27	158
SS-639-10	106	FR-052-12	123	FR-774-11	140	FR-906-10	159
SS-641-11	107	FR-054-12	124	FR-776-01	140	FR-907-11	159
SS-643-11	107	FR-054-15	124	FR-776-02	140	FR-908-10	159
SS-645-11	107	FR-250-12	124	FR-777-01	141	FR-910-09	160
SS-647-05	107	FR-250-15	124	FR-777-02	141	FR-912-08	160
SS-647-07	107	FR-250-18	124	FR-780-10	141	FR-912-12	160
SS-647-10	107	FR-252-12	125	FR-780-10 T	141	FR-914-09	160
SS-649-01	108	FR-253-12	125	FR-781-10	142	FR-916-09	161
SS-649-02	108	FR-256-12	125	FR-781-10 T	142	FR-918-09	161
SS-651-01	108	FR-256-15	125	FR-784-10	143	FR-919-09	161
SS-651-02	108	FR-630-09	126	FR-784-10 T	143	FR-920-09	161
SS-653-01	108	FR-632-10	126	FR-790-10	143	FR-921-09	161
SS-653-02	108	FR-680-20	127	FR-793-10	143	FR-922-01	162
SS-655-01	109	FR-680-23	127	FR-795-10	144	FR-923-01	162
SS-655-02	109	FR-682-20	127	FR-797-10	144	FR-923-02	162
SS-657-01	109	FR-682-23	127	FR-801-10	144	FR-925-10	162
SS-657-02	109	FR-685-20	127	FR-801-10 T	144	AF-062-03	163
SS-659-01	110	FR-687-20	127	FR-803-10	145	AF-063-03	163
SS-659-02	110	FR-689-23	127	FR-805-10	145	AF-140-10	163
SS-661-01	110	FR-690-23	127	FR-807-11	145	AF-141-10	163
SS-661-02	110	FR-692-20	128	FR-809-10	145	AF-160-12	163
SS-663-11	111	FR-695-23	128	FR-810-10	146	AF-161-12	163
SS-665-11	111	FR-700-01	129	FR-811-10	146	AF-170-10	163
SS-667-11	111	FR-700-02	129	FR-816-10	146	AF-171-10	163
SS-669-11	111	FR-700-03	129	FR-816-10 T	146	AF-172-10	163
SS-671-11	112	FR-702-01	129	FR-817-10	147	AF-173-10	163
SS-673-11	112	FR-702-02	129	FR-817-10 T	147	AF-761-05	164
SS-675-11	112	FR-702-03	129	FR-820-10	147	AF-761-09	164
SS-677-01	113	FR-710-07	130	FR-821-10	147	AF-773-08	164
SS-677-02	113	FR-712-07	130	FR-824-10	148	OT-751-02	165
SS-679-11	113	FR-714-07	130	FR-825-10	148	OT-751-03	165
SS-681-11	113	FR-716-07	130	FR-827-08	148	RH-305-03	165
SS-683-11	114	FR-718-07	130	R-827-08 T	148	RH-305-04	165
SS-685-11	114	FR-718-07 T	130	FR-830-10	149	RH-305-05	165
SS-688-11	114	FR-720-07	131	FR-834-12	149	RH-305-06	165
SS-689-11	114	FR-721-07	131	FR-834-17	149	RH-033-01	166
SS-690-11	115	FR-722-08	131	FR-836-10	149	RH-033-02	166
SS-691-11	115	FR-723-08	131	FR-838-10	149	RH-033-03	166
SS-694-08	115	FR-724-10	132	FR-840-10	150	RH-475-18	166
SS-695-08	115	FR-726-11	132	FR-842-10	150	RH-600-16	166
SS-696-08	116	FR-728-11	132	FR-844-11	150	RH-601-16	166
SS-697-08	116	FR-730-08	133	FR-847-08	150	RH-604-16	166
SS-698-08	116	FR-732-10	133	FR-846-10	150	RH-605-16	166
SS-699-08	116	FR-732-10 T	133	FR-848-11	151	SU-040-13	167
SS-701-08	117	FR-734-08	133	FR-850-11	151	SU-200-10	167
SS-703-08	117	FR-736-10	134	FR-855-11	151	SU-201-10	167
SS-705-08	117	FR-737-10	134	FR-857-11	151	SU-202-10	167
SS-706-10	118	FR-740-10	134	FR-859-11	151	SU-203-10	167
SS-707-10	118	FR-742-10	135	FR-861-11	152	SU-210-14	168
SS-709-10	118	FR-750-01	135	FR-863-11	152	SU-211-14	168
SS-710-11	119	FR-750-03	135	FR-867-11	152	SU-212-14	168
SS-711-11	119	FR-750-05	135	FR-869-10	152	SU-213-14	168
SS-713-11	119	FR-750-09	135	FR-871-11	153	SU-220-14	169
SS-715-11	119	FR-750-01 T	135	FR-875-11	153	SU-221-14	169
SS-720-20	120	FR-751-01	136	FR-877-10	153	SU-222-14	169
SS-721-20	120	FR-751-03	136	FR-879-10	154	SU-223-14	169
SS-721-23	120	FR-751-05	136	FR-882-10	154	SU-300-10	170
SS-721-30	120	FR-751-09	136	FR-883-10	154	SU-301-10	170
SS-721-45	120	FR-751-01 T	136	FR-885-10	155	SU-310-10	170
SS-723-23	120	FR-752-05	136	FR-887-10	155	SU-310-11	170
SS-723-60	120	FR-754-10	137	FR-889-10	155	SU-311-10	170
SS-723-80	120	FR-756-10	137	FR-891-12	156	SU-311-11	170
SS-725-23	121	FR-758-08	137	FR-893-12	156	SU-313-11	171
SS-726-12	122	FR-760-08	137	FR-895-11	156	SU-315-11	171
SS-727-12	122	FR-762-09	138	FR-897-12	156	SU-315-11 T	171
SS-729-12	122	FR-764-09	138	FR-899-12	157	SU-320-13	172

Numerical Index

Cat. No.	Page
SU-321-13	172
SU-322-13	172
SU-323-13	172
SU-331-14	173
SU-331-14 T	173
SU-333-14	173
SU-333-14 T	174
SU-341-12	174
SU-343-12	174
SU-500-13	175
SU-502-13	175
SU-510-14	176
SU-511-14	176
SU-512-14	177
SU-513-14	177
DI-850-20	178
DI-880-15	178
RT-048-01	178
RT-048-02	178
RT-048-03	178
RT-066-13	178
RT-068-13	178
RT-074-02	178
RT-074-05	178
RT-074-07	178
RT-074-10	178
RT-120-08	178
RT-120-10	178
RT-120-12	178
RT-120-14	178
OR-470-19	179
OR-933-01	179
OR-934-02	179
OR-937-16	179
PR-200-00	179
ST-110-01	180
ST-110-02	180
ST-110-03	180
ST-112-26	180



COPYRIGHT NOTICE

The copyright, trademarks, domain names and other intellectual property rights in all materials and information in this presentation belongs to Surtex Instruments Ltd. All rights are reserved. Technical alterations are possible. This presentation may be used for no other purpose than offering, buying and selling of Surtex Instruments products. No part may be reproduced or copied in any form. In the event of misuse we retain the right to recall our catalogues and take appropriate legal action.

DISCLAIMER

This document provides general information about Surtex Instruments Ltd, its products and services, and summarises general capabilities and offerings, which we deem to be relevant to our customers. Whilst Surtex Instruments has taken appropriate steps to ensure the accuracy of the information contained in the document Surtex Instruments gives no warranty regarding the accuracy or completeness of such information.

The presentation was prepared to provide information and does not constitute a contract. Surtex reserves the right, without prior notice, to change, delete, supplement, or otherwise amend at any time the information, images and offerings contained in the presentation or other documents. Surtex shall not be liable for any loss or damage of whatever nature (whether direct, indirect, consequential or other) may arise as a result of any third party relying on the information contained in this document.

For further detailed information, prices, terms and conditions, customer should contact Surtex Instruments at sales@surtex-instruments.co.uk.

© 2021 SURTEX INSTRUMENTS. ALL RIGHTS RESERVED. PRINTED IN UNITED KINGDOM. BT1-01 Rev 1.0



Unit 202 Kingspark Business Centre
152-178 Kingston Road, New Malden
Surrey KT3 3ST United Kingdom
Phone: +44 (0) 208 123 3513

E-mail: sales@surtex-instruments.co.uk
Website: www.surtex-instruments.co.uk

