




SURTEX
Maxview
RETRACTOR SYSTEM



SURTEX
INSTRUMENTS

Crafted to Excellence



SURTEX
INSTRUMENTS

Crafted to Excellence

SURTEX
Maxview
RETRACTOR SYSTEM

Surtex Maxview Retractor System

Our **Maxview** Retractor is an interdisciplinary surgical retractor system:

The system is modular in construction and retaining arms can be positioned according to specified operating field.

Maxview Blade Holder can be fixed to any area of retaining system.


Because of our innovative blade holding locking system, it's possible to change the blade - without dismantling the holder itself - within a few seconds (Press-Turn-Locking/PTL).

Same blade can be user-defined utilized as swivel blade or as fixated blade by turning the rotating wheel.

The self retaining modular construction makes the Surtex **Maxview** Retractor system efficient and optimally cost-effective to use.

Areas of use and access:

- Abdominal
- Thyroid/Carotid
- Axilla
- Upper abdomen/Rochard
- Vaginal
- Perineal
- Spine
- Paediatric Surgery
- Transplant Surgery
- Endoprosthetic






Socket Clamp

Catalog #	Size
RT-1000-01	25 x 10 mm (Euro Rail)
RT-1000-02	28.5 x 9.5 mm (US Rail)



Socket Clamp Universal

Catalog #	Size
RT-1000-03	Standard



Pillar

Catalog #	Size
RT-1070-35	35 cm, incl. hole Ø 12 mm



Distance Holder

Catalog #	Size
RT-1072-00	22 cm, Ø 18 mm, 25x10 mm



Support Arm

Catalog #	Size, Profile
RT-1076-43	43 cm, Straight
RT-1076-50	50 cm, Straight





Crossbar

Catalog #	Size
RT-1080-22	22 cm



Mini Frame Arm

Catalog #	Size, Profile
RT-1084-30	30 cm, Angled 2 Fold



Ortho Frame Arm

Catalog #	Size, Profile
RT-1084-40	40 cm, Angled 2 Fold



Maxi Frame Arm

Catalog #	Size, Profile
RT-1084-50	50 cm, Angled 2 Fold



Stack

Catalog #	Size
RT-1086-15	15 cm
RT-1086-45	45 cm






Extension Arm

Catalog #	Size, Profile
RT-1088-35	Ø 12 mm, 2X, 35 cm
RT-1088-48	Ø 12 mm, 2X, 48 cm



P-Clamp

Catalog #	Size, Profile
RT-1090-00	Two Diameter (Ø 16 + 12 mm)




S-Clamp

Catalog #	Size, Profile
RT-1094-00	Two Diameter (Ø 12 + 8 mm)



Q-Clamp

Catalog #	Profile
RT-1098-00	With Lock



Blade Holder

Catalog #	Size, Profile
RT-1106-20	20 cm, Press-Turn-Locking





Blade Holder Swivel
With Toothed Bar

Catalog #	Size
RT-1112-21	21 cm



Adaptor For Blade Holder
With Toothed Bar

Catalog #	Size
RT-1100-00	Standard



Rochard Traction Device

Catalog #	Size
RT-1102-22	22 cm



Holder W. Ball and Socket Joint

Catalog #	Size
RT-1114-23	23 cm Ø 8 mm





Extension Arm

Catalog #	Size, Profile
RT-1088-48	Ø 12 mm, 2 x 48 cm



Holder W.Hexagonal Snap Closing

Catalog #	Size
RT-1110-20	20 cm



Holder With Ball & Socket Joint

Catalog #	Size
RT-1120-00	Standard



Fixing Adapter W. Ball F.Hohmann-Retractor

Catalog #	Size
RT-1122-00	Standard



Fixing Adapter For Ball F.Hohmann-Retractor

Catalog #	Size
RT-1124-00	Standard





Ortho Frame Arm
Angled



Catalog #	Size, Profile
RT-1084-40	40 cm, 2 Fold

Holder W.Hexagonal



Catalog #	Size, Profile
RT-1110-20	20 cm, Snap Closing

Fixing Adapter For Ball/Hohmann Retractor




Catalog #	Size
RT-1124-00	Standard

Handle For Blade Holder



Catalog #	Size
RT-1108-00	Standard

Optoclip
Endoscope Holder



Catalog #	Size
RT-1002-00	Standard



Bushing For Endoscopes

Catalog #	Size
RT-1002-04	Ø 4 mm
RT-1002-05	Ø 5 mm
RT-1002-08	Ø 8 mm
RT-1002-10	Ø 10 mm
RT-1002-12	Ø 12 mm



Kirshner Body Wall Retractor Blade

Catalog #	Size
RT-1010-40	55 mm x 40 mm
RT-1010-50	65 mm x 50 mm
RT-1010-60	85 mm x 60 mm
RT-1010-78	65 mm x 78 mm
RT-1010-80	90 mm x 80 mm





Langenbeck Retractor Blade

Catalog #	Size
RT-1020-30	11 mm x 30 mm
RT-1020-35	11 mm x 35 mm
RT-1020-40	11 mm x 40 mm
RT-1020-45	11 mm x 45 mm
RT-1020-50	11 mm x 50 mm
RT-1020-55	11 mm x 55 mm
RT-1020-60	11 mm x 60 mm
RT1020-65	11 mm x 65 mm
RT-1022-04	20 mm x 40 mm
RT-1022-06	20 mm x 60 mm
RT-1022-08	20 mm x 80 mm
RT-1022-10	20 mm x 100 mm
RT-1022-12	20 mm x 120 mm
RT-1022-14	20 mm x 140 mm
RT-1022-16	20 mm x 160 mm
RT-1024-09	40 mm x 90 mm
RT-1024-10	40 mm x 100 mm
RT-1024-12	40 mm x 125 mm
RT-1024-15	40 mm x 150 mm
RT-1026-11	50 mm x 115 mm
RT-1026-15	50 mm x 150 mm
RT-1026-20	50 mm x 200 mm



Langenbeck Retractor Blade

Catalog #	Size
RT-1022-04 T	20 mm x 40 mm
RT-1022-06 T	20 mm x 60 mm
RT-1022-08 T	20 mm x 80 mm
RT-1022-10 T	20 mm x 100 mm
RT-1022-12 T	20 mm x 120 mm
RT-1022-14 T	20 mm x 140 mm
RT-1022-16 T	20 mm x 160 mm





Mini-Langenbeck Retractor Blade

Catalog #	Size
RT-1014-25	6 mm x 25 mm



Retractor Blade, Malleable

Catalog #	Size
RT-1030-50	50 cm x 160 mm
RT-1030-75	75 cm x 160 mm



Deaver Retractor Blade

Catalog #	Size
RT-1036-07	25 mm x 75 mm
RT-1036-09	25 mm x 90 mm
RT-1036-10	25 mm x 105 mm
RT-1036-12	25 mm x 120 mm
RT-1036-15	25 mm x 150 mm
RT-1036-20	25 mm x 200 mm
RT-1038-12	50 mm x 125 mm
RT-1038-15	50 mm x 150 mm
RT-1040-11	70 mm x 115 mm
RT-1040-20	70 mm x 200 mm
RT-1042-10	75 mm x 100 mm



Slotted Retractor Blade Malleable


Catalog #	Size
RT-1044-12	50 mm x 125 mm
RT-1044-16	50 mm x 165 mm
RT-1046-12	80 mm x 125 mm
RT-1046-16	80 mm x 165 mm





Finger Hook
Adjustable Blade

Catalog #	Size
RT-1048-12	90 mm x 125 mm
RT-1048-14	90 mm x 140 mm
RT-1048-16	90 mm x 165 mm



Doyen Retractor Blade

Catalog #	Size
RT-1050-09	60 mm x 92 mm
RT-1050-12	60 mm x 122 mm
RT-1050-16	60 mm x 160 mm



Harrington-Mini Retractor
Malleable

Catalog #	Size
RT-1052-50	25 mm x 50 mm
RT-1052-70	25 mm x 70 mm



Harrington Retractor
Malleable

Catalog #	Size
RT-1054-12	45 mm x 125 mm
RT-1054-17	45 mm x 175 mm
RT-1056-12	65 mm x 125 mm
RT-1056-17	65 mm x 175 mm



**Breisky Specula**

Catalog #	Size
RT-1060-09	25 mm x 90 mm
RT-1060-11	25 mm x 110 mm
RT-1060-13	25 mm x 130 mm
RT-1060-15	25 mm x 150 mm
RT-1060-20	25 mm x 200 mm
RT-1060-25	25 mm x 250 mm
RT-1062-09	35 mm x 90 mm
RT-1062-11	35 mm x 110 mm
RT-1062-13	35 mm x 130 mm
RT-1062-15	35 mm x 150 mm
RT-1062-20	35 mm x 200 mm

**Fence**

Catalog #	Size
RT-1064-12	100 mm x 120 mm
RT-1064-15	100 mm x 150 mm
RT-1064-20	100 mm x 200 mm
RT-1064-25	100 mm X 250 mm

**Richardson Retractor Blade**

Catalog #	Size
RT-1066-18	15 mm x 18 mm
RT-1066-25	21 mm x 25 mm
RT-1066-30	29 mm x 30 mm
RT-1066-37	38 mm x 37 mm





Kelly Retractor Blade

Catalog #	Size
RT-1074-64	52 mm x 64 mm
RT-1074-80	64 mm x 80 mm



Roux Retractor Blade

Catalog #	Size
RT-1078-01	11 mm x 17 mm
RT-1078-02	11 mm x 38 mm
RT-1078-03	16 mm x 17 mm
RT-1078-04	18 mm x 22 mm
RT-1078-05	23 mm x 28 mm



Volkman Retractor Blade

Sharp

Catalog #	Size
RT-1082-02	8 mm x 12 mm, 2 Teeth
RT-1082-03	16 mm x 20 mm, 3 Teeth
RT-1082-04	22 mm x 20 mm, 4 Teeth
RT-1082-06	35 mm x 20 mm, 6 Teeth



Volkman Retractor Blade

Blunt

Catalog #	Size
RT-1083-02	8 mm x 12 mm, 2 Teeth
RT-1083-03	16 mm x 20 mm, 3 Teeth
RT-1083-04	22 mm x 20 mm, 4 Teeth
RT-1083-06	35 mm x 20 mm, 6 Teeth



**Israel Retractor Blade**

Blunt

Catalog #	Size
RT-1092-05	41 mm x 51 mm, 5 Teeth
RT-1092-06	62 mm x 53 mm, 6 Teeth

**Prostate Retractor**

Malleable

Catalog #	Size
RT-1096-50	50 mm x 125 mm

**Bladder-Retractor Blade**

Malleable

Catalog #	Size
RT-1096-60	60 mm x 180 mm
RT-1096-60	80 mm x 150 mm

**Bladder-Retractor Blade**

Malleable

Catalog #	Size
RT-1104-25	25 mm x 100 mm

**St. Marks Pelvis Retractor**

(Kidney Vein)

Catalog #	Size
RT-1116-17	60/25 mm x 174 mm



St. Marks Pelvis Retractor Blade

Catalog #	Size
RT-1118-17	60/45 mm x 178 mm



St. Marks Pelvis Retractor Blade

Catalog #	Size
RT-1118-12	60/45 mm x 124 mm
RT-1118-27	60/45 mm x 270 mm



Retractor Blade For Axilla

Catalog #	Size
RT-1126-13	35 mm x 130 mm
RT-1126-17	25 mm x 175 mm



Alan Parks Retractor Blade

Catalog #	Size
RT-1128-75	24 mm x 75 mm
RT-1128-95	24 mm x 95 mm





Heavy Duty Rail Clamp

Catalog #	Size
RT-1150-00	Standard



Holding Arm Single With Hexagonal Quick Coupling

Catalog #	Size
RT-1152-00	Standard



Swan Neck Fiber Optic Light Cable

Catalog #	Size
RT-1160-40	40 cm, Ø 5 mm



Cushing Retractor Blade

Catalog #	Size
RT-1162-08	8 mm x 10 mm
RT-1162-10	10 mm x 11 mm
RT-1162-12	12 mm x 12 mm
RT-1162-14	14 mm x 13 mm
RT-1162-16	16 mm x 14 mm
RT-1162-18	18 mm x 16 mm



Frame Arm Angled 45°

Catalog #	Size
RT-1170-12	12 mm x 10 cm



Basket Including Inlay
For Basic Set


Catalog #	Size
RT-1199-01	Standard

Basket Including Inlay
For Blades & Accessories


Catalog #	Size
RT-1199-02	Standard

Blade Holder
Flex Swivel, With Spring


Catalog #	Size
RT-1172-00	Standard

Can be used for small incisions, such as resection of the thyroid gland, Head-Neck Dissection

Screw Terminal


Catalog #	Size
RT-1174-00	Standard



Frame Clamp

Catalog #	Size
RT-1176-00	Standard



Standard Instruments Clamp

Catalog #	Size
RT-1178-00	Standard



Isolated Optic/Endoscope Clamp

Catalog #	Size
RT-1180-00	Standard



Blade Holder With Ratchet Tip

Catalog #	Size
RT-1182-00	Standard



Liver Retractor Hook

Catalog #	Size, Profile
RT-1185-01	Ø 5 mm, Small
RT-1185-02	Ø 5 mm, Medium
RT-1185-03	Ø 5 mm, Large





COPYRIGHT NOTICE

The copyright, trademarks, domain names and other intellectual property rights in all materials and information in this presentation belongs to Surtex Instruments Ltd. All rights are reserved. Technical alterations are possible. This presentation may be used for no other purpose than offering, buying and selling of Surtex Instruments products. No part may be reproduced or copied in any form. In the event of misuse we retain the right to recall our catalogues and take appropriate legal action.

DISCLAIMER

This documents provides general information about Surtex Instruments Ltd, its products and services, and summarises general capabilities and offerings, which we deem to be relevant to our customers. Whilst Surtex Instruments has taken appropriate steps to ensure the accuracy of the information contained in the document Surtex Instruments gives no warranty regarding the accuracy or completeness of such information.

The presentation was prepared to provide information and does not constitute a contract. Surtex reserves the right, without prior notice, to change, delete, supplement, or otherwise amend at any time the information, images and offerings contained in the presentation or other documents. Surtex shall not be liable for any loss or damage of whatever nature (whether direct, indirect, consequential or other) may arise as a result of any third party relying on the information contained in this document.

For further detailed information, prices, terms and conditions, customer should contact Surtex Instruments at sales@surtex-instruments.co.uk.

© 2020 SURTEX INSTRUMENTS. ALL RIGHTS RESERVED. PRINTED IN UNITED KINGDOM. BT1-01 Rev 1.0



Unit 202 Kingspark Business Centre
152-178 Kingston Road, New Malden
Surrey KT3 3ST United Kingdom
Phone: +44 (0) 208 123 3513

E-mail: sales@surtex-instruments.co.uk
Website: www.surtex-instruments.co.uk



0005